



Lovibond® Colour Measurement

Tintometer® Group



**Understanding Colour
Communication**

www.lovibondcolour.com

A Guide to Understanding Colour Communication

Contents	Page
What is Colour?	1
Communication of Colour	1
How to Describe Colour	2
Hue	
Chroma	
Lightness	
Ways to Measure Colour	3
Reflected Light	4
Spherical	
0°:45° (or 45°:0°)	
Multi-Angle	
Non-Contact	
Tristimulus Colorimeter	
Transmitted Light	7
Regular	
Diffuse	
Visual Comparators vs. Automatic Measurements	8
Colour Consistency	9
Accuracy of Measurements	10

Correct Record Keeping 11

The Supply Chain 11

Methods of Quantifying Colour – Colour Scales and Spaces 12

- Lovibond® RYBN Colour
- Neutral Filters
- Colour Nomenclature
- The Munsell Scale
- CIE Colour Systems
- CIE XYZ
- Expressing Colour Uniformly
- CIELAB L*a*b*
- CIELAB L*C*h

Metamerism 18

Colour Differences, Notation and Tolerancing 19

- Delta CIELAB L*a*b* and CIELAB L*C*h
- Delta (Δ) E* colour difference
- CIE Colour Space Notation
- Visual Colour and Tolerancing
- Delta E* Tolerancing
- CIELAB L*a*b* Tolerancing
- CIELAB L*C*h Tolerancing
- CMC Tolerancing
- CIE2000 Tolerancing
- Choosing the Right Tolerance

Colour Scales 24

- ASTM Colour (ASTM D 1500, ASTM D 6045, ISO 2049, IP196)
- EBC (European Brewing Convention)
- European Pharmacopoeia (EP) Colour
- Gardner Colour (ASTM D 1544, ASTM D 6166, AOCS Td 1a, MS 817 Part 10)
- Platinum-Cobalt/Hazen/APHA Colour (ASTM D 1209)
- Saybolt Colour (ASTM D 156, ASTM 6045)
- White and Yellow Indices

What is Colour?

Colour is a visual, perceptual property in human beings.

Colour is derived from the signals produced by three different types of light sensitive cell in the eye that respond to the spectrum of light (distribution of light energy versus wavelength). In our environment, materials are coloured depending on the wavelengths of light they reflect or transmit.

The visible colour spectrum runs from blue through to red wavelengths, approximately 360-720 nm.

Three things are necessary to see colour:

- A source of light
- An object
- An observer/processor

Communication of Colour

How would you describe the colour of this rose? Would you say it's yellow, lemon yellow or maybe a bright canary yellow?

Each person verbally describes and hence defines the colour of an object differently.

As a result, objectively communicating a particular colour to someone without some type of physical standard is difficult. Also describing in words the precise colour difference between two objects is very challenging.

Your perceptions and interpretations of colour and colour comparisons are highly subjective. Eye fatigue, age and other physiological factors can influence your colour perception.

But even without such physical considerations, each observer interprets colour based on their personal perspective, feelings, beliefs and expectations. For example you may convince yourself that a certain colour match is within tolerance if you are under pressure to declare a colour match as acceptable.

The solution to this dilemma is an instrument that explicitly identifies a colour by measuring it and comparing the colour to standards completely objectively and accurately each and every time. That is, an instrument that differentiates a colour from all others and assigns it a numeric value.

Before we measure a colour we need to establish a means of describing it.






Figure 1: Hue




Figure 2: Chromaticity

How to describe Colour?

Colour is typically described utilising three dimensions: hue, chroma and lightness. By describing a colour using these three attributes, you can accurately identify a particular colour and distinguish it from any other.

Hue

When asked to identify the colour of an object, you'll most likely speak first of its hue. Quite simply, hue is how we perceive an object's colour – red, orange, green, blue, etc.

The colour wheel in Figure 1 shows the continuum of colour from one hue to the next. As the wheel illustrates, if you were to mix blue and green paints, you would get blue-green. Add yellow to green for yellow-green, and so on.

Chroma

Chroma describes the vividness or dullness of a colour – in other words, how close the colour is to either grey or the pure hue. For example, think of the appearance of a tomato and a radish. The red of the tomato is vivid, while the radish appears duller.

Figure 2 shows how chroma changes as we move from the centre to the perimeter. Colours in the centre appear grey (dull) and they become more saturated (vivid) as we move toward the perimeter. Chroma also is also known as saturation.

Lightness

The luminous intensity of a colour can be described by its lightness. Colours can be classified as light or dark when comparing their lightness.

For example, when a tomato and a radish are placed side by side, the red of the tomato appears to be much lighter. In contrast, the radish has a darker red lightness. In Figure 3, the lightness dimension is represented on the vertical axis.




Figure 3: Lightness

Ways to Measure Colour

Today, *colorimeters*, *spectrophotometers* and *spectrophotometric colorimeters* (sometimes referred to as *spectro-colorimeters*) are the most commonly used instruments for measuring colour worldwide.

These technologies measure the proportion of reflected or transmitted light at many points in the visible spectrum. The points can be plotted graphically to form a spectral curve. Since the spectral curve of each colour is completely unique, like a signature or fingerprint, the curve is an excellent tool for identifying, specifying and matching colour.

When an object interacts with light, some of the wavelengths of light are absorbed and others are reflected or transmitted (in the case of a coloured but clear liquid). Therefore a red ball absorbs all wavelengths of light except for those in the red part of the spectrum, which it reflects.

A blue ball reflects only blue wavelengths. A glass of apple juice transmits green and yellow wavelengths. Fresh snow reflects most of the light that interacts with it and hence appears white. Black is due to the absence of reflected light.




Figure 4: Spectral curve from a measured sample

Due to the numerous potential applications for colour measurement and the wide range of different types of materials which need to be tested, there are many diverse instrumental options which have been optimised for specific purposes.

The first consideration is whether the colour being measured is from light being reflected from a sample or light being transmitted through a sample.

Reflected Light

Reflected light can be measured in a number of different ways: different instrument geometries affect the reading obtained but can be easily matched to your application.

The different instrument geometries do, however, share common elements. Light from a controlled light source and a sensitive light detector are always utilised.

Sphere Geometry

Instruments with sphere geometry have played a major role in colour quality control systems for over 50 years. The sphere has a white inner surface and a small circular aperture, against which the sample is placed to allow colour measurement. Additional apertures allow light to enter and leave the sphere at defined angles.

Specular Component Included (SCI) and Specular Component Excluded (SCE)

The Specular Component is the component of light that is reflected from a surface at an angle equal to the incident angle of the illumination. A high Gloss surface will reflect more light into the Specular direction (i.e. act like a mirror) and appear smooth and

shiny. A low Gloss surface will reflect less light and appear matt.

The light that is not reflected in a Specular direction, but scattered in many directions, is called Diffuse reflectance.


If two samples of identically coloured plastic, that differ only in surface effect (i.e. one shiny, the other matt), are measured using the same instrument in SCI and SCE modes respectively, the results will differ as follows.

In SCE mode, the Specular reflectance is excluded from the measurement and only the diffuse reflectance is measured. This produces a colour evaluation which correlates to the way the observer sees the colour of an object.

In SCI mode, the Specular reflectance is included with the diffuse reflectance during the measurement process. This type of colour evaluation measures the total appearance independent of surface conditions.

Therefore the two samples of plastic should provide values that are very similar when measured in the SCI mode and values that show a colour difference in the SCE mode. The difference between the two indicates the effect of the gloss on appearance.

Most sphere-based instruments use an angle close to the perpendicular, usually 8° , to give the capability of including or excluding the Specular reflectance while making a measurement. By opening a small port in the sphere, the Specular component can



Sphere

be excluded from the measurement. This is because light that would be reflected at an equal angle (on the opposite side of the perpendicular) to the sample viewing port is lost when the Specular port is open and is therefore not included in any measurement. In some instruments the same effect is achieved by using a light trap instead of an open port.


$0^\circ:45^\circ$ (or $45^\circ:0^\circ$) Geometry

The angles in this type of instrument refer to the relative angles of the detector and illuminating light source in the instrument. The first angle is the angle of the incident light from the light source, and the second

is the angle of the detector – a 0° angle being perpendicular to the surface of the sample being measured.

The detector position can be at a single point in a plane at 45° , or a number of detectors can be positioned at discrete points around a circumference to approximate an annular ring.

No instrument "sees" colour more like the human eye than that with $0^\circ:45^\circ$ geometry. This is because a viewer always tries to exclude the Specular Component or Specular Reflectance when judging colour. When we look at pictures in a glossy magazine for example, we usually hold the magazine such that the specular light does not reflect back to the eye.




$0^\circ / 45^\circ$

A $0^\circ:45^\circ$ instrument will remove specular light from the measurement and measure the appearance of the sample as the human eye would see it.

Multi-Angle

In the past 20 years or so, car makers have experimented with special effect colours. They use special additives such as mica, pearlescent materials, ground-up seashells, microscopically coated coloured pigments and interference pigments to produce different colours at different angles of view.

Large and expensive Goniometers can be used to measure these colours although battery-powered, hand-held, multi-angle instruments are also available and are now used by most auto manufacturers and their supply chains worldwide.



Multi-angle

measured in their natural, unaltered state, as the eye sees the sample, by using non-contact colour measurement technology.

Non Contact spectrophotometers are designed for the colour measurement of many types of wet and dry samples including powders, pastes, gels, plastics and paints.

Tristimulus Colorimeter

Tristimulus Colorimeters are not the same as spectrophotometers. They are tristimulus (three-channel) devices that make use of red, green and blue filters to emulate the response of the human eye to light and colour.



NC45

Non-Contact

Products that normally require protection from physical contact with measurement apertures, such as liquids and pastes, or in which the surface appearance is changed by the presentation method, such as when the sample is pressed behind glass, can now be

Transmitted Light

Transmitted light can be measured by analyzing the light from a controlled light source that is passed through the sample using a suitable detector.

Regular Geometry

This is the most common transmittance geometry where the sample is placed into a parallel beam of light between the light source and the detector system.

For meaningful results, the sample should give a sharp image when looking through it at other objects. Any level of turbidity will seriously interfere with the accuracy of the colour readings using regular transmitted light.

Diffuse

This geometry uses a sphere. Either the sample is illuminated by a parallel beam of light with a sphere used to collect the light transmitted through the sample or the sphere is used to illuminate the sample diffusely with measurement of the light coming perpendicularly through the sample.

This is the way translucent samples and those that scatter light, giving a cloudy or fuzzy image of objects viewed through them, should be measured. For example; fruit juice is frequently measured using this method.



Visual Comparators vs. Automatic Measurements

Ensuring colour accuracy every time is critical to manufacture products of a consistent, high quality.

Reliable and repeatable colour test results are the key to ensuring final product quality and also to minimise production costs. Simplicity of operation helps to reduce error and increase productivity.

The fundamental difference between the Lovibond® Visual and Automatic ranges is that the Visual instruments are based on subjective, visual comparison methods (relying to a high degree on the judgment and skill of the operator and hence their perspective, feelings, beliefs and expectations) while the Automatic ranges rely on automatic, non-subjective measurements and are, therefore, completely unaffected by the judgment of the operator.

With operators inexperienced in this area, visual comparison can be more time consuming and less precise than the fast automatic instruments. Visual agreement between different operators at one site or multiple sites cannot be guaranteed.

The skill/experience of operators, degree of acceptable error, sample preparation time, scale choice and required scale resolution should be carefully considered before making a decision on which instrument to purchase.

In addition, although Visual systems are of a lower initial cost, their limitations should be taken into account when selecting the correct instrument.



EComparator series

Colour Consistency

When comparing Visual (subjective) to Automatic (non-subjective) colour assessment, the fundamental differences between these methods need to be considered.

Basic steps can be taken to reduce colour communication problems:

- Do you have a systematic, consistent and reliable method for sample preparation and presentation?

For example: When measuring liquids are you using comparable, clean, cells?



Figure 5 – Top View

Figure 5 shows the same liquid sample viewed across a range of cell path-lengths. Figure 6 shows that, as path-length changes, the perceived colour of the samples will change significantly. The results of any visual or automatic methods would be influenced by this difference.




Figure 6 – Front View

For example, with **Lovibond® RYBN Colour**, it is advisable that the depth of colour should never be greater than that which may be matched by a total of 20 Lovibond® units. This is because slight differences are most easily perceived in intensities ranging between approximately 3 and 10 units of the predominant colour. To clarify:

(i) using a 5.25" cell, a sample gives **R** 1.7, **Y** 8, **B** 0 and **N** 0.1 = 9.8 in total: the path-length used is therefore correct.

(ii) using a 5.25" cell a sample gives **R** 5.7, **Y** 31, **B** 0 and **N** 0.1 = 36.8 in total: the path-length used is therefore incorrect since the Y is high. Use of a 1" cell should be considered.

The choice of cell path-length will impact on accuracy. Unless working to a particular standard or specification, the optical path-length of the cell used should be related to the colour intensity of the sample.

When comparing with other measurements, it is necessary to check that the cell path-length and type (Optical Glass, Borosilicate or Plastic) are identical and the cells used are clean and undamaged.

- Confirm the correct colour scale is selected on the Automatic instrument. With CIE $L^*a^*b^*$, other settings such as light source and observer (10° or 2°) also need to be defined and communicated.

Historically, a number of scales are available that report Red & Yellow values. This is a common source of error. For example; a standard Model F reports Lovibond® Red, Yellow, Blue and Neutral units (RYBN). An AF710 reports AOCS Tintometer® in terms of Red & Yellow (RY). An Automatic instrument may be configured to display both RYBN and RY.

Accuracy of Measurements

For peace of mind in production and quality control, it is crucial that the correct performance of instruments can easily be confirmed.

There are different options available:

1) Liquid Reference Standards:

High quality liquid samples with known values are used to check that an instrument is reporting the correct measurements. The range of Lovibond® colour reference standards includes AOCS-Tintometer®, ASTM, Gardner, Lovibond® RYBN, Pt-Co/Hazen/ APHA and Saybolt Colour. Each standard is shipped with a 12-month guarantee of colour stability.

2) Glass Filters: Conformance filter-sets allow quick and simple conformance checks on Lovibond® instruments. Each filter set is supplied with a Certificate of Conformity that confirms that they have been manufactured under the control of the ISO 9001: 2000 Quality Management System.



Liquid Reference Samples

3) Solid Reference Standards:

High quality physical standards with known values are used to check that a Lovibond® instrument is reporting the correct results.

Web based solutions are becoming increasingly popular. They allow remote testing and calibration of automatic instruments using conformance standards.

Bespoke reference liquid standards and conformance filters can be requested and the Tintometer® Group always endeavours to match the requested value as closely as possible.

On occasion, it may not be possible to match the requested value exactly. The value achieved and the expected performance tolerances will, of course, always be reported.



Glass Conformance Filters

Correct Record Keeping

Colour Consistency is simple to achieve if care is taken to eliminate possible causes of variation between instruments. Measurement procedures should be documented, circulated and complied with.

Ideally, the following information should be recorded for each colour measurement:

- Instrument model used
- Serial number (if multiple instruments are available)
- Date
- Time (if multiple samples are measured on the same day)
- Name of observer/operator
- Description of sample
- Sample reference number (if relevant)
- Temperature of sample if heated

- Any comments on condition of sample (e.g. turbid or dirty etc)
- Path-length of optical cell
- Colour space or colour scale used
- Light source and observer (if using CIE systems)
- Individual colour values
- Any other information required by the organisation
- Any additional comments relevant to the colour match

Special care should be taken with regards to colour scales utilising Red and Yellow values to avoid miscommunication.

Instruments should be maintained correctly and regularly checked using third party standards.

The Supply Chain

Defining the colour of a product and ensuring colour accuracy every time is critical to long-term success and accurate communication within a supply chain.

The clear and correct description of colour standards and tolerances is critical when:

- Specifying materials when sourcing
- Communicating colour within the wider supply chain
- Inspecting incoming materials
- Conducting continual production Quality Control
- Inspecting final/outgoing products
- Guaranteeing compliance with National and International Standards

Reliable and repeatable colour test results are the key to ensuring final product quality, and also to minimising production costs. Speed of analysis can also be vital for efficient process control. Simplicity of operation helps to reduce error and increase productivity.

Methods of Quantifying Colour – Colour Scales and Spaces

There are many different types of colour systems available. Some are applicable to any type of sample, whereas some are specific to opaque materials or transparent materials. The use of specific colour scales or parameters varies from one industry to another depending on standards and requirements.

Lovibond® RYBN Colour

The Lovibond® RYBN colour scale is optimised for the colour measurement of clear (but coloured) liquids. In the 1890's, Joseph Lovibond, the founder of The Tintometer Ltd, developed the original Lovibond® Scale, based on a calibrated series of red, yellow, blue and neutral glasses.

The Lovibond® Scale is based on 84 calibrated glass colour standards of different densities of magenta (red), yellow, blue and

neutral glasses, that graduate from de-saturated to fully saturated. Sample colours are matched by a suitable combination of the three primary colours together with neutral filters, resulting in a set of Lovibond® RYBN units that define the colour.

Since several million combinations are available, it is possible to match the colour of almost any sample. The scale is particularly popular for measuring the colour of oils and fats, chemicals, pharmaceuticals and syrups. After more than a century, the Tintometer® Group still manufactures and grades the glass filters used for visual colour measurement in terms of Lovibond® units. It is this unparalleled knowledge and experience that has enabled the company to accurately replicate the scale in its automatic instruments.



Model F with Racks

The scale quoted by others as the Lovibond® scale does not guarantee validated Lovibond® Colour readings and may not conform to any visual instrument for Lovibond® Colour.

Neutral Filters

If, for any reason, an operator alters the method of use or changes any convention, it is important that they should give details when recording results, otherwise confusion could ensue. For example, observers employ neutral filters to dull a bright sample but omit to report the fact. In other cases they endeavour to make the best possible match without stating neutral values although they were needed, or use different colours in combination only in a fixed ratio according to some arbitrary convention.

Colour Nomenclature

The Lovibond® Scale provides its own simple language of colour which can fully describe the appearance of any colour in the least possible number of words and figures to avoid language difficulties. For convenience of laboratory records, or in communicating readings between laboratories, many industries record their results on a three colour basis, quoting the Red, Yellow and Blue instrumental values.

Some industries find it more convenient to simplify these terms by using the six divisions of the spectrum.

Red

Orange – combination of red and yellow.

Yellow

Green – combination of yellow and blue.

Blue

Violet – combination of red and blue.

These six terms are used in combination with “bright” and “dull”.

A sample is described as being bright when the nearest possible match appears dull in comparison. When this occurs, neutral values are introduced and recorded as sample brightness.

A sample is described as being dull when red, yellow and blue are required to make a match. When this occurs, the colour with the lowest value is expressed as dullness.

The Munsell Scale

In 1905, artist Albert H. Munsell originated a colour order system – or colour scale – which is still used today. The Munsell System of Colour Notation is significant from a historical perspective because it is based on human perception. Moreover, it was devised before instrumentation was available for measuring and specifying colour. The Munsell System assigns numerical values to three properties of colour known as Munsell Hue, Munsell Value and Munsell Chroma. Adjacent colour samples in the system represent equal intervals of visual perception.

The model in Figure 7 depicts the Munsell Colour Tree which is a representation of the Munsell Book of Color which provides physical samples for judging visual colour. Today's colour systems rely on instruments that utilise mathematics to help us judge colour.

As noted earlier, there are three things necessary to see colour:

- A source of light (also known as the illuminant)
- An object (also known as the sample)
- An observer/processor

CIE Colour Systems

The CIE, or Commission Internationale de l'Eclairage (translated as the International Commission on Illumination), is the body responsible for international recommendations for photometry and colorimetry. In 1931, the CIE first recommended a colour measurement system by specifying the light source (or illuminant), the observer and the methodology used to derive values that describe colour. As time has passed the system has been updated and added to.

The CIE Colour Systems utilise three coordinates to locate a colour in a 'colour space'. These colour spaces include-

- CIE XYZ
- CIE L*a*b*
- CIE L*C*h.

Colours can be quantified using these colour spaces by different calculations based upon specification of light source type and defining 'standardised' observers (these two parameters are accounted for using different numerical data).

Instruments quantify colour by gathering and filtering the wavelengths of light transmitted through, or reflected from, an object. The instrument detects the different intensities of different light wavelengths and these intensity values are recorded as points across the visible spectrum. This Spectral data is represented as a spectral reflectance or transmittance curve. This curve is the colour's fingerprint (Figure 8).




Figure 7: The Munsell Colour Tree




Figure 8: Daylight (standard illuminant D65)

Once we obtain the transmittance or reflectance curve of a colour, we can apply mathematics to map the colour into a colour space.

To do this, we take the reflectance or transmittance curve and multiply the data by a CIE standard illuminant or other illuminant. The illuminant is a theoretical representation of the light source under which the samples are viewed. Each light source has a power distribution that affects how we see the colour. Examples of different illuminants are Illuminant A – incandescent, Illuminant D65 – daylight (Figure 8) and Illuminant F2 – cool white fluorescent.

We multiply the result of this calculation by a CIE standard observer. The CIE commissioned work in 1931 and 1964 to derive the concept of a standard observer, which is based on the average human response to wavelengths of light (Figure 9). In short, a standard observer represents how an average person sees colour across the visible spectrum when using a defined area of the

retina in the eye. These values can now be used to identify a colour numerically.

CIE XYZ

The three co-ordinates required to define a colour in the CIE XYZ colour space are known as the Tristimulus values (XYZ). These are calculated using the standard illuminant, the sample's spectral curve and a standard observer. Unfortunately Tristimulus values have limited use as colour specifications because they correlate poorly with visual attributes. While Y relates to lightness, X and Z do not correlate to any visual attributes.

As a result, when the 1931 CIE standard observer was established, the CIE defined the x,y chromaticity coordinates which can be correlated to chroma and hue. These are derived from XYZ. The x,y coordinates are used to form the chromaticity diagram in Figure 11. The notation Y,x,y specifies colours by identifying lightness (Y) and the colour as viewed in the chromaticity diagram (x,y).




Figure 9: CIE 2° and 10° Standard Observers




Figure 10: Tristimulus Values

As Figure 12 shows, hue is represented at all points around the perimeter of the chromaticity diagram. Chroma, or saturation, is represented by a movement from the central white (neutral) area out toward the perimeter of the diagram, where 100% saturation defines the pure hue.




Figure 12: Chromaticity diagram




Figure 11: CIE 1931 (x,y) chromaticity diagram

Expressing Colour Uniformly

One issue with the x,y chromaticity diagram is that the different colours are not uniformly distributed. In an attempt to solve this problem, the CIE recommended a more uniform chromaticity diagram u',v' . In addition, CIE recommended two alternate, uniform colour scales: CIE 1976 ($L^*a^*b^*$) or CIELAB, and CIE 1976 ($L^*u^*v^*$) or CIELUV.

These colour scales are based on the opponent-colours theory of colour vision, which says that a colour cannot appear to be a mixture of both green and red at the same time, or both blue and yellow at the same time.

As a result, single values can be used to describe the red/green and the yellow/blue attributes.

CIELAB $L^*a^*b^*$

When a colour is expressed in CIELAB, L^* defines lightness, a^* denotes the red/green value and b^* the yellow/blue value.

Figures 13 and 14 show the colour-plotting diagrams for $L^*a^*b^*$. The a^* axis runs from left to right. A colour measurement movement in the $+a^*$ direction depicts a shift toward red. Along the b^* axis, a $+b^*$ movement represents a shift toward yellow. The centre L^* axis shows $L^* = 0$ (black or total absorption) at the bottom and $L^* = 100$ or white at the top. In between are greys. All colours on this axis can be considered as neutrals as they are not coloured in any particular direction. To demonstrate how the $L^*a^*b^*$ values represent the specific colours of Flowers A and B, they are plotted on the CIELAB Colour Chart in Figure 13.





Figure 13: CIELAB colour chart



The L^* value is represented on the centre axis. The a^* and b^* axes appear on the horizontal plane


Figure 14

These points specify each flower's hue (colour) and chroma (vividness/dullness).


When their L^* values (lightness) are added in Figure 14, the final colour of each flower is obtained.

CIELAB L^*C^*h

While CIELAB $L^*a^*b^*$ uses Cartesian coordinates to represent a colour in a colour space, CIELAB L^*C^*h uses polar coordinates. The C^* and h values can be derived from CIELAB a^* and b^* . The L^* , as before, defines lightness, C^* specifies chroma and h denotes hue angle, an angular measurement of hue.



Flower A



Flower B

Metamerism

Coloured materials can sometimes exhibit metamerism, a phenomenon whereby a pair of colours that match under one illuminant does not match under a second different illuminant. To limit metameric effects, it is always advisable to match colours under a Primary, Secondary and Tertiary illuminant.

Colour Differences, Notation and Tolerancing

Delta CIELAB L*a*b* and CIELAB L*C*h

Assessment of colour is often more than a numeric expression. Usually it is an assessment of the colour difference (delta) of a sample relative to a known standard or reference. Difference in CIELAB (L*a*b*) and CIELAB (L*C*h) values can both be used to compare the colours of two objects.

The expressions for these colour differences are ΔL^* , Δa^* , Δb^* , or ΔL^* ΔC^* ΔH^* (Δ symbolises "delta," which indicates difference). Given ΔL^* , Δa^* and Δb^* , or ΔL^* , ΔC^* and ΔH^* , the total difference or distance in CIELAB space can be stated as a single value, known as ΔE^* .

Delta (Δ) E* colour difference

$$\Delta E^*_{ab} = [(\Delta L^*)^2 + (\Delta a^*)^2 + (\Delta b^*)^2]^{1/2}$$

$$\Delta E^*_{ab} = [(\Delta L^*)^2 + (\Delta C^*)^2 + (\Delta H^*)^2]^{1/2}$$

Let us compare the colour of Flower A to Flower C. Separately, each would be classified as a yellow rose. But what is their relationship when set side by side? How do the colours differ?

Using the equation for ΔL^* , Δa^* , Δb^* , the colour difference between Flower A and Flower C can be expressed as:

$$\Delta L^* = +11.10$$

$$\Delta a^* = -6.10$$

$$\Delta b^* = -5.25$$

The total colour difference can therefore be expressed as $\Delta E^* = 13.71$

The values for Flowers A and C are shown at the bottom of this page. On the a* axis, a reading of $\Delta a^* = -6.10$ indicates greener or less red. On the b* axis, a reading of $\Delta b^* = +5.25$ indicates bluer or less yellow. On the L* plane, the measurement difference of $\Delta L^* = +11.10$ shows that Flower C is lighter than Flower A.


If the same two flowers were compared using CIELAB L*C*h, the colour differences would be expressed as:

$$\Delta L^* = +11.10$$


$$\Delta C^* = -5.88$$

$$\Delta H^* = +5.49$$

Referring again to the flowers shown below, the ΔC^* value of -5.88 indicates that Flower C is less chromatic, or less saturated. The ΔH^* value of +05.49 indicates that Flower C is greener in hue than Flower A. The ΔL^* values are identical for CIELAB L*C*h and CIELAB L*a*b*.



Flower A



Flower C

CIE Colour Space Notation

ΔL^* = difference in lightness/darkness value
(+ve = lighter, -ve = darker)

Δa^* = difference on red/green axis
(+ve = redder, -ve = greener)

Δb^* = difference on yellow/blue axis
(+ve = yellower, -ve = bluer)

ΔC^* = difference in chroma
(+ve = brighter, -ve = duller)

ΔH^* = difference in hue

ΔE^* = total colour difference value

Refer to Figure 13 for visualisation.

Visual Colour and Tolerancing

Tolerances for an acceptable colour match typically consist of a three-dimensional boundary with varying limits for lightness, hue and chroma, and must agree with visual assessment. CIELAB can be used to create those boundaries. The simplest method to create a spherical tolerance is Delta E* tolerancing. Additional tolerancing formulas, known as CMC and CIE94, produce ellipsoidal tolerances.

Delta E* Tolerancing

Delta E* is the total colour difference computed with a colour difference equation as defined above:

$$\Delta E^*_{ab} = [(\Delta L)^2 + (\Delta a)^2 + (\Delta b)^2]^{1/2}$$

$$\Delta E^*_{ab} = [(\Delta L^*)^2 + (\Delta C^*)^2 + (\Delta H^*)^2]^{1/2}$$

It is used widely as a single value method of identifying if a colour is within or outside a specified tolerance – ie an acceptable Pass or Fail. Delta E* should not be confused with ΔE CMC or other single value tolerancing techniques.

CIELAB L*a*b* Tolerancing

When tolerancing with L*a*b*, you may choose a difference limit for ΔL^* (lightness), Δa^* (red/green), and Δb^* (yellow/blue). These limits create a tolerance cuboid box around the standard (Figure 15).




Figure 15: CIELAB tolerance box

When comparing this tolerance box with an ellipsoid/sphere tolerance, some problems emerge. A box-shaped tolerance around the ellipsoid/sphere can give acceptable values for unacceptable colour differences (Figure 16). If the tolerance box is made small enough to fit within the ellipsoid/sphere, it is possible to get unacceptable values for visually acceptable colour differences.




Figure 16: Numerically correct vs. Visually acceptable




Figure 17: CIELAB (L*C*h) tolerance wedge

CIELAB L*C*h Tolerancing

CIELAB L*C*h users must choose a difference limit for ΔL^* (lightness), ΔC^* (chroma) and ΔH^* (hue). This creates a wedge-shaped box around the standard. Since CIELAB L*C*h is a polar-coordinate system, the tolerance box can be rotated in orientation with the hue angle (Figure 17).

When this tolerance is compared with the ellipsoid, we can see that it more closely matches human perception. This reduces the amount of disagreement between the observer and the instrumental values (Figure 18).




Figure 18: CIELAB (L*C*h) tolerance ellipsoids

CMC Tolerancing

As the eye does not detect differences in hue (red, yellow, green, blue, etc.), chroma (saturation) or lightness equally, the average observer will have variable sensitivity to hue, chroma and lightness differences: the eye

is usually more sensitive to changes in hue, then chroma and lastly lightness. Because of this, visual acceptability is best represented by an ellipsoid (Figure 19).




Figure 19: Tolerance Ellipsoid

CMC is not a colour space but rather a tolerancing system. CMC tolerancing is based on CIELAB L^*C^*h and provides better agreement between visual assessment and measured colour difference. CMC tolerancing was developed by the Colour Measurement Committee of the Society of Dyers and Colourists in Great Britain and is now recognised by a British Standard (BS 6923).

The CMC calculation mathematically defines an ellipsoid around the standard colour with semi-axis corresponding to hue, chroma and lightness. The ellipsoid represents the volume of acceptable colour and automatically varies in size and shape depending on the position of the colour in colour space.

Figure 20 shows the variation of the ellipsoids throughout colour space. The ellipsoids in the orange region of colour space are longer and narrower than the broader and rounder ones in the green region. The size and shape of the ellipsoids also change as the colour varies in chroma and/or lightness.

The CMC equation allows you to vary the overall size of the ellipsoid to better match what is visually acceptable. By varying the commercial factor (cf), the ellipsoid can be made as large or small as necessary to match visual assessment. The cf value is the tolerance, which means that if $cf=1.0$, then a value of ΔE CMC less than 1.0 would pass, but more than 1.0 would fail (see Figure 21, *top of the next page*).

Since the eye will generally accept larger differences in lightness (l) than in chroma (c), a default ratio for (l:c) is 2:1 which will allow twice as much difference in lightness as in chroma. The CMC equation allows this ratio to be adjusted to achieve better agreement with visual assessment in specific industries.




Figure 20: Tolerance ellipsoids in colour space

- Tolerance ellipsoids are tightly packed in the orange region
- Tolerance ellipsoids are larger in the green area




Figure 21: Commercial Factor (CF) of tolerances




Figure 22: CMC Tolerance Ellipsoids

CIE2000 Tolerancing

This is similar to CMC Tolerancing but based on a new colour-difference equation recommended by CIE in 2000.

Choosing the Right Tolerance


When deciding which colour difference calculation to use, consider these five rules;

1. Select a single method of calculation and use it consistently.
2. Always specify exactly how the calculations are made.
3. Never attempt to convert between colour differences calculated by different equations through the use of average factors.
4. Use calculated colour differences only as a first approximation in setting tolerances, until they can be confirmed by visual judgments.
5. Always remember that nobody accepts or rejects colour because of numbers – it is the way it looks that counts.

Colour Scales

Although the CIE colour notation systems can be used for any colour control application, it is often simpler to use an industry specific colour scale for routine grading of many types of material.

Grading techniques are widely used to assess product colour by comparison with a representative series of fixed colour standards. For many product types, a characteristic set of standards has been agreed and adopted to aid colour control and the communication of colour specifications; the result is a selection of traditional colour grading scales that have been adopted as industry standards and are still in common use today.



Example Colour Scale

The following Scales are examples of those currently most widely used by industry. Please note: when measuring liquids, cell path-length is critical as the same material measured in different path-length cells will appear a different colour.

The specific path-length of cell must be stated for accurate colour communication. The measurement standards issued for some colour scales require specific path-length cells. More details and the latest information can be found within the Scales Section at www.lovibondcolour.com

ASTM Colour (ASTM D1500, ASTM D6045, ISO 2049, IP 196)

The ASTM Colour Scale is widely utilised for the grading of petroleum products such as lubricating oils, heating oils and diesel fuel oils.

Mineral oils are constantly checked for colour during processing in order to establish when they have been refined to the required grade. Colour is also used as a means of confirming that the correct oil or fuel is being used for its intended purpose and that no contamination or degradation of quality has occurred.

ASTM D1500 is a single number, one dimensional, colour scale ranging from a pale straw colour through to a deep red colour in sixteen steps (0.5 - 8.0 units, in increments of 0.5 units). Visual Lovibond® Comparators can achieve a resolution of 0.5 units; Automatic instruments can achieve a resolution of 0.1 units. To comply with specifications, a 33mm path-length cell must be used. ASTM D1500 superseded the 12-step D155 NPA (National Petroleum Association) scale in 1960.

Other petroleum products that do not fall within the scope of ASTM D1500, such as undyed gasoline, white spirit, petroleum wax and kerosene, may be graded using the Saybolt test ASTM D156 or IP (Institute of Petroleum) 17.

EBC (European Brewing Convention)

The EBC colour scale, developed by the Institute of Brewing and the European Brewing Convention, is a recognised method for the colour grading of beers, malts and caramel solutions as well as similarly coloured liquids. It has a range of 2 to 27 visual units; yellower pale worts and lagers at the low end of the scale and the amber of dark worts, beers and caramels at the upper end of the scale.

If the sample falls outside this range (e.g. concentrates, syrups) then sample dilution

and/or a different path-length cell can be used to bring the reading within the EBC range.

European Pharmacopoeia (EP) Colour

The EP Colour Standards were originally visual colour standards intended to improve colour communication between sites by defining a sample colour as being close to a physical liquid standard (“near EP Y2”) rather than using the words “light yellow”.

EP consists of 3 primary standard colour solutions (yellow, red, blue) that are combined with hydrochloric acid to make 5 standard solutions that, when further diluted with hydrochloric acid (10 mg/l), make 37 reference EP standards; Red (R1 - R7); Yellow (Y1 - Y7); Brown (B1 - B9); Brown/Yellow (BY1 - BY7); Green/Yellow (GY1 - GY7).

Automatic instruments can now measure EP Colour. Measurements are based on specific wavelengths for each of the European Pharmacopoeia colour scales (US and Chinese Pharmacopoeia scales are local alternatives).

Gardner Colour (ASTM D1544, ASTM D6166, AOCS Td 1a, MS 817 Part 10)

The Gardner Colour scale as specified in ASTM D1544 is a single number, one dimensional, colour scale for grading the colour of similarly coloured liquids such as resins, varnishes, lacquers, drying oils, fatty acids, lecithins, sunflower oil and linseed oil.

The scale ranges from pale yellow to red in shade and is described in terms of the values 1-18. The glass standards used with the Lovibond® Comparator can achieve a

resolution of 1 unit. Automatic instruments can achieve a resolution of 0.1 units. To comply with the specifications, a 10mm path-length must be used.

The light yellow Gardner colour numbers (1 to 8) are based on potassium chloroplatinate solutions, numbers 9 to 18 on solutions of ferric chloride, cobaltous chloride and hydrochloric acid.

In 1958, The Tintometer Ltd was instrumental in the development of the master glass standards that were utilised when the current Gardner scale was specified in 1963; the earlier 1933 and 1953 versions are available upon request in the form of Lovibond® Comparator discs.

Platinum-Cobalt/Hazen/APHA Colour (ASTM D1209)

This scale can be referred to as Pt-Co, Platinum-Cobalt, Hazen or APHA Colour. All terms are interchangeable and equally valid. It is used to measure clear to dark amber liquids.

The scale was originally defined by specified dilutions of a platinum-cobalt stock solution, ranging from 0 at the light end of the scale to 500 at the dark end. The scale is now available in a digital format on Automatic instruments. The scale is used extensively in the water industry but also for clear oils, chemicals and petrochemicals such as glycerine, plasticisers, solvents, carbon tetrachloride and petroleum spirits.

Saybolt Colour (ASTM D156, ASTM D6045)

The Saybolt Colour scale is used for grading light coloured petroleum products including

aviation fuels, kerosene, naphthas, white mineral oils, hydrocarbon solvents and petroleum waxes.

The colour range of the Saybolt-scale is similar to that of the Platinum-Cobalt/Hazen/APHA Colour (ASTM D1209) scale and is therefore employed for the measurement of clear water and colourless to slightly yellowish products. The faintest coloration is Saybolt-colour number +30; the strongest evaluable Saybolt coloration value is -16.

White and Yellow Indices

Certain industries, such as paint, textiles and paper manufacturing, evaluate their materials and products based on standards of whiteness. Typically, this whiteness index is a preference rating for how white a material should appear, be it photographic and printing paper or plastics.

Therefore the Whiteness Index is a measurement which correlates the visual rating of whiteness for certain white and near-white samples.

The Yellowness Index is a number calculated from spectrophotometric data that describes the change in colour of a test sample from colourless through to yellow.

ASTM has defined whiteness and yellowness indices. The E313 whiteness index is used for measuring near-white, opaque materials such as paper, paint and plastic.

Tintometer GmbH

Lovibond® Colour Measurement
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
Fax: +49 (0)231/94510-20
sales@tintometer.de

Germany

The Tintometer Limited

Lovibond House
Sun Rise Way, Amesbury
SP4 7GR
Tel.: +44 (0)1980 664800
Fax: +44 (0)1980 625412
sales@tintometer.com

UK

Tintometer AG

Hauptstraße 2
5212 Hausen AG
Tel.: +41 (0)56/4422829
Fax: +41 (0)56/4424121
info@tintometer.ch

Switzerland

Tintometer Inc

6456 Parkland Drive
Sarasota
FL 34243
Tel.: +1 941 756 6410
Fax: +1 941 727 9654
sales@tintometer.us

USA

Tintometer China

Room 1001, China Life Tower,
16 Chaoyangmenwai Avenue,
Beijing, 100020
Tel: +86 10 85251111 ext.330
Fax: +86 10 85251001
chinaoffice@tintometer.com

China

Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,
Lebuh Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com

Malaysia

Tintometer Brazil

Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@tintometer.com.br

Brazil

Tintometer India Pvt. Ltd.

B-91, A.P.I.E. Sanath Nagar,
Hyderabad, 500018
Tel.: +91 (0)40 4647 9911
Toll Free: 1 800 102 3891
indiaoffice@tintometer.com

India

www.lovibondcolour.com

Lovibond® & Tintometer® are Registered Trade Marks of the Tintometer® Group. All translations and transliterations of LOVIBOND® & TINTOMETER® are asserted as Trade Marks of the Tintometer® Group. Registered Office: Lovibond House, UK.

Registered in England No. 45024 –
Errors and Omissions Excepted –
Content subject to alterations without notice.
920173_V3_02/16