


# Lovibond® Water Testing

Tintometer® Group



## Manual of Methods

MD50 • MD150





KS4.3 T / 20


Method name

Method number

Bar code for the detection of the methods

$K_{S4.3 T}$ 
20

0.1 - 4 mmol/l  $K_{S4.3}$ 
S:4.3

Acid / Indicator

Measuring range

Chemical Method

Display in the MD 100 / MD 110 / MD 200

### Instrument specific information

The test can be performed on the following devices. In addition, the required cuvette and the absorption range of the photometer are indicated.

Instrument Type	Cuvette	$\lambda$	Measuring Range
MD 200, MD 600, MD 610, MD 640, MultiDirect, PM 620, PM 630	ø 24 mm	610 nm	0.1 - 4 mmol/l $K_{S4.3}$
SpectroDirect, XD 7000, XD 7500	ø 24 mm	615 nm	0.1 - 4 mmol/l $K_{S4.3}$

### Material

Required material (partly optional):

Reagents	Packaging Unit	Part Number
Alka-M-Photometer	Tablet / 100	513210BT
Alka-M-Photometer	Tablet / 250	513211BT

### Application List

- Waste Water Treatment
- Drinking Water Treatment
- Raw Water Treatment

### Notes

1. The terms Alkalinity-m, m-Value, total alkalinity and Acid demand to  $K_{S4.3}$  are identical.
2. For accurate results, exactly 10 ml of water sample must be used for the test.

Language codes ISO 639-1

Revision status

EN Handbook of Methods 01/20

Performing test procedure

### Implementation of the provision Acid capacity $K_{S_{4.3}}$ with Tablet

Select the method on the device

For this method, no ZERO measurements are to be carried out with the following devices: XD 7000, XD 7500



Fill 24 mm vial with **10 ml sample**.



Close vial(s).



Place **sample vial** in the sample chamber. • Pay attention to the positioning.

• • •



Dissolve tablet(s) by inverting.




Place **sample vial** in the sample chamber. • Pay attention to the positioning.



Press the **TEST** (XD: **START**) button.

The result in Acid Capacity  $K_{S_{4.3}}$  appears on the display.



**Suspended solids 24**

**M384**

**10 - 750 mg/L TSS**

**SuS**

**Turbidity / Attenuated Radiation Method**

EN

## Material

Required material (partly optional):

Reagents	Packaging Unit	Part Number
no reagent required		

## Sampling

1. Measure the water sample as soon as possible after sampling. It is possible to store the sample at 4 °C for 7 days in plastic or glass containers. The measurement should be at the same temperature as the sample. Temperature differences between measurement and sampling can change the result of the measurement.

## Notes

1. The photometric determination of Suspended Solids is based on a gravimetric method. In a laboratory this is usually done by evaporation of the filter residue of a filtrated water sample in a furnace at 103 °C – 105 °C and weighing of the dried residue.
2. When higher accuracy is required perform a gravimetric determination of a water sample. The result can be used to calibrate the photometer with the same water sample.
3. The estimated detection limit is 20 mg/L TSS.

## Determination of Total suspended solids


Select the method on the device.

For this method, a ZERO measurement does not have to be carried out every time on the following devices: XD 7000, XD 7500


Homogenize mL of the water sample in a blender on high speed for minutes




Fill 24 mm vial with **10 mL deionised water** .




Close vial(s).



Place **sample vial** in the sample chamber. Pay attention to the positioning.




Press the **ZERO** button.




Remove the vial from the sample chamber.


For devices that require **no ZERO measurement** , start here.






Mix homogenised water sample thoroughly.



Pre-rinse vial with water sample.



Fill 24 mm vial with **10 mL prepared sample** .



EN

Place **sample vial** in the sample chamber. Pay attention to the positioning.

Press the **TEST** (XD: **START**) button.

The result in mg/L TSS (Total Suspended Solids) appears on the display.

## Chemical Method

Turbidity / Attenuated Radiation Method

## Appendix

### Interferences

#### Persistent Interferences

- Colour interferes if light is absorbed at 660 nm.

#### Removeable Interferences

- Air bubbles interfere and can be removed by swirling the vial gently.

### Method Validation

<b>Limit of Detection</b>	10 mg/L
<b>Limit of Quantification</b>	30 mg/L
<b>End of Measuring Range</b>	750 mg/L
<b>Sensitivity</b>	550 mg/L / Abs
<b>Confidence Intervall</b>	4.24 mg/L
<b>Standard Deviation</b>	1.79 mg/L
<b>Variation Coefficient</b>	0.47 %

#### Derived from

EN 872:2005











**Tintometer GmbH**

Lovibond® Water Testing  
Schleefstraße 8-12  
44287 Dortmund  
Tel.: +49 (0)231/94510-0  
sales@lovibond.com  
www.lovibond.com  
Germany

**Tintometer South East Asia**

Unit B-3-12, BBT One Boulevard,  
Lebuh Nilam 2, Bandar Bukit Tinggi,  
Klang, 41200, Selangor D.E  
Tel.: +60 (0)3 3325 2285/6  
Fax: +60 (0)3 3325 2287  
lovibond.asia@tintometer.com  
www.lovibond.com  
Malaysia

**Tintometer India Pvt. Ltd.**

Door No: 7-2-C-14, 2<sup>nd</sup>, 3<sup>rd</sup> & 4<sup>th</sup> Floor  
Sanathnagar Industrial Estate,  
Hyderabad, 500018  
Telangana  
Tel: +91 (0) 40 23883300  
Toll Free: 1 800 599 3891/ 3892  
indiaoffice@lovibond.in  
www.lovibondwater.in  
India

**The Tintometer Limited**

Lovibond House  
Sun Rise Way  
Amesbury, SP4 7GR  
Tel.: +44 (0)1980 664800  
Fax: +44 (0)1980 625412  
sales@lovibond.uk  
www.lovibond.com  
UK

**Tintometer Brazil**

Caixa Postal: 271  
CEP: 13201-970  
Jundiaí – SP  
Tel.: +55 (11) 3230-6410  
sales@lovibond.us  
www.lovibond.com.br  
Brazil

**Tintometer Spain**

Postbox: 24047  
08080 Barcelona  
Tel.: +34 661 606 770  
sales@tintometer.es  
www.lovibond.com  
Spain

**Tintometer China**

9F, SOHO II C.  
No.9 Guanghualu,  
Chaoyang District,  
Beijing, 100020  
Customer Care China Tel.: 4009021628  
Tel.: +86 10 85251111 Ext. 330  
Fax: +86 10 85251001  
chinaoffice@tintometer.com  
www.lovibond.com  
China

**Tintometer Inc.**

6456 Parkland Drive  
Sarasota, FL 34243  
Tel: 941.756.6410  
Fax: 941.727.9654  
sales@lovibond.us  
www.lovibond.us  
USA



Technical changes without notice  
Printed in Germany 08/24

No.: 00386766

Lovibond® and Tintometer® are Trademarks of  
the Tintometer Group of Companies

