# Lovibond® Water Testing

Tintometer® Group

- Ⓜ GB Alkalinity P
- Ⓜ DE Alkalität P-Wert
- Ⓜ FR Alcalinité P
- Ⓜ IT Valore Alcalinità P
- Ⓜ ES Valor P de alcalinidad
- Ⓜ PT Alcalinidade/ Valor P
- Ⓜ DK Alcalitet P-værdi
- Ⓜ NL Alkaliteit P-waarde

## MINIKIT

20 - 500 mg/l  
CaCO<sub>3</sub>  
(1,12 - 28°dH)

41 41 40

## **GB** Introduction

The Alkalinity P Minikit is for determination of alkalinity where the pH of the water is over 8.5. The test range is 20 - 500 mg/l  $\text{CaCO}_3$  and based on the tablet count method. Alkalinity P is a measure of the water up to the pH-value 8.3 (colour change of phenolphthaleine) resultant from the presence of hydroxides, carbonates, bicarbonates and a minor part from silicates and phosphates.

Material Safety Data Sheets: [www.lovibond.com](http://www.lovibond.com)

### **Contents**

calibrated sample container, 100 ml  
calibrated sample container, 30 ml  
Alkalinity P tablets, 100 pieces  
operating manual in 8 languages

### **Instructions**

Range 20 - 500 mg/l  $\text{CaCO}_3$                       Sample Size 100 ml

1. Take the sample of water into the sample container provided, filling to the 100 ml mark.
2. Add one Alkalinity P tablet and shake the container until the tablet disintegrates.

Continue adding Alkalinity P tablets one at a time in this manner until the colour of the solution changes from

**BLUE to YELLOW.**

3. Note the number of tablets used until the colour of the solution changed and calculate the result from the formula below appropriate to the sample volume taken:

$$\text{Alkalinity P (mg/l CaCO}_3\text{)} = (\text{No. of tablets} \times 40) - 20$$

### **Example**

No. of tablets used = 6

$$(6 \times 40) - 20 = 220 \text{ mg/l CaCO}_3$$

### **Cleaning**

Thoroughly rinse out the sample container after each use.

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are Trademarks of the
The Tintometer Ltd.
Tintometer GmbH