DE

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GB

Test Instructions

1. Use the photometer as described in the instructions in the “Abs”/“Measure absorbance“-Mode.
2. Clean the outside of the vial with a clean tissue.
3. Zero the instrument according to instrument instructions using the “Verification Zero”.
4. Use the Verification Standard for the required wavelength instead of the test sample and take the test reading.
The reading should correspond to the stated calibration value within the given tolerances.

Notes

1. If an instrument and the Verification Standards show Non-Compliance:
 - the cell chamber should be checked for dirtiness, if so, cleaned and the measurement should be repeated
 - the instrument and the Verification Standards should be returned to Tintometer GmbH, Technical Service Department, or an authorized distributor, for servicing
2. Place standards in the case after use.
Store in the dark at room temperature (15 - 25°C).
Do not use the standards after the stated expiry date.
3. Verification Standards should be used under laboratory conditions.
Avoid working in direct sunlight and carry out calibrations at a temperature of $20 \pm 3^\circ\text{C}$.