# Lovibond® Water Testing

Tintometer® Group

Ⓜ GB Nitrite

Ⓜ DE Nitrit

Ⓜ FR Nitrite

Ⓜ IT Nitrito

Ⓜ ES Nitrito

Ⓜ PT Nitrito

Ⓜ DK Nitrit

Ⓜ NL Nitriet

## MINIKIT

70 - 1540 mg/l  
NaNO<sub>2</sub>

41 42 70

## **GB** Introduction

The Nitrite Minikit provides a simple method of checking nitrite levels in water over the range 70 - 1540 mg/l  $\text{NaNO}_2$  and it is based on the tablet count method.  
Material Safety Data Sheets: [www.lovibond.com](http://www.lovibond.com)

### **Contents**

Calibrated sample container, 100 ml  
Nitrite No.1 tablets, 40 pieces in foil  
Nitrite No.2 tablets, 100 pieces in foil  
Instruction manual in 8 languages

### **Instructions\***

Range 140 - 1540 mg/l (70-770 mg/l)                      Sample Size 10 ml (20 ml)

1. Take a 10 ml (20 ml) sample size of the water under test in the sample container. Make up to the 50 ml mark using distilled or nitrite free water.
2. Add two (four) Nitrite No.1 tablets, stopper the container and shake until the tablets have disintegrated.
3. Add one Nitrite No.2 tablet, stopper the container and shake until the tablet disintegrates. Continue adding Nitrite No.2 tablets in this manner one at a time until a PINK colour persists for one minute.

4. Note the number of Nitrite No.2 tablets used and calculate the result from the formula below: Nitrite (as  $\text{NaNO}_2$ ) = No. of Nitrite No.2 tablets x 140 (70)

### **Example**

Nitrite No.2 tablets used = 5

$5 \times 140 (70) = 700 \text{ mg/l (350 mg/l) NaNO}_2$

### **Cleaning**

Thoroughly rinse out the sample container and the stopper after each use.

### **Notes**

1. If the cooling water sample is cloudy, filter as necessary to obtain a clear solution prior to carrying out the test.
2. Note that the Nitrite No.1 tablets are to condition the sample only.  
Do not include these tablets in the calculation of nitrite level.

\*Note: 10 ml or 20 ml sample volume can be used. For 20 ml sample volume the values in brackets apply.

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are Trademarks of the
The Tintometer Ltd.
Tintometer GmbH