Lovibond® Water Testing

Tintometer® Group

Photometer-System

MD 6x0 • PM 6x0

(EN) Instruction manual..... 3	(DE) Gebrauchsanleitung..... 36
(ES) Manual de instrucciones..... 74	(FR) Mode d'emploi..... 110
(PT) Manual de instruções..... 148	(IT) Manuale d'istruzione..... 183
(NL) Gebruiksaanwijzing..... 220	(TR) Kullanım Talimatları..... 256
(RU) Инструкции по эксплуатации..... 292	(ZH) 使用说明..... 330

Table of Contents

	Page
Photometer Series	
1 Introduction	5
1.1 General Information	5
1.1.1 Read instructions before use	5
1.1.2 Purpose of the Product	5
1.1.3 Authorized use	5
1.1.4 Requirements for safe use	5
1.1.5 Keep instructions	6
1.1.6 User qualification	6
1.1.7 Handling of hazardous chemicals	6
1.1.8 Disposal notes	6
1.2 List of all used signs in the document	6
2 Product overview	8
2.1 Guide to symbols	8
2.2 Features	8
2.3 Product description	8
2.3.1 Instrument view	9
2.3.2 Interface description	9
3 Commissioning	10
3.1 Operating environment	10
3.2 Contents of delivery	10
3.3 Inserting and replacing the batteries	11
3.4 Mounting	12
3.4.1 Mechanical installation	12
3.4.2 Internet Updates	12
4 Operation	13
4.1 First start up	13
4.2 General operation principles	13
4.3 Control elements	14
4.3.1 List of control elements and their function	14
4.4 Description of the operation procedure or the individual sections of the user interface	15
4.4.1 Mode Selection	15
4.5 Description of advanced operation	20
4.5.1 Calibration special methods - Calcium Hardness Method 191 – Calibration of a method blank	20
4.5.2 Calibration special methods - Fluoride Method 170	21
4.5.3 Calibration special methods - PTSA Method 500	22
4.5.4 Calibration special methods - PTSA 2P Method 501	22
4.5.5 Calibration special methods - Fluorescein Method 510	23

	Page
4.5.6 Calibration special methods - Fluorescein 2P Method 511	23
4.5.7 User Calibration	24
4.5.8 User Methods	25
User Concentration Methods	25
User Polynomials	26
4.5.9 One Time Zero (OTZ)	27
4.5.10 Profi-Mode	27
5 Troubleshooting	28
5.1 Error messages and what to do	28
5.2 Problem solving strategies	29
6 Accessories & Replacement parts	31
6.1 List of Accessories	31
7 Specifications	33
8 Appendix	35
8.1 Literature	35
8.2 Copyright and Trademark Notice	35

1 Introduction

1.1 General Information

1.1.1 Read instructions before use

This manual provides important information on the safe operation of the product. Please read this manual carefully and familiarize yourself with the product before use.

1.1.2 Purpose of the Product

Lovibond photometers are designed to measure different parameters in aqueous samples in a variety of applications such as e.g. drinking water, waste water, industrial processing water, pool water and science & research.

Lovibond photometers are suitable to be used in laboratory settings as benchtop as well as in the field as portable instrument. For portability instruments are delivered in robust carrying cases with required accessories.

Performance of the instruments can be impacted by exposure to extreme light and temperature that is why the instruments should always be used under recommended environmental conditions.

1.1.3 Authorized use

The manufacturer's liability and warranty for damage is voided with improper use, failure to follow this manual, use by unqualified personnel, or unauthorized changes to the product.

The manufacturer is not liable for costs or damages that arise from the user or third parties due to the use of this product, especially in cases of improper use of the product or misuse or faults in the connection of the product.

The manufacturer assumes no liability for print errors.

1.1.4 Requirements for safe use

Note the following points for safe use:

- If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.
- The product may only be used according to the authorized use specified above.
- The product may only be supplied with power by the energy sources mentioned in this operating manual.
- The product may only be used under the environmental conditions mentioned in this operating manual.
- The product must not be opened or modified.

The product must not be used if:

- it is visibly damaged (e.g. after being transported)
- it was stored under adverse conditions for a lengthy period of time (storage conditions, see chapter "Specifications")

1.1.5 Keep instructions

The manual must be kept in the vicinity of the product so you can always find the information you need.

1.1.6 User qualification

The operating personnel must be able to understand and correctly implement the safety labels and safety instructions on the packages and inserts of the products.

The user must be capable and able to read and understand this manual in order to familiarize themselves with the handling and to ensure safe use.

1.1.7 Handling of hazardous chemicals

Chemical and/or biological hazards may exist where this product is used. Abide by all governing laws, regulations and protocols when using this product.

For the development of products, Lovibond® pays close attention to safety. Some hazards from dangerous substances cannot be avoided. If self-produced tests or solutions are used, the responsibility concerning any risks caused by those tests or solutions lies with the user (personal responsibility).

1.1.8 Disposal notes

Dispose of the batteries and electrical devices at a suitable facility according to local legal requirements.

It is illegal to dispose of the batteries with household waste.

Within the European Union, the batteries are removed at a specialized treatment center at the instrument's end of life.

Instruments marked with this symbol must not be disposed of in normal domestic waste.

1.2 List of all used signs in the document

The following symbols are used in this manual to identify sections that require special attention:

Danger!

A hazard exists that will result in death or severe injury if not avoided.

Warning!

Improper handling of certain reagents can cause damage to your health. In any case follow the safety labels on the packing, the safety instructions of the package insert and available SDS. Protective measures specified there have to be followed exactly.

Caution!

A hazard exists that may result in minor or moderate injury.

Notice!

Important information or specific instructions need to be strictly followed.

EN

2 Product overview

2.1 Guide to symbols

Labels attached to the product should be strictly observed to avoid personal injury or damage to the product. Refer to this chapter for information regarding the nature of the danger or risk before taking any action where such label is present.

For professional users in the European Union:

If you wish to discard electrical and electronic equipment (EEE), please contact your dealer or supplier for further information.

For disposal in countries outside of the European Union:

This symbol is only valid in the European Union (EU). If you wish to discard this product please contact your local authorities or dealer and ask for the correct method of disposal.

2.2 Features


Lovibond MD6x0 and PM6x0 photometers provide high level of accuracy and efficiency with user friendly interface.

- Robust design, water & dust proof
- 5 nm interference filter for higher precision
- Updateable to include the latest methods and firmware version
- Multilingual UI for global use
- Memory capacity up to 1000 data sets for standard and 500 sets for Bluetooth versions
- Data transfer capability
- Up to 35 User programmable methods
- PTSA and Fluorescein capability (MD 640)

2.3 Product description

Lovibond's MD & PM 6xx series are multiparameter and multiwavelength photometers designed specifically to include all important parameters in their respective application. With high quality and long term stable LEDs as light source, instruments provide high accuracy in measurement results and operational efficiency. Depending on variant, the instruments measure at 3 to 6 wavelengths of 430, 530, 560, 580, 610 and 660 nm. To insure reproducible results photometers use high precision interference filters with a FWHM of 5 nm. Instruments are operated by 4 AA batteries and have a compact design which makes the perfect handheld instruments ideal in the fields. MD6x0 series offer up to 120 and PM6x0 series offer up to 32 preprogrammed methods and with the help of locally produced reagents and test standards, they provide a complete solution to meet requirements in the water analysis industry.

2.3.1 Instrument view

Front Side

1 - Sample Chamber

2 - Display

3 - Model Number

4 - Power Button

EN

2.3.2 Interface description

MD 610 / MD 640 / PM 630

The MD 610 / MD 640 / PM 630 has a Bluetooth® 4.0 interface which enables the wireless transmission of data. Now it is possible to transmit current readings automatically and manually. Stored results can also be shared manually. Bluetooth® 4.0 is also known as Bluetooth® Smart or Bluetooth® LE (Low Energy). Data is transmitted from the photometer as a .csv file. Details on how information is transmitted from the photometer can be found on www.lovibond.com. To receive the data, there are several options on offer from the Tintometer® Group.

The App, AquaLX®, is available for mobile devices such as Smartphones and Tablets and enables the user to manage and graphically chart the received data. Both the data and charts can then be shared via email. AquaLX® can be downloaded free of charge from the iTunes Store® for iOS® and from Google Play™ Store for Android™.

A software tool is available for PCs to receive data stored on the photometer. The data can be exported to an Excel® spreadsheet which enables users to process the information according to their usual practice. If Excel® is not available, the data can be stored as a .txt file for processing at a later date. A Bluetooth® dongle is required to receive the data. This is included in the standard shipment.

Bluetooth® Module - Specifications:

- Modul: BLE 113-A
- Bluetooth® 4.0 LE
- FCC ID: QQQBT113
- IC: 5123A-BGTBLE113

3 Commissioning

3.1 Operating environment

The instrument can be used in any indoor or outdoor environment that is deemed safe for an operator to perform the analysis. However, it is recommended not to use the instrument in bright sunlight to avoid errors caused by stray light. This environment includes a temperature range from 5 to 40 degrees Celsius with a non-condensing humidity up to 90%. Large temperature differences between the instrument and the environment can lead to errors – e.g. due to the formation of condensation in the area of the lens or on the vial. The instrument should be placed on a level surface to achieve the best measurement performance.

3.2 Contents of delivery


Carefully inspect all items to ensure that every part of the list below is present and no visible damage has occurred during shipment. If there is any damage or something is missing, please contact your local distributor immediately.

Content	MD 600	MD 610	MD 640	PM 600	PM 620	PM 630
Instrument in case with inserts	✓	✓	✓	✓	✓	✓
4 Batteries (AA)	✓	✓	✓	✓	✓	✓
3 Vials 24 mm Ø	✓	✓	✓	✓	✓	✓
3 Vials 16 mm Ø	✓	✓	✓			
1 Adapter each (16 mm Ø and 13 mm Ø cuvettes)	✓	✓	✓			
Plastic stirring rod 13 cm	✓	✓	✓	✓	✓	✓
Brush 11 cm	✓	✓	✓	✓	✓	✓
Plastic beaker 100 ml				✓	✓	✓
Syringe			✓	✓	✓	✓
Screw driver	✓	✓	✓			
Instruction manual	✓	✓	✓	✓	✓	✓
Methods manual (print form)				✓	✓	✓
Methods manual (USB)	✓	✓	✓			
Certificate of Compliance	✓	✓	✓	✓	✓	✓
Warranty Information	✓	✓	✓	✓	✓	✓

Reagents for:	MD 600	MD 610	MD 640	PM 600	PM 610	PM 630
Chlorine (free, combined, total)				✓	✓	✓
pH value				✓	✓	✓

Reagents for:	MD 600	MD 610	MD 640	PM 600	PM 610	PM 630
Calcium Hardness				√	√	√
Alkalinity-M				√	√	√

3.3 Inserting and replacing the batteries

- (1) Screwdriver
- (2) Instrument back
- (3) Battery
- (4) Notch
- (5) Screw
- (6) Seal ring
- (7) Battery compartment cover
- (8) RJ 45 Connector for updates

EN

1. Switch the instrument off.
2. If necessary remove vial from the sample chamber.
3. Place the instrument upside down on a clean and even surface.
4. Unscrew the four screws (5) of the battery compartment cover (7).
5. Lift off battery compartment cover (7) at the notch (4).
6. Remove old batteries (3).
7. Place 4 new batteries.
Ensuring the correct polarity!
8. Replace the battery compartment cover (7).
Check the seal ring (6) of the notch to make sure it is tight-fitting
9. Tighten the screws (5) carefully.

**Caution!**

To ensure that the instrument is water proof:

- Seal ring must be in position
- Battery compartment cover must be fixed with all screws

**Caution!**


Dispose of used batteries in accordance with all federal, state and local regulations.

**Caution!**

The batteries save data (stored results and photometer setting). During battery change the data in the photometer is saved for 2 minutes. If the change time exceeds 2 minutes all stored data and settings are lost. Recommendation: for replacement a screwdriver and new batteries must be available.

3.4 Mounting**3.4.1 Mechanical installation**

Insertion of the adapter:

**3.4.2 Internet Updates**

To connect the instrument to the serial interface of a computer the optional connection cable with integrated electronic system is required.

It is possible to update instrument firmware via the internet. Please find detailed information about current firmware under the download section of each instrument version.

How to open and close the battery compartment cover see chapter: Inserting and replacing the batteries.

**Notice!**

To prevent loss of stored test results store or print them out before performing an update. If the update procedure is interrupted (eg. interruption of connection, LoBat., etc.) the instrument isn't able to work (no display). The instrument will only work again after completing the data transfer.

4 Operation

4.1 First start up


Before using the photometer perform the following settings in the Mode-Menu:

- MODE 10: Select language
- MODE 12: Set date and time
- MODE 34: Perform „Delete data“
- MODE 69: Perform “User m. init” to initialise the userpolynomial system

Please refer to chapter "Mode Selection" for instructions.

EN

4.2 General operation principles





1. Interference filter
2. Detector
3. Display
4. Sample compartment
5. Light source LED
6. Light beam

The photometer has pre-programmed methods which are derived from standard analytical procedures. To ensure simplified and error free analysis, necessary calibration curves with reagents, reaction times and sequences are all programmed into these methods.

The optical setup of photometer consists of light source, interference filters and photo sensors. Lovibond photometers use LEDs as light source which are highly energy efficient and have a long life to last as long as the life of the photometer itself. High quality interference filters are used to limit the wavelength to ensure high precision in measurement results. Photo sensors collect the signal and transfer it to the microprocessor which digitally calculates concentration and displays results in respective units.

4.3 Control elements

4.3.1 List of control elements and their function

Keys	Description / Function	How to execute
	Time and date, display timer setup	Press key once and the display shows date and time. Press again and the display shows last used user countdown, press Enter to use existing countdown or enter a new one. Press Enter to start countdown.
	Backlight on/off	Press Shift + F1 keys to turn backlight on or off. It is switched off automatically during measurement.
 or 	Method Selection	The display shows methods list after start. Choose method by scrolling with Arrow keys and confirming with Enter key or enter method number by pressing Shift + method Number keys. To measure absorbance please select method associated with each wavelength at the end of method list.
	Method Information	Press F1 key to see method information such as name, range, reagent form, vial type and size and reagents used. Press F1 key again to revert back to methods list.
	Chemical Species Information	Press F2 key to see chemical and corresponding range. Press F2 key again to revert back to methods list. You can change the chemical species of some methods by scrolling with arrow keys once the result is displayed.
	Differentiation Selection	Press Arrow keys to select required determination and press Enter key to confirm.
	Perform Zero	Prepare a vial according to method description in methods manual. Place vial in the sample chamber making sure that the marks of instrument and vial are aligned. Press Zero key to perform zero.
	Test measurement	Prepare a vial according to method description in methods manual. Place vial in the sample chamber making sure that the marks of instrument and vial are aligned. Press Test key to perform measurement.

Keys	Description / Function	How to execute
 and 	Execute integrated countdowns	Press Enter key to start first countdown before placing vial in the sample chamber. Place vial in the sample chamber and press Test key to start 2nd countdown. The instrument will perform measurement after the countdown is finished. It is possible to cancel countdown and start direct measurement, by pressing Enter key again.
	Store test results	Press Store key when results are displayed. Enter 6-digit code and press Enter key to confirm. The user is prompted when memory capacity is near it's limit.
	Print test results	When printer is connected through IRiM module (optional) press F3 key to print results. The printed data set contains, date, time, code no, method and test results.

Other keys and their functions

	Mode menu	Opens mode functions list.
	Power on/off	Turns Instrument on or off.
	Shift key	Must be pressed to activate secondary / numeric keypad.
	Escape key	Cancels actions and takes to the previous menu.

4.4 Description of the operation procedure or the individual sections of the user interface

4.4.1 Mode Selection

- Press Mode key to open the mode functions list.
- Select required mode function from list or press Shift + mode Number keys then press Enter key to confirm.

Mode Number	Description / Display	How to Execute function
10	Language	<ul style="list-style-type: none"> • Press Arrow keys to select desired language. • Press Enter key to confirm.
11	Key beep	<ul style="list-style-type: none"> • Press Shift + 0 keys to turn off key beep. • Press Shift + 1 keys to turn on key beep. • Press Enter key to confirm.
12	Date/Time	<ul style="list-style-type: none"> • Enter/Adjust date and time according to the format on display. • While confirming date and time the seconds are adjusted to zero automatically.
13	Countdown	<ul style="list-style-type: none"> • Press Shift + 0 keys to switch off countdown. • Press Shift + 1 keys to switch on countdown. • Press Enter key to confirm. <p>Note: It is possible to skip countdown by pressing Enter key. If the countdown function is switched off, the operator is responsible for ensuring the necessary reaction period. Non-compliance with reaction periods leads to incorrect test results.</p>
14	Signal beep	<ul style="list-style-type: none"> • Press Shift + 0 keys to switch off signal beep. • Press Shift + 1 keys to switch on signal beep. • Press Enter key to confirm.
18	Bluetooth *	<ul style="list-style-type: none"> • Press Shift + 0 keys to switch off Bluetooth connection. • Press Shift + 1 keys to switch on Bluetooth connection. • Press Enter key to confirm.
19	Auto transfer*	<p>Autotransfer enables user to transfer measurement results automatically to the App or PC without storing.</p> <ul style="list-style-type: none"> • Press Shift + 0 keys to switch off auto transfer. • Press Shift + 1 keys to switch on auto transfer. • Press Enter key to confirm.
20	Print / Transfer	<ul style="list-style-type: none"> • Press Enter key to start print/transfer all stored results, press Esc key to cancel mode.
21	Print / Transfer, date	<ul style="list-style-type: none"> • Enter start date and end date to print/transfer data from specified time period and press Enter key to start print / transfer.

Mode Number	Description / Display	How to Execute function
22	Print / Transfer, code no.	<ul style="list-style-type: none"> Specify code numbers range by entering first and last code number. Press Enter key to start print/transfer. Enter same code in both fields to print/transfer one data set only. Enter 0 in code no. twice to print/transfer results without code.
23	Print / Transfer, method	<ul style="list-style-type: none"> Choose method from list or enter method number. In case of differentiated methods select the required determination and press Enter key to confirm. To start press Enter key.
29	Printing parameters	<ul style="list-style-type: none"> Display shows Baud rate, press Shift + 2 keys to enter Baud rate selection menu. Select Baud rate from list and press Enter key to save. Press Enter key again to execute.
30	Storage	<ul style="list-style-type: none"> Press Enter key to display all results in chronological order starting from latest result. Press Arrow keys to go to next or previous result. Press F3 key to print/transfer displayed result only. Press F2 key to print/transfer all stored results.
31	Stor., date	<ul style="list-style-type: none"> Enter start date and end date to specify time period and press Enter key to view results of selected time period. Press F3 key to print/transfer displayed result only. Press F2 key print/transfer all stored results.
32	Stor., code	<ul style="list-style-type: none"> Enter first required code no. and last required code no. to specify code no. range and press Enter key to view results of selected code no. range. Press F3 key to print/transfer displayed result only. Press F2 key print/transfer all stored results.
33	Stor., method	<ul style="list-style-type: none"> Select method from list or enter method number, press Enter key to view results of selected method. Press F3 key to print/transfer displayed result only. Press F2 key print/transfer all stored results.
34	Delete data	<ul style="list-style-type: none"> Press Shift + 0 keys to retain data. Press Shift + 1 keys to delete data. Press Enter key to confirm. <p>Attention: All stored results are deleted. Press Esc key to cancel without deleting data.</p>

Mode Number	Description / Display	How to Execute function
40	Calibration	<ul style="list-style-type: none"> • Select method from list for calibration and press Enter key to confirm. • Refer to chapter 4.5. for detailed procedures. • Press Shift + 2 keys to reset method blank calibration of Method 191 to factory Calibration.
45	User Calibration	<ul style="list-style-type: none"> • Select mode when test results are displayed. • Press Arrow keys to increase or decrease test results to match the corresponding Standard values. • Press Enter key to save the results. • Refer to chapter 4.5. for more information.
46	Clear Calibration	<ul style="list-style-type: none"> • Select mode before zero measurement of required method. • Press Shift + 0 keys to retain user calibration. • Press Shift + 1 keys to delete user calibration.
50	Profi Mode	<ul style="list-style-type: none"> • Press Shift + 0 keys to switch off profi mode. • Press Shift + 1 keys to switch on profi mode. • Press Enter key to confirm. <p>Refer to chapter 4.5. for more details.</p>
55	One Time Zero	<ul style="list-style-type: none"> • Press Shift + 0 keys to switch off OTZ. • Press Shift + 1 keys to switch on OTZ. • Press Enter key to confirm. <p>Refer to chapter 4.5. for more details.</p>
60	Method list	<p>This setting allows user to create a user defined list of favorite methods. The program structure requires that the list must have at least one active method. For this reason, it is necessary to activate all methods and then deactivate the ones which are not required.</p> <ul style="list-style-type: none"> • Press Enter key to display favourite methods. • Press Arrow keys to select desired method. • Press F2 key to activate or inactivate specific method and press Enter key to confirm. • Methods with [*] behind method number will be displayed only.
61	M list all on	<ul style="list-style-type: none"> • Press Shift + 1 keys to switch all user methods on. • Press Shift + 0 keys to keep valid methods list.
62	M list all off	<ul style="list-style-type: none"> • Press Shift + 1 keys to display only one method on the list. • Press Shift + 0 keys to keep valid methods list.
64	User concentration	<ul style="list-style-type: none"> • Create a user method by entering concentration values. • Refer to chapter 4.5. for detailed instructions.

Mode Number	Description / Display	How to Execute function
65	User polynoms	<ul style="list-style-type: none"> • Create a user method by entering polynomial data. • Refer to chapter 4.5. for detailed instructions.
66	User methods clear	<ul style="list-style-type: none"> • Enter user method number which is to be deleted. • Press Shift + 1 keys to delete user method. • Press Shift + 0 keys to keep user method.
67	User methods print / transfer	<ul style="list-style-type: none"> • Press Enter key to start print/transfer all user methods data.
69	User methods init	<ul style="list-style-type: none"> • Press Enter key to confirm user method intializing. • Press Shift + 1 keys to start user method. • Press Shift + 0 keys to cancel without initialisation. <p>Attention: All stored user methods are deleted with initialisation.</p>
70	Langelier	<ul style="list-style-type: none"> • Enter the temperature value (T) in the range between 3 and 53 °C or 37 and 128 °F if °F was selected. • Enter the value for Calcium hardness in the range between 50 and 1000 mg/l CaCO₃. • Enter the value for Total Alkalinity in the range between 5 and 800 mg/l CaCO₃. • Enter the value for TDS in the range between 0 and 6000 mg/l. • Enter the pH-value in the range between 0 and 12. • The display shows the Langelier Saturation Index. • Press Enter key to start new. • Press ESC key to return to Mode Menu.
71	Temperature	<ul style="list-style-type: none"> • Press Shift + 1 keys to select degree Celsius. • Press Shift + 2 keys to select degree Fahrenheit.

Mode Number	Description / Display	How to Execute function
80	LCD Contrast	<ul style="list-style-type: none"> • Press Arrow keys to increase and decrease contrast. • Press Store key to increase contrast by 10 units and Test key to decrease by 10 units. • Press Enter key to confirm.
81	LCD	<ul style="list-style-type: none"> • Press Arrow keys to increase and decrease brightness. • Press Store key to increase brightness by 10 units and Test key to decrease by 10 units. • The brightness can be selected between 0 and 254 units. • Press Enter key to confirm.
91	System info	<ul style="list-style-type: none"> • Display shows software version. • Press Arrow key to view number of performed tests and free memory capacity.

* Only MD610, MD640 and PM630

4.5 Description of advanced operation

4.5.1 Calibration special methods - Calcium Hardness Method 191 – Calibration of a method blank

- Method Selection
 - After selection of calibration mode 40, select method 191 by pressing shift + 1 key or by scrolling with arrow keys in case of MD 640.
 - Press enter key to confirm.
- Zero
 - Perform zero with exactly 10 ml of deionised water in a 24 mm Ø vial.
- Sample Preparation
 - Pipette 100 ml of water free calcium into a beaker.
 - Add 10 CALCIO H No. 1 tablets straight from the foil and dissolve them completely by crushing with stirring rod.
 - Add 10 CALCIO H No. 2 tablets straight from the foil to the same water, dissolve them completely by crushing with stirring rod.
- Sample
 - Press Enter key to start the test countdown.
 - After countdown is finished rinse the vial (24 mm Ø) with the coloured sample from the beaker and fill with 10 ml of the sample.
 - Press Test key.
 - The batch related method blank is saved.
 - Press Enter key, to go back to Mode menu.

5. Notes
- If a new batch of CALCIO tablets is used a calibration of the method blank has to be performed to optimise the results.
 - Use deionised or tap water.
 - If no water free Calcium is available these ions can be masked by using EDTA.
Preparation: Add 50 mg (a spatula tip full) EDTA to 100 ml water and dissolve.
 - To achieve the most accurate method blank it is important to adhere exactly to the sample volume of 100 ml.

4.5.2 Calibration special methods - Fluoride Method 170

1. Method Selection
- After selection of calibration Mode 40, select method 170 by pressing Shift + 3 or by scrolling with Arrow keys in case of MD 640.
 - Press Enter key to confirm.
-
2. Zero
- Perform zero with exactly 10 ml of deionised water in a 24 mm Ø vial.
-
3. Sample Preparation & Test
- Add exactly 2 ml SPADNS reagent solution to the water sample.
Caution: Vial is filled up to the top!
 - Close the vial tightly with the cap and swirl gently several times to mix the contents.
 - Place the vial in sample chamber and press Test key.
-
4. Standard Preparation & Test
- Remove the vial, empty it, rinse vial and cap several times and then fill the vial with exactly 10 ml Fluoride standard (Concentration 1 mg/l F).
 - Add exactly 2 ml SPADNS reagent solution to the Fluoride standard. Caution: Vial is filled up to the top!
 - Place the vial in the sample chamber and Press Test key.
 - Press Enter key to confirm the displayed results.
-
5. Calibration Check
- Press Esc key to exit the calibration mode.
 - Select Fluoride Method from methods list and if error (Error, absorbance T2>T1) appears, repeat the adjustment.
-
6. Notes
- The same batch of SPADNS reagent solution must be used for adjustment and test. The adjustment process needs to be performed for each new batch of SPADNS reagent solution (see Standard methods 20th, 1998, APHA, AWWA, WEF 4500 F D., S. 4-82).
 - As the test result is highly dependent on exact sample and reagent volumes, the sample and reagent volumes should always be metered by using a 10 ml resp. 2 ml volumetric pipette (class A).

4.5.3 Calibration special methods - PTSA Method 500

1. Method Selection
 - After selection of calibration mode 40, select PTSA method 500 by scrolling with Arrow keys.
 - Press Enter key to confirm.

-
2. Procedure

Perform following tests:

 - **T1 - Blank**
 - Place the 0 ppb standard vial (blank) from the PTSA calibration set in the sample chamber and press Test key.
 - The display shows "T1 accepted".
 - **T2 - 200 ppb**
 - Place the 200 ppb standard vial from PTSA calibration set in the sample chamber and press Test key.
 - The display shows "T2 accepted".
 - **T3 - 1000 ppb**
 - Place 1000 ppb standard vial from PTSA calibration set in the sample chamber and press Test key.
 - The display shows "Calibration accepted".
 - Confirm with Enter key to save the calibration.

-
3. Reagents

Reagent: PTSA calibration set
Form of reagent/Quantity: 3 x vials (0, 200, 1000 ppb)
Order-No.: 461245

4.5.4 Calibration special methods - PTSA 2P Method 501

1. Method Selection
 - After selection of calibration mode 40, select PTSA 2P method 501 by scrolling with Arrow keys.
 - Press Enter key to confirm.

-
2. Procedure

Perform following tests:

 - **T1 - Blank**
 - Fill a clean vial (24 mm Ø) with 10 ml of deionised water or 0 ppb PTSA standard solution, close tightly with a black cap and place it in the sample chamber.
 - Press Test key.
 - The display shows "T1 accepted".
 - **T2 - (50 ... 400) ppb**
 - Press Enter key to save predetermined concentration or enter a concentration in the range from 50 to 400, press Enter key again to confirm.
 - Place the vial with standard solution of set concentration and press Test key.
 - The display shows "Calibration accepted".
 - Press Enter key to save the calibration.

4.5.5 Calibration special methods - Fluorescein Method 510

- | | |
|---------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1. Method Selection | <ul style="list-style-type: none"> • After selection of calibration mode 40, select Fluorescein method 510 by scrolling with Arrow keys. • Press Enter key to confirm. |
| 2. Procedure | <p>Perform following tests:</p> <ul style="list-style-type: none"> • T1 - Blank • Fill a clean vial (24 mm Ø) with 10 ml of deionised water or 0 ppb Fluorescein standard solution (blank), close tightly with a black cap and place it in the sample chamber. • Press Test key. • The display shows "T1 accepted". • T2 - 75 ppb • Remove and empty vial completely then fill with 10 ml 75 ppb Fluorescein standard solution. • Place the vial in the sample chamber and press Test key. • The display shows "T2 accepted". • T3 - 400 ppb • Remove and empty vial completely then fill with 10 ml 400 ppb Fluorescein standard solution. • Place the vial in the sample chamber and press Test key. • The display shows "Calibration accepted". • Confirm with Enter key to save the calibration. |
| 3. Reagents | <ol style="list-style-type: none"> 1. Reagent: Fluorescein standard addition solution 400 ppb form of reagent/Quantity: solution / 50 ml
Order-No.: 461230 2. Reagent: Fluorescein calibration set form of reagent/Quantity: 2 x 50 ml 0 ppb, 2 x 50 ml 75 ppb, 1 x 50 ml 400 ppb
Order-No.: 461240 |

4.5.6 Calibration special methods - Fluorescein 2P Method 511

- | | |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1. Method Selection | <ul style="list-style-type: none"> • After selection of calibration Mode 40, select Fluorescein 2P Method 511 by scrolling with Arrow keys. • Press Enter key to confirm. |
| 2. Procedure | <p>Perform following tests:</p> <ul style="list-style-type: none"> • T1 - Blank • Fill a clean vial (24 mm Ø) with 10 ml of deionised water or 0 ppb Fluorescein standard solution, close tightly with a black cap and place it in the sample chamber. • Press Test key. • The display shows "T1 accepted". |

- **T2 - (010 ... 300) ppb**
- Press Enter key to save predetermined concentration or enter a concentration in the range from 10 to 300, press Enter key again to confirm.
- Place the vial with standard solution of set concentration and press Test key.
- The display shows "Calibration accepted".
- Press Enter key to save the calibration.

EN

4.5.7 User Calibration

If a test method is user calibrated the method name is displayed inverse.

Procedure:

- Prepare a standard of known concentration and use this standard instead of the sample according to the test procedure.
- It is recommended to use well known standards which are formulated according to DIN EN, ASTM or other international norms or to use certified standards which are commercially available.
- After measuring this standard solution it is possible to change the displayed results to the required value.
- If a method uses a mathematic equation for the calculation of the result, it is only possible to calibrate the basic tests since all the other tests use the same polynomial.
- The same applies for some test procedures which use a polynomial from another test procedure.

Return to factory calibration:

If the user calibration is deleted, the factory calibration is automatically activated.

Remarks:

The method "Fluoride" cannot be calibrated with Mode 45 since the test requires a calibration related to the batch of the liquid reagent (SPADNS) (Mode 40, chapter "Calibration special methods - Fluoride Method 170").

The recommended range for user calibration is between the lower end and the upper end of the measuring range (25% - 75%).

4.5.8 User Methods

User Concentration Methods

It is possible to enter and store up to 10 User Concentration Methods.

Therefore you need 2 to 14 standards of known concentration and one blank (deionised water or reagent blank value). The Standards should be measured with increasing concentrations and from the brightest to the darkest coloration.

The measuring range for "Underrange" and "Overrange" is defined with -2600 mAbs^* and $+2600 \text{ mAbs}^*$. After selection of a method the concentration of the lowest and highest used standard is displayed as measuring range. The operation range should be within this range to achieve best results.

*1000 mAbs = 1 Abs = 1 E (displayed)

EN

Method Entry Procedure:

1. Select Mode 64 and enter a method number in the range from 850 to 859, e.g.: Shift + 850 and confirm with Enter key.
Note: If the entered number has already been used to save a concentration the display shows the query to overwrite:
Press Shift + 0 or ESC keys to go back to method no. query.
Press Shift + 1 keys to overwrite.

2. Choose the required wavelength e.g.: Shift + 2 keys for 560 nm.

3. Press Arrow keys to select the required unit and confirm with Enter key.

4. Press the appropriate numerical key to select the required resolution, e.g.: Shift + 4 keys for [0][.][0][1].
Note: Please enter the required resolution according to the instrument pre-sets:

Range	Max. resolutions
0.000 ...9.999	0.001
10.00 ...99.99	0.01
100.0... 999.9	0.1
1000 ...9999	1

Measurement procedure with standards of known concentration:

1. Perform zero with deionised water or reagent blank.

2. Enter the concentration of the first standard; e.g.: Shift + [0][.][0][5] keys
 - Press ESC key to go one step back.
 - Press F1 key to reset numerical input.

3. Confirm with Enter key.

4. Prepare the first standard, place it in the vial chamber and press Test key.

5. The display shows the input value and the measured absorption value. Confirm with Enter key.

6. Enter the concentration of the second standard; e.g.: Shift + [0][.][1] keys

7. Prepare the second standard, place it in the vial chamber and press Test key.

8. The display shows the input value and the measured absorption value. Confirm with Enter key.
Note:
 - Perform as described above to measure further standards.
 - The minimum of measured standards is 2.
 - The maximum of measured standards is 14 (S1 to S14).

9. If all required standards or the maximum value of 14 standards are measured press Store key.

10. The concentration is stored and the instrument goes back to mode menu. Now the concentration is stored in the instrument and can be recalled by entering its method number or selecting it from the displayed method list.

User Polynomials

It is possible to enter and store up to 25 User Polynomials.

The program allows the user to apply a Polynomial up to the 5th degree:

$$y = A + Bx + Cx^2 + Dx^3 + Ex^4 + Fx^5$$

If only a Polynomial of a lower degree is necessary the other coefficients are specified as zero (0), e.g.: for the 2nd degree is D, E, F = 0.

The values of the coefficients A, B, C, D, E, F must be entered in an academic notation with maximal 6 decimal places, e.g.: 121,35673 = 1,213567E+02

Method Entry Procedure:

1. Select Mode 65 and enter a method number in the range from 800 to 824, e.g.: Shift + 800 and confirm with Enter key.
Note: If the entered number has already been used to save a polynomial, the display shows the query:
 - Press Shift + 0 or ESC keys to go back to method no. query.
 - Press Shift + 1 keys to start entry mode.

2. Choose the required wavelength (refer to step 2 in user concentration methods).

3. Enter data of the coefficient A including decimal point, e.g.: Shift + [1][.][3][2] keys
Note: Use Arrow keys to change between plus and minus sign.

4. Press F1 key to reset numerical input.

5. Confirm with Enter key.

6. Enter the exponent of the coefficient A, e.g.: Shift + 3, confirm with Enter key.

7. Successively the instrument queries the data for the other coefficients (B, C, D, E and F).
Note:
 - If zero [0] is entered for the value of the coefficient, the input of the exponent is omitted automatically.

8. Confirm every input with Enter key.

9. Enter measurement ranges from – 2600 to +2600 mAbs.

10. Enter the values in Absorbance (mAbs) for the upper limit (Max) and the lower limit (Min).

11. Confirm every input with Enter key.

12. Follow step 3 and 4 from concentration methods for entry of unit and resolution.

13. Once all values are entered the display shows “stored” and the instrument goes back to mode menu automatically. The method can be accessed from list of methods or by entering method number.

4.5.9 One Time Zero (OTZ)

One Time Zero is available for all methods where Zero is performed in a 24 mm Ø round vial with sample water.

One Time Zero can be used for different tests providing the tests are performed with the same sample water and under the same test conditions.

When the instrument is first being used for an OTZ compatible method and One Time Zero is activated, the instrument will request a new Zero with “prepare OT-Zero”. Perform Zero as described in the method. This Zero will be stored and used for all methods with OTZ function until the instrument is switched off.

If necessary, a new Zero can be performed by pressing Zero key at any time.

Note:

The specified accuracy is valid for all test results when Zero is performed for each test (One Time Zero function is switched off).

4.5.10 Profi-Mode

This function may be used for routine analyses with many samples of one method. The following information is always stored in the methods:

1. Method
2. Range
3. Date and time
4. Differentiation of results
5. Detailed operator instruction
6. Compliance with reaction periods

If the Profi-Mode is active, the photometer provides only a minimum of operator instructions. The criteria specified above in 4, 5, 6 are no longer included.

Note:

Storage of test results is possible. When results are stored the display also shows “Profi-Mode”. The selected settings are kept by the photometer even when it is switched off. To change photometer setting a new setting is required.

5 Troubleshooting

5.1 Error messages and what to do

Display	Possible Causes	Elimination
Overrange	Reading is exceeding the range.	If possible dilute sample or use other measuring range.
	Water sample is too cloudy.	Filtrate water sample.
	Too much light on the photo cell.	Seal on the cap? Repeat measurement with seal on the cap of the vial.
Underrange	Result is under the detection limit.	Indicate result with lower x mg/l; x = low end of measuring range; if necessary use other analytical method.
Storagesystem error use Mode 34	Main power fails or is not connected.	Insert or change battery. Delete data with Mode 34.
Battery warning 	Warning signal every 3 minutes.	Capacity of the battery is too low; change the batteries.
	Warning signal every 12 seconds.	
	Warning signal, the instrument switches itself off.	Change the batteries.
Jus Overrange E4 Jus Underrange E4	The user calibration is out of the accepted range.	Please check the standard, reaction time and other possible faults. Repeat the user calibration.
Overrange E1 Underrange E1	The concentration of the standard is too high/too low, so that during user calibration the limit of the range was exceeded.	Perform the test with a standard of higher/lower concentration.
E40 user calibration not possible	If the display shows Overrange/ Underrange for a test result a user calibration is not possible.	Perform the test with a standard of higher/lower concentration.
Zero not accepted	Light absorption is too great or too low.	Refer to chapter: Performing Zero. Clean sample chamber. Repeat zeroing.

Display	Possible Causes	Elimination
???	The calculation of a value (e.g. combined Chlorine) is not possible.	Test procedure correct? If not – repeat test.
Example 1: 0,60 mg/l free Cl ??? comb Cl 0,59 mg/l total Cl		Example 1: The readings for free and total Chlorine are different, but considering the tolerances of each reading they are the same. For this reason the combined Chlorine is most likely zero.
Example 2: Underrange ??? comb Cl 1,59 mg/l total Cl		Example 2: The reading for free Chlorine is under the detection limit. The instrument is not able to calculate the combined Chlorine. In this case the combined Chlorine is most likely the same as the total Chlorine.
Example 3: 0,60 mg/l free Cl ??? comb Cl Overrange		Example 3: The reading for total Chlorine is exceeding the range. The instrument is not able to calculate the combined Chlorine. The test should be repeated with a diluted sample.
Error absorbance e.g.: T2>T1	Fluoride calibration was not correct.	Repeat calibration.

5.2 Problem solving strategies

Finding	Possible Causes	Elimination
Test result deviates from the expected.	Chemical species not as required.	Press Arrow keys to select the required chemical species.
No differentiation: e.g. for the Chlorine test there is no selection between differentiated, free or total.	Profi-Mode is switched on.	Switch Profi-Mode off with Mode 50.
The preprogrammed countdown is not displayed.	Countdown is not activated and/or the Profi-Mode is activated.	Switch the countdown on with Mode 13 and/or switch the Profi-Mode off with Mode 50.

Finding	Possible Causes	Elimination
It seems that a method is not available.	Method is not activated in the user method list.	Activate the required method in the user method list with Mode 60.

6 Accessories & Replacement parts

6.1 List of Accessories

Title		Part Number
250 mL bottle, AF 631		375072
Adapter for round cuvettes 13 mm Ø	MD 600, MD 610, MD 640	19802192
Adapter for round cuvettes 16 mm ø		19802190
Batteries (AA), set of 4		1950025
Brush, 11 cm length	MD 600, MD 610, MD 640, PM 630	380230
Cleaning cloth		197635
Cuvette stand for 6 round cuvettes Ø 24 mm		418951
Cuvette stand for 10 round cuvettes Ø 16 mm	MD 600, MD 610, MD 640, PM 630	418957
Fixed price service package for MD600/ MD610	MD 600, MD 610	19802702
Fixed price service package for MD640	MD 640	19802703
Fixed price service package for PM600/ PM620/PM630	PM 600, PM 620, PM 630	19802704
Fluoresceine calibration set (0, 75, 400 ppb)	MD 640	461240
Fluoresceine standard addition solution, 400 ppb	MD 640	461230
Kit of BT data transfer software and Bluetooth Dongle	MD 610, MD 640, PM 630	2444480
Measuring beaker, 100 ml		384801
Membrane filter set for use when preparing samples, 25 membrane filters 0.45 µm, 2 syringes 20 ml	MD 600, MD 610, MD 640	366150
Mixing cylinder, 25 ml	MD 600, MD 610, MD 640	19802650
Multy cuvette-3, set of 12		197605
Plastic funnel with handle (white)		471007
Plastic syringe, 5 ml	PM 600, PM 620, PM 630	366120
PTSA calibration set (0, 200, 1000 ppb)	MD 640	461245
PTSA standard addition solution, 200 ppb	MD 640	461200
PTSA standard addition solution, 1000 ppb	MD 640	461210
Round cuvette 16 mm ø, set of 10		197665

Title		Part Number
Round cuvette 24 mm ø, set of 5		197629
Round cuvette 24 mm ø, set of 12		197620
Rubber seal cap	MD 600, MD 610, MD 640	19802223
Sample cuvettes with black lid, Height 48 mm, ø 24 mm, set of 12	MD 640	197657
Sealing ring for round vials 24 mm ø, set of 12		197626
Service plan - 3 years for MD600/MD610	MD 600, MD 610	19802802
Service plan - 3 years for MD640	MD 640	19802803
Service plan - 3 years for PM600/PM620/PM630	PM 600, PM 620, PM 630	19802804
Stirring rod, 10 cm length		364109
Stirring rod, 10 cm length, set of 10		364130
Stirring rod, 13 cm length		364100
Stirring rod, 13 cm length, set of 10		364120
Update cable		214030
Update cable set with USB/R232 adapter		214031
UV Pen Lamp, 254 nm	MD 600, MD 610, MD 640	400740
Water sampler with bottle 250 mL and lid, AF 631		170500

EN

7 Specifications

Optics	MD 6x0	LEDs, interference filters (IF) and photo sensor in transparent sample chamber. Wavelength range: 430 nm IF $\Delta\lambda = 5$ nm 530 nm IF $\Delta\lambda = 5$ nm 560 nm IF $\Delta\lambda = 5$ nm 580 nm IF $\Delta\lambda = 5$ nm 610 nm IF $\Delta\lambda = 6$ nm 660 nm IF $\Delta\lambda = 5$ nm IF = interference filter
	PM 6x0	LEDs, interference filters (IF) and photo sensor in transparent sample chamber. Wavelength range: 530 nm IF $\Delta\lambda = 5$ nm 560 nm IF $\Delta\lambda = 5$ nm 610 nm IF $\Delta\lambda = 6$ nm IF = interference filter
Wavelength Accuracy		± 1 nm
Photometric Accuracy		2 % FS (T = 20 °C – 25 °C)
Suitable Vials	MD 6x0	<ul style="list-style-type: none"> • Multi vial 10 mm • Round Cuvettes 13 mm • Round Cuvettes 16 mm • Round Cuvettes 24 mm
	PM 600	Round Cuvettes 24 mm
	PM 620	<ul style="list-style-type: none"> • Multi vial 10 mm
	PM 630	<ul style="list-style-type: none"> • Round Cuvettes 16 mm • Round Cuvettes 24 mm
Display		Graphic-display
Interfaces		Infrared
Operation		Acid and solvent resistant, touch-sensitive keypad with audible feedback via integrated beeper
Auto – OFF		Yes
Updates		Software updates via internet
Internal Storage	MD 600	Approx. 1000 data sets
	MD 610	approx. 500 data sets
	MD 640	
	PM 630	
	PM 600	approx. 1000 data sets
	PM 620	

Power Supply	4 batteries (Mignon AA/LR6)
Battery Life Time	approx. 26 h
Beeper	existing
Portability	Benchtop
Environmental Conditions	5 – 40 °C at max. 30 – 90 % Conditions rel. humidity (non condensing)
Protection Class	IP 68
Compliance	CE
Languages User Interface	<ul style="list-style-type: none"> • English • French • German • Indonesian • Italian • Polish • Portuguese • Spanish
Dimensions	95 x 45 x 210 mm
Weight	450 g

**Caution!**

Subject to technical modification!
To ensure maximum accuracy of test results, always use the reagent systems supplied by the instrument manufacturer.

8 Appendix

The precision of Lovibond® Reagent Systems (tablets, powder packs and tube tests) is identical to the precision specified in standard literature such as American Standards (AWWA), ISO etc.

Most of the data referred to in these standard methods relates to standard solutions.

Therefore they are not readily applicable to drinking-, boiler- or waste-water, since various interferences can have a major influence on the accuracy of the method.

For this reason we don't state such potentially misleading data.

Due to the fact that each sample is different, the only way to check the tolerances ('precision') is the Standard Additions Method.

According to this method, first the original sample is tested. Then further samples (2 to 4) are taken and small amounts of a standard solution are added, and further results are obtained. The amounts added range from approximately half, up to double the amount present in the sample itself.

These supplementary results make it possible to estimate the actual concentration of the original sample by comparison.

EN

8.1 Literature

The reagent formulations are based on internationally recognised test methods. Some are described in national and/or international guidelines.

- Deutsche Einheitsverfahren zur Wasser-, Abwasser- und Schlammuntersuchung
- Standard Methods for the Examination of Water and Wastewater; 18th Edition, 1992
- Photometrische Analysenverfahren, Schwedt, Wissenschaftliche Verlagsgesellschaft mbH, Stuttgart 1989
- Photometrische Analyse, Lange / Vejdelek, Verlag Chemie 1980
- Colorimetric Chemical Analytical Methods, 9th Edition, London
- Adapted from Merck, for more information see instructions delivered with the test

8.2 Copyright and Trademark Notice

The Bluetooth® word mark is a registered trademark owned by Bluetooth SIG, Inc. and any use by The Tintometer® Group is under license.

iOS® is a registered trademark of Cisco, Inc. and licensed to Apple, Inc.

iTunes Store® is a trademark of Apple, Inc., registered in the U.S. and other countries.

Android™ and Google Play™ are trademark of Google, Inc.

Excel® is a trademark of Microsoft Corp., registered in the U.S. and other countries.

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Trademarks of the Tintometer Group
Tintometer Inc.