# Lovibond® Water Testing

Tintometer® Group

## BOD-System BD 600

### BD 600

<b>(EN)</b> Instruction manual..... 3	<b>(DE)</b> Gebrauchsanleitung..... 39
<b>(ES)</b> Manual de instrucciones.....76	<b>(FR)</b> Mode d'emploi..... 114
<b>(PT)</b> Manual de instruções..... 152	<b>(IT)</b> Manuale d'istruzione..... 190
<b>(NL)</b> Gebruiksaanwijzing.....228	

# Table of Contents

	Page
BOD-System BD 600	
<b>1 Introduction</b>	<b>5</b>
1.1 General Information	5
1.1.1 Read instructions before use	5
1.1.2 Purpose of the Product	5
1.1.3 Authorized use	5
1.1.4 Requirements for safe use	5
1.1.5 Keep instructions	5
1.1.6 User qualification	5
1.1.7 Handling of hazardous chemicals	6
1.1.8 Disposal notes	6
1.2 List of all used signs in the document	6
<b>2 Product overview</b>	<b>8</b>
2.1 Guide to symbols	8
2.2 Certification	8
2.3 Features	8
2.4 Product description	9
2.4.1 Description of the keys and switches	9
2.4.2 Interface description	10
<b>3 Commissioning</b>	<b>11</b>
3.1 Operating environment	11
3.2 Contents of delivery	11
3.3 Inserting and replacing the batteries	12
<b>4 Operation</b>	<b>14</b>
4.1 First start up	14
4.2 General operation principles	14
4.3 List of control elements and their function	17
4.4 Main menu	18
4.5.1 Start test series	18
4.5.2 Show current values	19
4.5.3 Display test series	19
4.5.4 Export test series	21
4.5.5 Options	23
4.6 Advanced operations	24
4.7.1 Notice	26
4.7.2 Initial Steps	26
4.7.3 Selection of the sample volume	26
4.7.4 Sample Preparation	26
4.8 Description of setup and handling of important accessories	28

	<b>Page</b>
<b>5 Maintenance</b>	30
5.1 Calibration	30
5.2 Decommissioning	30
5.3 Maintenance and cleaning	30
<b>6 Troubleshooting</b>	31
6.1 Error Messages and Notices	31
6.2 Interpretation of BOD curve and Problem Solving Strategies	33
<b>7 Accessories &amp; Replacement parts</b>	34
7.1 List of Accessories	34
8.1 Specifications - Mains adapter	35
8.2 Specifications - Stirrer unit	35
8.3 Specifications - BD 600	36
<b>9 Appendix</b>	38
9.1 List of Trademarks	38

# 1 Introduction

## 1.1 General Information

### 1.1.1 Read instructions before use

This manual provides important information on the safe operation of the product. Please read this manual carefully and familiarize yourself with the product before use.

### 1.1.2 Purpose of the Product

The Lovibond® sensor system BD 600 is a 6 sample system that allows precise measurements of biological oxygen demand (BOD) based on the manometric principle.

### 1.1.3 Authorized use

The manufacturer's liability and warranty for damage is voided with improper use, failure to follow this manual, use by unqualified personnel, or unauthorized changes to the product.

The manufacturer is not liable for costs or damages that arise from the user or third parties due to the use of this product, especially in cases of improper use of the product or misuse or faults in the connection of the product.

The manufacturer assumes no liability for print errors.

### 1.1.4 Requirements for safe use

Note the following points for safe use:

- If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.
- The product may only be used according to the authorized use specified above.
- The product may only be supplied with power by the energy sources mentioned in this operating manual.
- The product may only be used under the environmental conditions mentioned in this operating manual.
- The product must not be opened or modified.

The product must not be used if:

- it is visibly damaged (e.g. after being transported)
- it was stored under adverse conditions for a lengthy period of time (storage conditions, see chapter "Specifications")

### 1.1.5 Keep instructions

The manual must be kept in the vicinity of the product so you can always find the information you need.

### 1.1.6 User qualification

The operating personnel must be able to understand and correctly implement the safety labels and safety instructions on the packages and inserts of the products.

The user must be capable and able to read and understand this manual in order to familiarize themselves with the handling and to ensure safe use.

### 1.1.7 Handling of hazardous chemicals

Chemical and/or biological hazards may exist where this product is used. Abide by all governing laws, regulations and protocols when using this product.

For the development of products, Lovibond® pays close attention to safety. Some hazards from dangerous substances cannot be avoided. If self-produced tests or solutions are used, the responsibility concerning any risks caused by those tests or solutions lies with the user (personal responsibility).

EN

### 1.1.8 Disposal notes

Dispose of the batteries and electrical devices at a suitable facility according to local legal requirements.

It is illegal to dispose of the batteries with household waste.

Within the European Union, the batteries are removed at a specialized treatment center at the instrument's end of life.

Instruments marked with this symbol must not be disposed of in normal domestic waste.

## 1.2 List of all used signs in the document

The following symbols are used in this manual to identify sections that require special attention:

**Danger!**

A hazard exists that will result in death or severe injury if not avoided.

**Warning!**

Improper handling of certain reagents can cause damage to your health. In any case follow the safety labels on the packing, the safety instructions of the package insert and available SDS. Protective measures specified there have to be followed exactly.

**Caution!**

A hazard exists that may result in minor or moderate injury.

**Notice!**

Important information or specific instructions need to be strictly followed.

## 2 Product overview

### 2.1 Guide to symbols

Labels attached to the product should be strictly observed to avoid personal injury or damage to the product. Refer to this chapter for information regarding the nature of the danger or risk before taking any action where such label is present.

**For professional users in the European Union:**

If you wish to discard electrical and electronic equipment (EEE), please contact your dealer or supplier for further information.

**For disposal in countries outside of the European Union:**

This symbol is only valid in the European Union (EU). If you wish to discard this product please contact your local authorities or dealer and ask for the correct method of disposal.

EN

### 2.2 Certification

Device/Accessories	Test approvals, EMC
BOD Measurement System	CE, EMC in accordance DIN EN 61326, basic requirements
Stirrer unit	CE, DIN EN 61326:2013
Mains adapter	CE, EN 55024

### 2.3 Features

The BD 600 is the ideal system for process control or as a secondary test to the dilution method. Using a respirometric system solves many of the issues associated with the dilution method for BOD analysis.

- Significantly reduced sample preparation - it is easy to begin collecting data - simply add sample to the bottle, add nitrification inhibitor and place the sensor head on the bottle. There is no need for sample dilutions or seeding the sample. Because this is a respirometric test, concerns regarding air in the sample are eliminated.
- Operators do not have to be present at completion of tests. Measurements are automatically taken and stored at regular time intervals and can be recalled at your convenience - no more weekend trips to the lab or plant!
- Easy and quick interpretation of measurement values - In addition to displaying results of a sample directly in mg/l BOD, the large, integrated, display is able to show the BOD curve which makes it is easy to see and understand measurement data points and trends. In addition, each sensor head has an integrated LED. This allows operators to visually confirm which sample the data is being displayed for and reduces potential confusion or errors over which sample data is presented.
- Consistent measurement conditions - The automatic start function ensures all measurements are taken at an equal sample temperature. When the measurement conditions are consistent across readings, operators do not have to work about temperature or pressure fluctuations skewing results.

- Sample identification and tracking - The BD 600 interface, makes it easy to assign an identification code or sample name when setting up a new sample. This helps track a sample from the time it is collected and ensures traceability during sample measurement.
- Use the remote control - If using an incubator with a glass door, operators can use the supplied remote control to display measurement data for a sample - without even opening the door! This further ensures that samples stay at a consistent temperature which is essential to accurate results.

## 2.4 Product description

Biochemical oxygen demand (BOD) is an expression for the quantity of oxygen required for biological degradation of organic matter in a waste water sample. The BD 600 measuring system records a measurement once every hour, independent of the length of the measuring period. This way the quality of the series of measurement can be evaluated at an early stage. Current values and stored values may be called up at any time. Stored values can be displayed numerically or graphically.


### 2.4.1 Description of the keys and switches

Key	Function
	<b>ON/OFF key</b> This key is used for switching the device on and off.
	<b>Function keys</b> The function keys have a different meaning in each menu. Info text above the keys indicates their meaning. If no text is shown above the key, it has no function.
	<b>Backspace</b> Characters which have been entered can be deleted with the backspace key.
	<b>Number keys</b> The number keys are used for entering the date, time, as well as file and head names.
	<b>Quick selection keys</b> The quick selection keys take you directly to the corresponding menu each submenu: Quick selection key 1 ,Start' --> ,Start measurement series' Quick selection key 2 ,List' --> ,Display current values' Quick selection key 3 ,Graph' --> ,Display measurement series'
	
	

Key	Function
	<b>Arrow keys</b> They are provided mainly for navigation through the menus.
	<b>Head keys</b> Individual measurement points and/or heads can be selected in the submenus with the head keys.

EN

## 2.4.2 Interface description

1. SD card holder
2. USB host interface

### NOTICE!

The USB host interface is only intended for USB sticks. USB hubs, external hard disks, and adapter sticks (e.g. USB SD card adapters) are not supported. Battery operation entails the limitation that the host interface can only provide 200 mA.

### 3. USB device interface

The USB device interface is located under the USB host interface on the side of the front section of the housing. It makes it possible to access data saved on the SD card with a PC.

For this purpose, connect the device to a PC via the appropriate USB cable. The device must be switched on. If necessary, deactivate the auto-off option. If there is an SD card in the card holder, the device appears as mass storage on the PC.

4. Connection jack for the mains adapter.
5. Window for the IR receiver for the remote control.

## 3 Commissioning

### 3.1 Operating environment

Temperature equalisation is essential prior to biological testing and reproducible BOD measurements, as temperature has a major effect on biological activity. BOD measurements, for example, are always performed in a thermostatically controlled cabinet at a temperature of 20 °C. For temperature equalisation, we recommend Lovibond® thermostatically controlled cabinets with a user-selectable temperature from 2 °C to 40 °C.

EN

### 3.2 Contents of delivery

#### CAUTION!

Inspect items to ensure no damage has occurred during shipment. If there is damage or something is missing, please contact local distributor immediately.


#### Scope of delivery

BD 600 or BD 600 GLP	BD 606	
1	2	complete unit with 6 sensors and control unit with batteries (BD 600 GLP with certificate)
1	2	power supply unit incl. Y-cable for common power supply of instrument and stirring unit
1	1	remote control (without batteries)
1	2	inductive stirring unit
6	12	sample bottles
6	12	rubber gaskets
6	12	magnetic stirring rods
1	1	overflow flask, 157 ml
1	1	overflow flask, 428 ml
1	1	bottle, 50 ml potassium hydroxide solution
1	1	bottle, 50 ml Nitrification inhibitor
1	1	instruction manual
1*	2	USB cable

\* BD 600 only

### 3.3 Inserting and replacing the batteries

The device can be supplied with power via the accompanying mains adapter or batteries. If the mains adapter is connected and batteries are inserted, the device is supplied via the mains adapter and the batteries are not used. If the mains adapter is removed during operation, a seamless transition to battery operation takes place automatically.

Remove all bottles from the rack in order to insert the batteries. Open the battery compartment (**A**) and insert three size C alkali-manganese cells (LR14 / ,baby cell') in the tubes provided for this purpose. The tube makes it easier to insert the batteries and prevent the batteries from jumping out again when you attempt to close the compartment.

Otherwise, if the battery polarity is reversed it can be damaged and leak and damage the device. Then insert the batteries with tubes into the battery compartment (ensure correct polarity).

**WARNING!**

Insert the batteries in the tube such that a positive pole always faces a negative pole.

**WARNING!**

The device is not designed for rechargeable batteries. Therefore, they must not be used. Rechargeable batteries can be damaged, leak, and damage the device.

**CAUTION!**

Do not handle depleted batteries or leaking material with bare hands! Wear protective gloves! Avoid contact with the eyes and skin!

**DANGER!**

Only handle the mains adapter with dry hands!

Protect the mains adapter from moisture.

Do NOT open the mains adapter!

Do not subject the mains adapter to great force!

Do NOT use the mains adapter if there is damage to the housing or the connection contacts!

Only use the mains adapter with an appropriate socket intended for this purpose!

Ensure that this socket is in fault-free condition!

The mains voltage and the mains frequency must always lie within the specified range of the mains adapter.

---


## 4 Operation

### 4.1 First start up

1. Place the instrument rack on the inductive stirring plate.
2. Use the Allen key to adjust the height between the rack and the inductive stirring plate. Or adjust the height later if continuous stirring is not achieved.
3. Use the Y-cable to connect the instrument and the inductive stirring plate to the power supply. Also, insert the batteries to prevent the unit from turning off in the event of a power failure.
4. Press the power button to start.
5. Select the desired language using the arrow keys and confirm with F2.
6. Set Date/Time using the arrow keys and confirm with F2.
7. The Lovibond logo appears briefly and will turn to the main menu. The instrument is now ready.

### 4.2 General operation principles

#### Measuring Principle

#### Biochemical Oxygen Demand (BOD)

The biochemical oxygen demand (BOD) in water (e.g. waste water, surface water) is the amount of oxygen that is consumed during the degradation of organic substances through biochemical processes.

#### Measuring principle

Respirometric methods provide direct measurements of the oxygen consumed by microorganisms in the waste water from an air or an oxygen-enriched environment in a closed vessel under conditions of constant temperature and agitation. Carbon dioxide produced metabolically by the bacteria is chemically bound by the potassium hydroxide solution contained in the sealed cup in the bottle. The result is a pressure drop in the system, which is directly proportional to the BOD value and is measured by the pressure sensor. The BOD level is then displayed directly in mg/l.

**General Operation Procedure (Brief Summary)**


- Estimate the measurement range of the sample to be tested and use the sample volume as indicated in section "Handling".
- If necessary, pre-treat the sample as indicated in section "Measuring Principle" (e.g. adjust pH value and temperature, filtration, etc.).
- Measure the sample volume precisely with the overflow volumetric flask and pour into the BOD bottle (use a funnel, if necessary).
- If necessary, add the nitrification inhibitor as indicated in section "Handling".
- Place the magnetic stir bar in the BOD bottle.
- Fill the seal cup with 3 - 4 drops of KOH solution and place the seal cup in the test bottle.
- Screw the BOD sensors on the test bottles.
- Place the sample bottle in the bottle rack.
- Start the test (see section "Handling").
- Incubate the sample according to specifications (e.g. BOD<sub>5</sub> at 20 °C).

## 4.3 List of control elements and their function

### Display overview


Description	Function
<b>1 Header</b>	Displays the name of the menu (left) and date/time (right)
<b>2 Display Area</b>	Area that displays application steps of the corresponding menu
<b>3 Footer</b>	Displays the meaning of function keys F1 & F2 (left & right). If nothing is displayed, the corresponding keys have no function. Displays the icon for data storage (SD/USB) and power supply (battery/plug)

### Icon overview

Icon	Description	Function
	SD Card*	Displays that SD card is recognized by the system.
	USB drive*	Displays that USB drive is recognized by the system.
	Battery power	Displays that system is powered by batteries and indicates battery life.
	Plug symbol	Displays that system is powered by mains adapter.

\*If the icons are faded, the components are not recognized or connected by the system.

## 4.4 Main menu


Use the up ▲ and down ▼ arrow keys to select a submenu.

Use the right ► arrow key to access a submenu. Alternatively, use the quick selection keys "Start", "List" and "Graphic" to access the first three submenus.

EN

### 4.5.1 Start test series

After sample preparation, enter the submenu "Start test series" to begin a measurement.


Use the up ▲ and down ▼ arrow keys to select a setting:

- Bottle position - Name
- Range - Sample volume / ATH
- Measurement duration / Time interval

Function	Description
<b>Bottle Position – Name</b>	When a sensor head is detected by the system, the bottle position and name are displayed on the screen. Use the left ◀ and right ▶ arrow keys or the head keys to select one of six possible sample bottles. The bottle selection can be tracked by a red LED on the top of the sensor head. A name for each bottle can be set in the "Options" submenu under "Name heads", otherwise no name will appear on the screen.
<b>Range – Sample Volume/ATH</b>	Select the expected BOD measuring range of your sample. The required sample volume and the number of drops for the nitrification inhibitor (ATH) are displayed. An overview of all BOD measuring ranges, required sample volumes and number of ATH drops is summarized under section "Handling".
<b>Measurement duration – Time interval</b>	Specify the measuring duration: 1 – 28 days. Standard storage interval: <ul style="list-style-type: none"> <li>• 1st day – measurement every hour</li> <li>• 2nd day – measurement every 2 hours</li> <li>• 3rd to 28th day – measurement every 24 hours</li> </ul>
<b>A</b>	The A-symbol indicates if the Autostart mode is activated in the "Options" submenu. If the background is grey, the Autostart mode is deactivated.
<b>Back</b>	Press function key F1 to return to the main menu.

Function	Description
<b>Start</b>	Press the function key F2 to start the test series. If data are already stored in the sensor head, a reminder message "The current test series is being overwritten!" appears on the screen. Press the function key F2 again to overwrite the current data and start the new test series.

#### 4.5.2 Show current values


The latest measurements can be tracked in the submenu "Show current value". All six bottle positions are displayed in the center of the screen. The latest value is displayed next to the sensor head if connected to the system.

EN

Function	Description
	Sensor head with a completed test series is recognized.
	Sensor head with an on-going test series is recognized.
	Sensor head is not connected or recognized.
<b>Back</b>	Press function key F1 to return to the main menu.
<b>Update</b>	Press function key F2 to update the current value. Alternatively, press the head keys to update the measurement value of a specific bottle position.

#### 4.5.3 Display test series


A schematic diagram of all bottle positions is displayed on the right side of the screen. The head keys can be used to select a bottle position. After, the measured data of the selected head are listed in mg/l with time stamp.

Use the up ▲ and down ▼ arrow keys to scroll through the measurements.


Use the left ◀ and right ▶ arrow keys to show additional information such as measurement parameters and head name.

Function	Description
	Sensor head is recognized and can be selected.
	Sensor head is recognized and selected.
	Sensor head is not connected or recognized.
x / y / z	<p>Example: "2 / 3 / 39"</p> <p>x: Selected data point in the test series; example "2": second measurement selected in the table.</p> <p>y: Number of data points already recorded during a test series; example "3": three measurements have already been recorded.</p> <p>z: Number of total measurements within the test series; example "39": thirty-nine measurements to be done</p>
Back	Press function key F1 to return to the main menu.
Graphics	Press the function key F2 to plot the measured data in a diagram (see below).

Use the head keys to select and deselect a sensor head to plot the measured data. A maximum of three plots can be displayed simultaneously for comparison.

## Possible Errors

**Underr.** Underrange appears instead of a measurement value, if the recorded measurement value is below the starting measurement value.

**Overr.** Overrange appears instead of a measurement value, if he recorded measurement value is outside the measuring range.


**Fail** Fail appears instead of a measurement value, if no measurement could be conducted for the provided time; e.g. the sample was removed before the end of the measurement.

The eye symbol indicates that no usable measurement value was determined in at least one measurement within the test series.


The warning symbol indicates that a test series has been ended, although not all measurements have been performed.

**(A)** Underrange is shown with zero.


**(B)** Overrange is shown with the maximum display value within the measuring range.

## 4.5.4 Export test series

A schematic diagram of all bottle positions is displayed on the left site of the screen. Use the head keys to select and deselect the sensor heads for data export. Only selected heads are considered for data export. Use the up ▲ and down ▼ arrow keys to select between SD-Card and USB-drive if both are connected. Use the right ▶ arrow keys to continue with the data export.

Function	Description
	Sensor head is recognized and can be selected.
	Sensor head is recognized and selected.
	Sensor head is not connected or recognized.
<b>M: (SD)</b>	Appears when SD-Card is connected to the instrument.
<b>U: (USB)</b>	Appears when USB-drive is connected to the instrument.
<b>Back</b>	Press function key F1 to return to the main menu.


After selecting a storage medium, existing files and folders are displayed. Use the up ▲ and down ▼ arrow keys to select a folder or a file.

Function	Description
 folder ▶	When a folder is selected, use the left ◀ and right ▶ arrow keys to enter the folder or return to the previous entry.
 bod 05.01.2021.csv ▶	When an existing file is selected, use the right ▶ arrow key to enter the status bar. The existing file can now be renamed using the number keys, the backspace key, the left ◀ and right ▶ arrow keys. Press F2 to save and overwrite the existing file. Use the up ▲ and down ▼ arrow keys to exit the status bar.
m:/bod 05.01.2021(1).csv	The destination folder is always displayed at the bottom of the list together with an automatically generated file name. The file name can be modified. Therefore, use the up ▲ and down ▼ arrow keys at the beginning or the end of the list.

Function	Description
	Folder and files names are displayed in lower-case letters. Folder and file names with special characters are not displayed, but are indicated by the warning symbol at the top right of the screen.
Back	Press function key F1 to return to the previous menu.
Save	Press function key F2 to save the measured data in csv-format onto the storage drive. Only data of selected sensor heads are saved.

### 4.5.5 Options

Instrument settings can be changed in the “Options” submenu. Use the up ▲ and down ▼ arrow keys to select a setting. Use the left ◀ and right ▶ arrow keys to change the selected setting. Use function key F2 to accept the new settings or function key F1 to return to the main menu.

Function	Description
Auto start	The auto-start function is used to equalise the sample temperature with the ambient temperature. Ideally, the measuring system is operated in an incubator with a constant temperature of 20 °C. If the sample temperature is cooler than the ambient temperature, a pressure increase is measured until the sample has reached the ambient temperature. The temperature adjustment is completed when the sensor head detects a first pressure drop. The first pressure drop is the starting value for the BOD measurement. If no pressure drop is detected, the measurement starts automatically after a maximum waiting time of 3 hours. Using the auto-start function does not replace proper temperature control of the sample (15 - 21 °C) before starting the measurement. BOD samples whose temperature deviates strongly from the target temperature lead to incorrect BOD values due to strong pressure differences.
Brightness	Adjust the background lighting of the screen in 10 % increments from 0 % to 100 %.
Contrast	Adjust the contrast of the screen in 10 % increments from 0 % to 100 %.
Name heads	Name the sensor heads
Language	Change the display language
Date/Time	Change date and time

Function	Description
Date (format)	Change the representation of the date in the header. The following options are available: <ul style="list-style-type: none"> <li>• DD.MM.YYYY / 31.12.2021</li> <li>• MM.DD.YYYY / 12.31.2021</li> <li>• YYYY.MM.DD / 2021.12.31</li> </ul>
Time (format)	Change the representation of the time in the header. The following options are available: <ul style="list-style-type: none"> <li>• 12 h</li> <li>• 24 h</li> </ul>
Auto-Off	Use the auto-off function to save power. Select between 3, 5 or 10 minutes for the unit to switch off automatically at the selected time. When the Auto-Off is deactivated, the unit runs continuously.
Remote Control	Select "on" or "mains supply" to activate the IR interface and enable the use of the remote control. If "mains supply" is selected, the IR interface is only active when the instrument is powered by DC.
Device ID	Assign a Device ID to operate the unit using the remote control
Update	Enter update to see the latest version of the firmware or to update the firmware of the instrument. Firmware Update: Transfer the ‚bod.hex‘ update file to the root directory of an USB drive or a SD card and insert the storage medium to the unit. If both USB drive and SD card are inserted, the file on the SD card is used. A software update is only possible in mains operation. Press function key F2 to start the update. After completion, the instrument restarts. Firmware updates can be found on our website <a href="http://www.lovibond.com">www.lovibond.com</a> .


**NOTICE!**

To prevent loss of stored test results store or print them out before performing an Update. If the update procedure is interrupted (eg. interruption of connection, LoBat., etc.) the instrument isn't able to work (no display). The instrument will only work again after completing the data transfer.

## 4.6 Advanced operations


### Name Heads

Enter the setting "Name heads" under "Options" to give the sensor heads a name. All six bottle positions are displayed in the center of the screen. The head keys can be used to select a bottle position. The bottle selection can be tracked by a red LED on the top of the sensor head. Use the number keys, the backspace key and the left ◀ and right ▶ arrow keys to enter a name for a selected sensor head.

## Remote control

### Purpose

With the remote control, the submenus "Show current values", "Display test series" and "Export test series" can be accessed during operation or when the unit is switched off. This enables remote monitoring of a test series, e.g. when the BD 600 measuring system is placed in an incubator for temperature control.

### Initial Set-Up

1. The remote control is delivered without a battery. Before use, insert a lithium button cell (e.g. CR2025).
2. Activate the remote control in the "Options" menu.
3. Assign a Device ID to the instrument in the "Options" menu. If multiple instruments are used, assign different Device IDs.

Device ID 13:00  
05.01.2021

01

SD USB F5

### How to use the remote control

1. Point the remote control at the instrument. Press the On key to switch on the instrument. The Device ID appears on the screen. (If other instruments are nearby and appropriately configured, their screen will also display the Device ID)
2. Use the number keys to enter the Device ID of the desired instrument. (If no or an incorrect ID is entered, the instrument will return to its previous status)
3. The submenu "Show test series" will appear on the screen.
4. The remote control can now be used to operate the unit in the same way as the keypad.

Press the Off key to switch off the instrument or press the On key to access the Device

ID screen in order to continue with another instrument.

#### 4.7.1 Notice

##### Notice!

The safety datasheets of the chemicals comprise all instructions on safe handling, occurring hazards, preventive actions and actions to take in hazardous situations.

EN

#### 4.7.2 Initial Steps

Place the BOD system, consisting of rack and inductive stirring plate, in an incubator for temperature control. Set the temperature of the incubator to  $20 \pm 1$  °C as recommended for the determination of BOD according to EN 1899. Connect the system to the mains supply.

#### 4.7.3 Selection of the sample volume

The expected BOD value of the sample defines the necessary sample volume for the experiment. Select a measuring range that results in a BOD value in the upper half of the measuring range, e.g. for a BOD of 150 mg/L select the measuring range 0 - 200 mg/L. If the expected BOD value is unknown, an estimate can be made on the basis of the COD value (COD = chemical oxygen demand), e.g. for domestic waste water the  $BOD_5$  is approx. 80 % of the COD value.

BOD range in mg/L	Sample volume in mL	Nitrification inhibitor ATH dosage
0 – 40	428	10 drops
0 – 80	360	10 drops
0 – 200	244	5 drops
0 – 400	157	5 drops
0 – 800	94	3 drops
0 – 2000	56	3 drops
0 – 4000	21.7	1 drops

#### 4.7.4 Sample Preparation

1. Use clean equipment for sampling.
2. Ensure the use of representative samples.
3. Test the pH value of the sample. The ideal pH value must be in the physiological range of the microorganisms, which is usually between pH 6.5 and 7.5. Any greater deviations provide a lower BOD value. If the sample pH is too high, use diluted hydrochloric acid (1 M) or diluted sulphuric acid (1 M) to adjust the pH. If the pH value is too low, use a sodium hydroxide solution (1 M) to adjust the pH.

4. Prepare the sample in accordance with the respective normative requirements so that it can be used as settled, stirred, filtered or even homogenized with a blender. It is recommended to test each sample two to three times.
5. The sample must be brought to the desired measuring temperature  $\pm 1$  °C.
6. Add a clean magnetic stir bar to the BOD sample bottle for continuous agitation.
7. Measure the required sample volume depending on the expected measuring range (see "4.7.3 Selection of the sample volume") with the corresponding overflow flask and fill the exact sample volume into the BOD sample bottle. Use a funnel if necessary.
8. Depending on the sample volume, add the recommended number of drops of nitrification inhibitor (ATH) to the BOD sample bottle (see "4.7.3 Selection of the sample volume").
9. Fill the dry rubber gasket with 3 - 4 drops of 45 % potassium hydroxide solution. Then place the rubber gasket on the BOD sample bottle.
10. Screw on the sensor head hand-tight. The rubber gasket ensures the necessary sealing. (Do not use grease or other lubricants!).
11. Place the prepared BOD sample bottle into the pre-tempered rack inside the incubator at  $20 \pm 1$  °C (see "4.7.2 Initial Steps"). (Note: The BD 600 has an optional autostart function which allows the use of samples with a temperature of 15 to 21 °C. With the auto start function switched on, the system checks within the first three hours if a pressure drop is recognized by the sensor head. If so, the measurement will start immediately or after the three hours.)
12. Check if the sensor head is recognized by the system and if the sample is continuously stirred. If necessary, adjust the 4 ball screws on the rack.
13. Start the test series (see "4.5.1 Start test series")
14. Incubate the sample in accordance with the respective normative requirements (e.g. BOD<sub>5</sub> at 20 °C for 5 days).

#### **Additional Information:**

#### **Sample water**

- Domestic wastewater usually contains sufficient nutrient salts and suitable microorganisms and no toxic or retarding substances. Therefore, the samples can be used undiluted without the need for additional nutrients or seeding by microorganisms.
- Industrial wastewater may contain inappropriate microorganisms, toxic and inhibitory substances as well as insufficient nutrient content. Such samples must be diluted with dilution water containing a sufficient amount of nutrients to achieve an adequate BOD:N:P ratio and to eliminate the inhibitory effect of toxic substances. Furthermore, an additional seeding by a microbial source (inoculum) is often needed. An application paper on strongly-loaded organic waste water can be downloaded on [www.lovibond.com](http://www.lovibond.com).

#### **Potassium Hydroxide**

When oxygen is consumed by the microorganisms to form carbon dioxide during the degradation of organic compounds, there is no direct change in pressure. Potassium hydroxide in the rubber gasket and carbon dioxide react chemically to form potassium carbonate:  $2 \text{KOH} + \text{CO}_2 \rightarrow \text{K}_2\text{CO}_3 + \text{H}_2\text{O}$


This removes the carbon dioxide formed from the gas phase creating a measurable negative pressure which correlates with the oxygen consumption of the microorganism and ultimately with the measured BOD.

### Nitrification Inhibitor

Nitrifying bacteria also consume oxygen. This consumption can already occur within the first five days, especially in samples with low BOD values. In wastewater monitoring, often only the BOD from the degradation of carbon compounds is determined (carbon BOD), with the addition of allylthiourea (ATH) nitrification is inhibited to prevent oxygen consumption from the degradation of nitrogen compounds. If the oxygen consumption in the course of the nitrification (N-BOD) should be determined, a comparison of the sample with and without nitrification inhibitor can be done. The difference between both BOD values will correspond to the oxygen requirement of nitrifying bacteria.

## 4.8 Description of setup and handling of important accessories

### Inductive stirring system

1. Mains adapter
2. DC connection cable
3. Stirrer drive
4. Stirring position
5. Y-cable

### Device description and functional description

The inductive stirring system is designed for the stirring of liquids in special BOD bottles. It comprises a super-flat stirrer drive with 6 stirring positions and the mains adapter. The inductive stirrer drive does not have a motor and is thus wear-free. It is especially well-suited for stirring during continuous operation in thermocabinets. It is completely protected from dripping water due to its hermetically sealed encapsulation. It can be used in harsh environmental conditions.

The wide recess between the stirring positions also ensures greater ventilation of the containers in the thermocabinet. The stirrer drive is supplied with power via the mains

adapter and has integrated control electronics. When switched on, a soft-start at reduced speed assures an even start-up and centring of the magnetic stir bars. The electronically controlled automatic monitoring unit reduces the speed about every 40 seconds. In the process, each magnetic stir bar is re-centred in the bottle for a few seconds. Therefore, you can easily exchange individual bottles while the unit is switched on.

Due to the synchronous operation, mutual interference of the magnetic stir bars can be practically ruled out altogether.

## Magnetic stir bars

Use the PTFE covered magnetic stir bars included in the scope of delivery.

### **NOTICE!**

The length of the magnetic stir bars should not exceed 40 mm.

---

### **NOTICE!**

#### **Stirring**

Do not place hot containers on the stirrer drive; maximum temperature: 56 °C.

Fill the BOD bottles.

Place a magnetic stir bar in each bottle.

Place the BOD bottles in the rack.

---

### **Warning due to magnetism!**

The magnetic fields can influence parts that are sensitive to magnetic fields, magnetic parts, or metal parts (e.g. data carrying media, pacemakers, watches, etc.). Keep these parts away from the stirrer drive and magnetic stir bars.

---

## 5 Maintenance

### 5.1 Calibration

A test kit (Art. no. 2418328) is available for testing the BD 600.

The test kit enables the testing of all components. It comprises special reagent tablets which generate a defined vacuum in the closed BOD bottle.

### 5.2 Decommissioning

Disconnect the external power supply from the unit in order to decommission the system. Disconnect the mains adapter plug from the mains supply. Remove batteries from the measurement system. Remove seal caps from the test bottles and clean in the correct manner. Empty and clean test bottles correctly. Clean sensor heads. Store the stirrer drive and magnetic stir bars in a manner such that no hazards arise for parts which are sensitive to magnetic fields.

**CAUTION!**

Observe the SAFETY INSTRUCTIONS at the beginning of the manual.  
Observe the local legal regulations for all emptying and cleaning.

### 5.3 Maintenance and cleaning

- High-quality metal alloys are used for the contacts between the BOD sensors and the bottle rack. Carefully clean the contacts as necessary with a soft cloth. To smooth out any unevenness, for optimisation of the contact between the BOD sensor and bottle rack, and for optimisation of the stirrer position, the BOD base unit has 4 adjusting screws on the bottom. Carefully clean the BOD base unit (including bottle rack) and the BOD sensors as necessary with a dry cloth. The parts coming into contact with the samples (BOD bottle, seal cup, magnetic stir bar) must be carefully cleaned after each test. Empty the bottles after the testing is complete – observe local regulations in the process – and rinse out repeatedly with hot water. Rinse thoroughly after using cleansers! Residue from cleansers can destroy the BOD test.
- The stirrer drive is maintenance-free. The magnetic drive is installed inside the housing to ensure that it is waterproof. Clean the surface of the stirrer drive regularly. The stirrer drive can be cleaned with cleansers or di-sinfecting solutions which are suitable for PVC. Wipe off the surface of the mains adapter with a dry cloth.

**CAUTION!**

The device may only be opened by an authorised service location in the case of a repair. Disconnect the device from the mains before opening!

## 6 Troubleshooting

### 6.1 Error Messages and Notices


Message in the display	Meaning
RTC Error! Please set the date and time again.	The message appears when the device is switched on and the RTC (real-time clock) has forgotten its time. Adjust the time in the following menu and then check all pending measurements.
LOBAT!	This appears just before the device switches off on its own. Replace the batteries.
USB-Host overcurrent! Remove the USB device.	This message appears if a short-circuit has occurred at the USB host port. Remove the USB device.
Incorrect head ID! The measurement series was not started!	This appears at the start of a new measurement series when a different head was found at the selected measurement position than for the reading of measurement parameters.
Communication error! The measurement series was not started!	This appears when a new measurement series is started, if there is no head available at the selected measurement position.
The battery level is low!	This appears when a new measurement series is started, if the batteries are nearly depleted. Replace the batteries or connect the mains adapter. The measurement series can still be started.
The current measurement series will be overwritten!	Appears as a warning when a new measurement series is started. The current measurement series will be overwritten!
Internal processing error! The measurement series was not started!	This appears when a new measurement series is started, if the displayed measurement parameters cannot be correctly interpreted. Change the measurement parameters back to the previous parameters. Then restart measurement series.
Start error! Please restart the measurement series!	This appears when a new measurement series is started, if the measurement could not be started. Check the contacts at the head and the measurement position. Restart the measurement series.
The measurement series was started.	Note. This appears if a measurement series was started successfully.
Updates are only possible in mains operation.	This appears if an update is started from the Options menu and the mains adapter is not connected.
The update file ,bod.hex' was not found.	This appears if an update was started from the Options menu and the file ,bod.hex' cannot be found on the SD card or the USB stick.
Do you really want to start the update?	Security query. This appears if an update is started from the Options menu and no errors have occurred.

<b>Message in the display</b>	<b>Meaning</b>
No head was selected!	This appears on data export, if the export is started and no head was selected.
The medium is write-protected!	This appears on data export to the SD card, if it is write-protected. Remove the write protection.
The path or file name is invalid!	This appears on data export if the path or file name is invalid. Enter a new path and/or file name.
The file already exists. Would you like to replace it?	This appears on data export when an attempt is made to overwrite an existing file.
The file was saved.	Note. This appears after a successful data export.
Communication error! Error when saving!	This appears after an unsuccessful data export. Carry out the export again. Ensure that there is sufficient memory available on the medium and that the device may write to the folder.
Incorrect head ID!	This appears when saving the head name, if a different head was found at the selected measurement position than for the reading of measurement parameters.
Communication error! The name was not saved!	This appears when saving the head name, if there is no head available at the selected measurement position.
The name was saved.	Note. This appears if the head name was saved successfully.
Adjusting the time while measurement series are in progress can falsify the correct measurement times!	This appears as a warning if you attempt to adjust the date and time in the options, because errors can occur when measurement series are in progress.
Communication error! Error when saving the options!	This appears when options are adopted, if they could not be saved successfully. Switch off the device, wait at least 10 seconds, and then switch on the device again. Repeat the settings and attempt to adopt the options. If the message appears again, there is a defect in the device. In this case, contact customer service.

## 6.2 Interpretation of BOD curve and Problem Solving Strategies

Example: Measuring Range: 0 – 40 mg/L

Graph	Description	Action
A	Ideal BOD curve	
B	Additional BOD due to nitrification	Add nitrification inhibitor
C	The BOD value is above the measuring range	Increase measuring range or dilute sample
D	Measuring system is leaking	Check BOD bottle for leakage and replace accessories (e.g. rubber gasket) if necessary
E	<ul style="list-style-type: none"> <li>• Not sufficient microorganisms</li> <li>• Not temperature adjusted</li> </ul>	<ul style="list-style-type: none"> <li>• Inoculate sample</li> <li>• Adjust temperature</li> </ul>

## 7 Accessories & Replacement parts

### 7.1 List of Accessories

Title	Part Number
6 BOD Sample flasks, brown glass	418645
BOD Sample flask, brown glass, 500 ml	418644
BOD sensor	2444470
Complete Set Overflow measuring flasks	418654
Inductive Stirring Unit	2444456
Magnetic stirring rod	418633
Magnetic stirring rod remover	418638
Nitrification inhibitor	2418642
Overflow measuring flask 21.7 ml	418664
Overflow measuring flask 56 ml	418655
Overflow measuring flask 94 ml	418656
Overflow measuring flask 157 ml	418657
Overflow measuring flask 244 ml	418658
Overflow measuring flask 360 ml	418659
Overflow measuring flask 428 ml	418660
Potassium hydroxide solution 45%	2418634
Power supply unit for inductive stirring unit	444454
Remote control	2444481
Rubber gasket, 4.5 cm	418636
Test set for checking the system (10 tablets)	2418328
USB cable 3 m	2444482
Y-cable	2444475

EN

## 8.1 Specifications - Mains adapter

Type	SRB1502300P
Design	Mains adapter plug, switching power supply
Input voltage, frequency	100 – 240 V ± 10 %, 50 / 60Hz
Input current	1000 mA
Primary adapter	Europe, UK, Australia, USA
Protection class	II
Output voltage, frequency	15 V, DC
Max. output current	2300 mA
Output protection	Short-circuit-proof
DC cable length approx.	1800 mm
Environmental conditions	0 °C to 40 °C, with max. 93 % relative air humidity, non-condensing, maximum altitude 2000 m, pollution degree I
Energy efficiency class	VI
Safety standard	EN 60950, EN 62368-1
EMC	IEC 61204-3
Dimensions [mm]	62 x 84 x 53
Weight, with EU adapter	258 g
Approval, EMC	CE, EN 55024

### Caution!

Subject to technical modification!  
To ensure maximum accuracy of test results, always use the reagent systems supplied by the instrument manufacturer.

## 8.2 Specifications - Stirrer unit

Type	Inductive stirring system
Number of positions	6
Stirring output	7W
Speed	320 rpm, brief centring phase every 40 s
Dimensions (W x D x H) (mm)	270 x 180 x 25
Stirring positions spacing	88 mm
Weight (stirrer drive)	1204 g

Type	Inductive stirring system
Housing material	PVC
Environmental conditions	-10°C to 56°C, 95 % relative air humidity, non-condensing, maximum altitude 2000m, pollution degree I
Test approvals, EMC	CE, DIN EN 61326

**Caution!**

Subject to technical modification!  
To ensure maximum accuracy of test results, always use the reagent systems supplied by the instrument manufacturer.

### 8.3 Specifications - BD 600

Measuring Principle	Manometric; mercury-free; electronic pressure sensor
Measuring Range	0-40, 0-80, 0-200, 0-400, 0-800, 0-2000, 0-4000 mg/l
Display	Backlit
Interfaces	<ul style="list-style-type: none"> <li>• SD-Card</li> <li>• USB</li> </ul>
Auto – OFF	No
External Storage	USB / SD-Card / Cartão USB / SDUSB / carte SDUSB / tarjeta SDUSB / scheda SDUSB / SD-kaar
Measurement Time	User-selectable, between 1 and 28 days
Power Supply	<ul style="list-style-type: none"> <li>• 3 alkali-manganese batteries (Baby cells/ size C)</li> <li>• 100 - 240 V / 50-60 Hz</li> </ul>
Clock	Real Time Clock and Date
Portability	Benchtop
Autostart	Yes
Measuring Stations	6
Storage Interval	- hourly (1st day) - every 2 hours (2nd day) - 1x daily (3rd - 28th day)
Protection Class	IP 53
Compliance	<ul style="list-style-type: none"> <li>• CE</li> <li>• EMC according to DIN EN 61326</li> </ul>

Dimensions	181 x 230 x 375 mm
Weight	4.1 kg (4100 g, unit with bottles & batteries 5775 g, complete with stirring unit)
Weight with Packaging	(4100 g, unit with bottles & batteries 5775 g, complete with stirring unit)

**Caution!**

Subject to technical modification!  
To ensure maximum accuracy of test results, always use the reagent systems supplied by the instrument manufacturer.

## 9 Appendix

### 9.1 List of Trademarks

Lovibond® and Tintometer® are registered trademarks of the Tintometer group of companies. All translations and transliterations of Lovibond® and Tintometer® are asserted as trademarks of The Tintometer® Group.

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- Consistent measurement conditions - The automatic start function ensures all measurements are taken at an equal sample temperature. When the measurement conditions are consistent across readings, operators do not have to work about temperature or pressure fluctuations skewing results.
- Sample identification and tracking - The BD 600 interface, makes it easy to assign an identification code or sample name when setting up a new sample. This helps track a sample from the time it is collected and ensures traceability during sample measurement.
- Use the remote control - If using an incubator with a glass door, operators can use the supplied remote control to display measurement data for a sample - without even opening the door! This further ensures that samples stay at a consistent temperature which is essential to accurate results.

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The wide recess between the stirring positions also ensures greater ventilation of the containers in the thermocabinet. The stirrer drive is supplied with power via the mains adapter and has integrated control electronics. When switched on, a soft-start at reduced speed assures an even start-up and centring of the magnetic stir bars.

The electronically controlled automatic monitoring unit reduces the speed about every 40 seconds. In the process, each magnetic stir bar is re-centred in the bottle for a few seconds. Therefore, you can easily exchange individual bottles while the unit is switched on.

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### 9.1 Hinweis zu Copyright und Warenzeichen
Trademarks of the Tintometer Group
**Tintometer Inc.**