
Lovibond® Water Testing

Tintometer® Group

Turbidimeter

TB 350

(EN) Instruction manual.....	3	(DE) Gebrauchsanleitung.....	42
(ES) Manual de instrucciones.....	82	(FR) Mode d'emploi.....	122
(PT) Manual de instruções.....	164	(IT) Manuale d'istruzione.....	204
(NL) Gebruiksaanwijzing.....	244	(TR) Kullanım Talimatları.....	284
(RU) Инструкции по эксплуатации.....	324	(ZH) 使用说明.....	365

Table of Contents

	Page
1 Introduction	5
1.1 General Information	5
1.1.1 Read instructions before use	5
1.1.2 Purpose of the Product	5
1.1.3 Authorized use	5
1.1.4 Requirements for safe use	5
1.1.5 Keep instructions	6
1.1.6 User qualification	6
1.1.7 Handling of hazardous chemicals	6
1.1.8 Disposal notes	6
1.2 List of all used signs in the document	7
1.3 Abbreviations	7
2 Product overview	8
2.1 Certification	8
2.2 Features	9
2.3 Product description	9
2.4 Instrument view	9
3 Commissioning	11
3.1 Operating environment	11
3.2 Contents of delivery	11
3.3 Inserting and replacing the batteries	12
4 Operation	14
4.1 First start up	14
4.2 General operation principles	15
4.3 List of control elements and their function	16
4.3.1 Display Overview	16
4.3.2 Status Bar	17
4.3.3 Home Bar	19
4.4 Measurement	20
4.4.1 Measurement Overview	20
4.4.2 Performing a Measurement	21
4.4.3 Measurement Tips	22
4.4.4 Measurement Settings	23
4.5 Calibration	24
4.5.1 Calibration Overview	24
4.5.2 Performing a Calibration	24
4.5.3 Calibration Settings	25

	Page
4.6 Verification	26
4.6.1 Verification Overview	26
4.6.2 Performing a Verification	27
4.6.3 Verification Settings	27
4.7 Logs	28
4.7.1 Selecting a Log Type	28
4.7.2 Retrieving a Log	28
4.8 Instrument Settings	29
5 Maintenance	32
5.1 Recommended maintenance schedule	32
6 Troubleshooting	34
6.1 Error messages and what to do	34
7 Accessories & Replacement parts	36
7.1 List of Accessories	36
7.2 Replacement Parts	36
8 Specifications	37
8.1 Specifications - Mains adapter	39
9 Appendix	41
9.1 List of Trademarks	41

1 Introduction

1.1 General Information

1.1.1 Read instructions before use

This manual provides important information on the safe operation of the product. Please read this manual carefully and familiarize yourself with the product before use.

1.1.2 Purpose of the Product

The TB Series Portable Turbidimeter is designed to analyze aqueous samples that are collected into a glass vial that is then capped. The vial is prepared for measurement using a scrupulous technique to minimize interferences from dust and scratches. Once prepared, the vial is placed in the instrument for measurement.

The instrument can be used in the laboratory or taken to the field to perform measurements on samples as they are collected.

The handling of this turbidimeter should be done with care during any transport. If the instrument is taken into the field, the carrying case should be used. If used for regulatory applications, consultation with the local regulatory authorities is advised to ensure the correct model of instrument is selected.

1.1.3 Authorized use

The manufacturer's liability and warranty for damage is voided with improper use, failure to follow this manual, use by unqualified personnel, or unauthorized changes to the product.

The manufacturer is not liable for costs or damages that arise from the user or third parties due to the use of this product, especially in cases of improper use of the product or misuse or faults in the connection of the product.

The manufacturer assumes no liability for print errors.

1.1.4 Requirements for safe use

Note the following points for safe use of the unit:

- If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.
- The product may only be used according to the authorized use specified above.
- The product may only be supplied with power by the energy sources mentioned in this operating manual.
- The product may only be used under the environmental conditions mentioned in this operating manual.
- The product must not be opened or modified.

The product must not be used if:

- it is visibly damaged (e.g. after being transported)
- it was stored under adverse conditions for a lengthy period of time (storage conditions, see Chapter "Specifications")

Prerequisites for safe use of the plug-in power supply unit

Observe the following points for safe use:

- Check for visible damage before use. If damaged, do not put into operation.
- Use only in dry rooms.
- Handle only with dry hands.
- Only the power supply recommended by the manufacturer may be used.
- The product must not be opened or modified.
- Only the USB-C cable recommended by the manufacturer may be used to power the device.

1.1.5 Keep instructions

The manual must be kept in the vicinity of the product so you can always find the information you need.

1.1.6 User qualification

The operating personnel must be able to understand and correctly implement the safety labels and safety instructions on the packages and inserts of the products.

The user must be capable and able to read and understand this manual in order to familiarize themselves with the handling and to ensure safe use.

1.1.7 Handling of hazardous chemicals

Chemical and/or biological hazards may exist where this product is used. Abide by all governing laws, regulations and protocols when using this product.

For the development of products, Lovibond® pays close attention to safety. Some hazards from dangerous substances cannot be avoided. If self-produced tests or solutions are used, the responsibility concerning any risks caused by those tests or solutions lies with the user (personal responsibility).

For your own protection, for turbidity standards greater than 5 NTU, wearing protective gloves or goggles is recommended. Observe the relevant safety data sheet.

The safety datasheets of the chemicals comprise all instructions on safe handling, occurring hazards, preventive actions and actions to take in hazardous situations.

1.1.8 Disposal notes

Dispose of the batteries and electrical devices at a suitable facility according to local legal requirements.

It is illegal to dispose of the batteries with household waste.

Within the European Union, the batteries are removed at a specialized treatment center at the instrument's end of life.

Instruments marked with this symbol must not be disposed of in normal domestic waste.

1.2 List of all used signs in the document

The following symbols are used in this manual to identify sections that require special attention:

 Danger!	<p>A hazard exists that will result in death or severe injury if not avoided.</p>
 Warning!	<p>Improper handling of certain reagents can cause damage to your health. In any case follow the safety labels on the packing, the safety instructions of the package insert and available SDS. Protective measures specified there have to be followed exactly.</p>
 Caution!	<p>A hazard exists that may result in minor or moderate injury.</p>

1.3 Abbreviations

Abbreviation	Definition
NTU	Nephelometric Turbidity Unit
FNU	Formazin Nephelometric Units
mg/l	Milligram per liter

2 Product overview

2.1 Certification

EMC Conducted and radiated emissions	CISPR 11 (Class A Limits)	CE mark
EMC Immunity	EN 61326-1 (Industrial limits)	CE mark
Safety	EN 61010-1	TÜV safety mark
FCC	FCC Class A	FCC mark

FCC Class A Notice

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

- This device may not cause harmful interference.
- This device must accept any interference received, including interference that may cause undesired operation.

Note: This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at their expense.

Shielded Cables

Connections between the system and its peripherals must be made using shielded cables in order to maintain compliance with FCC radio frequency emission limits.

Caution!

Modifications

Any modifications made to this device that are not approved by Tintometer may void the authority granted to the user by the FCC to operate this equipment.

Caution!

DOC Class A Notice - Avis DOC, Class A

This Class A digital apparatus meets all requirements of the Canadian Interference-Causing Equipment Regulations.

2.2 Features

Laboratory accuracy meets portability



The Lovibond® TB Series eliminates the complexities in turbidity measurement. This new technology provides the highest level of accuracy, operational and regulatory efficiencies.


- Touchscreen interface for ease of operation
- Guided routine procedures to eliminate errors
- Straightforward data management protocols to ensure data integrity

2.3 Product description

Turbidity is an expression of the optical properties of a sample that causes light rays to be scattered and absorbed rather than transmitted in straight lines through a fluidic sample. In water, turbidity is often caused by the presence of matter such as clay, silt, finely divided organic matter, plankton, other microscopic organisms, organic acids, and pigments. Turbidity measurement is used to monitor processes that remove or alter the concentrations of particles in solutions, such as in the production of drinking water. These Lovibond® turbidimeters are designed to combine high accuracy and high sensitivity that are expected of a laboratory turbidimeter with portability.

2.4 Instrument view

Position	Description	Function
	1. Lid	When lifted, allows the user to place the sample cell in the instrument for a reading.
	2. Power button	Powers the instrument ON and OFF.
	3. Touch Screen	Displays the user interface of the instrument.
	4. Model Number	Identifies instrument version.
	5. USB-C Port	Supplies power to the instrument from an external source or will re-charge the battery pack accessory (if installed).
	6. USB-A Port	Supports external USB-A device for functions such as data transfer or backups.
	7	
	5	
	6	

Position	Description	Function
	7. Plug	Used to cover USB-C and USB ports when connections are not in use.

EN

3 Commissioning

3.1 Operating environment

The instrument can be used in any indoor or outdoor environment that is deemed safe for an operator to perform the analysis.

The instrument can accept aqueous samples from 0 to 70 °C.

The instrument should be level $\pm 15^\circ$ when performing a measurement.

The instrument should be placed on a level surface to achieve the best measurement performance.

EN


3.2 Contents of delivery

Inspect items to ensure no damage has occurred during shipment. If there is damage or something is missing, please contact local distributor immediately.

TB	350 IR	350 WL
Instrument	X	X
Silicone Oil, 15 ml	X	X
Microfiber Cloth	X	X
2 x Sample Cells	X	X
Screwdriver	X	X
Brush	X	X
Case with insert	X	X
Quick Start Guide	X	X
Certificate of Conformity	X	X
Warranty Booklet	X	X
4 x AA Batteries	X	X
Rechargeable Battery Pack		
USB-C Cable		
AC Power Supply		
Adapter for US, UK and EU Plugs		
T-CAL [®] Standard, 5 NTU	X	X
T-CAL [®] Standard, 20 NTU	X	X
T-CAL [®] Standard, 800 NTU	X	X
T-CAL [®] Standard, 2000 NTU		X
T-CAL [®] Standard, 4000 NTU	X	X

3.3 Inserting and replacing the batteries

When using batteries

EN


Step	Description
1	Use a screwdriver to remove the four screws on the back of the instrument and remove the panel.
2	Insert four AA batteries in the correct orientation.
3	Attach the panel back onto the instrument.

Caution!

- Only use the approved type: alkaline manganese battery, LR6, with 1.5V voltage.
- Remove empty batteries from the unit as soon as possible.
- If the unit is not used for a long time, remove the batteries from the unit.

When using rechargeable battery packs

Step	Description
1	Use a screwdriver to remove the four screws on the back of the instrument, remove the panel and the battery tray.
2	Insert the rechargeable battery pack in the correct orientation.
3	Attach the panel back onto the instrument.

Warning!

The rechargeable battery pack (19820-011) is only compatible with instruments with the serial number 23041000 and higher.

Using the battery pack for instruments with serial numbers prior to 23041000 may cause the device to overheat during charging.

This overheating can lead to incorrect readings and damage to the device and must therefore be avoided.

For devices with serial numbers prior to 23041000, please contact Lovibond® for details on the use of this accessory.

4 Operation

4.1 First start up

1. Welcome Screen

Initial screen user sees upon first power on.
Press > button

2. Select Language

Select desired language to set.
Press > button

3. Set Date and Time

Choose between 12- and 24-hour format.
Change values using up and down arrows.
Press > button

4. Password

Set a user-selected password to implement security features on instrument.
Press save button.
If no password is desired, select "Skip".
Press > button

5. Default Units

Select desired unit type to set.
Press > button

6. Contact Details

Enter instrument owner information.
Press save button.
If no details are desired, select "Skip".
Press > button


7. Summary Screen

Review of initial setting.
Press < to change setting before proceeding.
Press > to proceed with the Home Screen.
User can select task to perform with instrument.

4.2 General operation principles

The designs of the TB Series Portable Turbidimeters are broken down into their respective functional components that are discussed in detail. These include the light source, the 90-degree scattered light trap, the scattered light detector, turbidimeter sample vial and how these key components work together to deliver the turbidity result.

The simplified optical configuration for the TB Series of Portable Turbidimeters. This figure shows the incident light path and the view of the 90° scattered light detector.

- 1 Illuminator (LED and Optics)
- 2 Incident Light Path
- 3 90° Stray Light Trap
- 4 Sample Vial
- 5 Scattered Light (10°-20° Detector View)
- 6 Scattered Light Detector and Optics


4.3 List of control elements and their function

4.3.1 Display Overview


Description	Function
1 Status Bar	Displays the icons and labels that are universal including Time, Bluetooth®, USB, Security, Light Source & Power.
2 Action Area	Area that displays application steps and interaction.
3 Home Bar	Used primarily for navigation, to start or save or dismiss actions & return Home.

Icon	Description
	Measurement Mode
	Calibration Mode
	Verification Mode
	Data Logs

Icon	Description
	Bluetooth connection to PTV process turbidimeter (not available for TB 350)
	Instrument settings

4.3.2 Status Bar

Icon	Description	Function
	Application Name	Title of the app that the user is currently using.
	Time	Time set by user in initial setup or by accessing the instrument settings. 12 hour or 24 hour format is available.
	USB	If displayed, indicates a USB module is attached to the unit.
	Bluetooth Icon - Plain ²⁾	If displayed, Bluetooth® module is attached.
	Bluetooth Icon – With arrows ²⁾	If displayed, Bluetooth® module is attached and actively communicating with a PTV or AquaLX App.
	Red Light Source – Filled	Instrument contains a red LED light source; light is switched on for measurement.
	Red Light Source – Outline	Instrument contains a red LED light source; no measurement is active.
	White Light Source – Filled	Instrument contains a white LED light source; light is switched on for measurement.
	White Light Source – Outline	Instrument contains a white LED light source; no measurement is active.

Icon	Description	Function
	Infrared Light Source – On	Instrument contains an infrared LED light source; light is switched on for measurement.
	Infrared Light Source – Off	Instrument contains an Infrared LED light source; no measurement is active.
	Laser Light Source – Filled ¹⁾	Instrument contains a laser LED light source; light is switched on for measurement.
	Laser Light Source – Outline ¹⁾	Instrument contains a laser LED light source; no measurement is active.
	AC Battery Charging	AC Battery is charging.
	AC Power Icon	Indicates instrument is powered via AC power.
	Battery Power – Green	Indicates 100% life of battery remains.
	Battery Power – Yellow	Indicates 75% of battery life remains.
	Battery Power – Orange	Indicates 50% of battery life remains.
	Battery Power - Red	Indicates 25% of battery life remains.
	Sample chamber cover open	Indicates open Sample chamber cover.
	Sample chamber cover closed	Indicates closed Sample chamber cover.

¹⁾ only TB 600

²⁾ not TB 350

4.3.3 Home Bar


Icon	Description	Function
	Home Button	Returns user to home screen.
	Home Button - Locked	Security is active; password has been enabled and is required to perform certain functions.
	Home Button - Unlocked	Security is inactive; no functions are blocked to the user.
	Save Button	Saves the content displayed in the active area.
	Settings Button	Opens the setting menu of an application.
	Export Button	Exports the log(s) displayed in the active area.
	Measurement Button	Starts or Begins or Initiates a measurement.
	Measurement Button – with Bluetooth ^{®2)}	Takes a measurement while connected to a PTV unit.
	Bluetooth [®] Refresh Button ²⁾	Refreshes the list of available PTV units within proximity of the instrument.

Icon	Description	Function
	Cancel Button	Gives option to cancel current activity in the active area.
	Right Chevron	Advance to the next step within a process.
	Left Chevron	Return to a previous step within a process.

²⁾ not TB350

4.4 Measurement


4.4.1 Measurement Overview

Description	Function
1 Application Bar	Navigation path to help the user through the process steps.
2 Measurement Value	Display of measured sample value with unit type.
3 Settings	Settings applied to measurement are displayed.
4 Log Button	Creates a Log when pressed.
5 Measurement Button	Starts a measurement when pressed.
6 Settings Button	Used to change settings.
7 Home Button	Returns user to home screen.

4.4.2 Performing a Measurement

Please prepare the sample cells included in the scope of delivery before measurement. Make sure that the sample cells are free of stains, scratches or blemishes.

1. Collect sample

Ensure sample is representative (well mixed).
Glassware should be free of stains, scratches or other blemishes.

Rinse collection vessel at least two times with sample prior to sample collection.

2. Prepare sample cell

Using a lint free cloth, dry the exterior of the sample vial.


Use silicone oil to minimize stray light reflections.

3. Insert vial into instrument

Ensure sample is bubble free.

Invert slowly to mix. Do not shake.

4. Close lid

Measurement will not be allowed without lid in place.

5. Check settings

Ensure the measurement settings are correct.

Settings for measurement cannot be changed after measurement occurs.

6. Press measure button

Starts Measurement.

7. Data Log

Enter Data Log to save measurement.

8. Save measurement

If desired, press Save Button in order to create a measurement log

9. Pause

Pauses the measurement in "continuous" measurement mode.

- Do not bring any objects into the measuring chamber except the cuvettes intended for the measurements.
- Do not bring any liquids into the chamber.

Warning!

4.4.3 Measurement Tips

Indexing the vial: This is best performed with a vial that is filled with water that is filtered through a 0.2 µm filter.

- Press Measurement button.
- Rotate the sample. While rotating the vial, identify the lowest reading.
- When position is determined, stop the measurement. Mark the vial.
- Before proceeding with additional measurements, ensure measurement mode is updated, as necessary.
- Indexing the sample vial is recommended if there is a noticeable scratch or pit in the glassware.

Use of Silicone Oil

- Take care to handle the vial by the top third only.
- Wipe the vial with the supplied cleaning cloth to remove water, dust and fingerprints.
- Oil the outside of the vial with silicone oil.
 - Place a very thin bead of silicone oil down the side of the vial. Evenly spread the bead of oil over the bottom 2/3's of the vial.
 - Wipe the silicone from the vial through rotation of the vial against a clean, dust-free cloth.
- Oiling cells is critical and recommended when measuring samples with an expected turbidity that is below 0.5 NTU.

Samples containing excess bubbles

- Allow sample to stand for a minimum of five minutes.
- Gently invert sample two to three times to resuspend the turbidity.
- Place in instrument and proceed as instructed.

Other Notes

- Use of ultrasonic baths is not recommended. Ultrasonic energy can change particle distribution.
- Use Fast Setting Mode for samples that have high concentrations of sands or particles of high density.

4.4.4 Measurement Settings


Press the Settings Button to access the Measurement Settings Mode.

Icon	Description	Option	Function
	Mode	Single	Push to Read.
		Continuous ¹⁾	Measurement updates every X seconds until stopped.
		Fast Settling	Recommended for samples with rapid settling characteristics.
	Location	Default	Used if no Custom ID is set.
		Custom	Custom description of measurement location and measurement limits.
	Signal Avg.	On	If on, can be set up to 255 measurements. Function should be utilized if measurements are unstable.
		Off	deactivated
	Sample ID	Default	Used if no Custom ID is set.
		Custom	Unique identifier of sample; user selectable.

¹⁾not TB350

4.5 Calibration

4.5.1 Calibration Overview

Description	Function
1 Application Bar	Breadcrumb to help the user through the process steps.
2 Calibration Curve	Specifies calibration type.
3 Standard Set	Unique identifier of standard set.
4 Settings Button	Used to change settings.
5 Home Button	Returns user to home screen.
6 Next Button	Progress to next screen.


4.5.2 Performing a Calibration

Please prepare the T-CAL[®] standards included in the scope of delivery before calibration. Make sure that the T-CAL[®] vials are free of stains, scratches or blemishes.

- T-CAL[®] standards for TB 350 IR: 20 / 800 / 4000 NTU
- T-CAL[®] standards for TB 350 WL: 20 / 800 / 2000 / 4000 NTU

- 1. Zero Adjustment:**
Sets zero standard.

2. Prepare Vial:
Oil sample vial if necessary.

3. Insert Standard:
Insert standard value specified on screen.
User will insert series of various turbidity standards to complete calibration.

4. Modal:
Represents active process.

5. Log Summary:
Press Save Button to save log.
Press Cancel Button to return to beginning of calibration process.

4.5.3 Calibration Settings

Press the Settings Button to access the settings mode.

Icon	Description	Option	Function
	Cal Curve	T-CAL®	Selects prepared stabilized formazine standard at pre-defined concentrations
		User Formazin	Selects user prepared formazine standard at pre-defined concentrations
	Standard Set	Default	Used if no unique identifier is set.
		Unique	Unique identifier of standard set.
	Reminder ¹⁾	On	Reminds user that calibration is due. Option to force Calibration before additional measurements occur.
		Off	Default setting.
	Default Cal		Option to return to factory calibration setting.
	Post Verify	On	Requires user to perform a Verification after calibration is complete.
		Off	Default setting.


Cal Checker

On (Default)	During the calibration procedure, if an individual standard is out of tolerance, the user must correct the standard to proceed or cancel the calibration procedure. Feature prevents user from a calibration error due to an incorrect standard. The log screen will have a line item for Calibration Checker and it will show an "On" or "Enabled" state.
Off	During the calibration procedure, if an individual standard is out of tolerance, user will not be notified that a calibrator is out of Tintometers predefined range. The log screen will have a line item for Calibration Checker and it will show an "Off" or "Disabled" state.

¹⁾not TB350

4.6 Verification


4.6.1 Verification Overview

Description	Function
1 Application Bar	Breadcrumb to help the user through the process steps.
2 Verification Standard	Value of Standard to be measured.
3 Pass / Fail Criteria	Tolerance used for determining if a verification passes or fails.
4 Standard Set	Unique identifier of standard set.
5 Settings Button	Used to change settings.

Description	Function
1 Application Bar	Breadcrumb to help the user through the process steps.
6 Home Button	Returns user to home screen.
7 Next Button	Progress to next screen.

4.6.2 Performing a Verification

1. Prepare Vial

Oil sample vial if necessary.

3. Insert Standard

Insert standard value specified on screen.

4. Modal

Represents active process.

5. Log Summary

Press Save button to save log.

Press Cancel Button to return to beginning of Verification process.

4.6.3 Verification Settings

Press the Settings button to enter the settings mode.

Icon	Description	Option	Function
	Standard	T-CAL 5.0 NTU	Selects prepared stabilized formazin standard at 5.0 NTU concentration.
		User Prepared	Selects user prepared standard between 1 and 1000 NTU.

Pass/Fail criteria	$\pm 10\%$	Tolerance used for determining if a verification passes or fails.
	$\pm 5\%$	
	± 0.05 NTU	
	± 0.10 NTU	
	Custom	User entered value. Tolerance of $\pm 20\%$ or $\pm 0.05 - 0.2$ NTU allowed.

Reminder¹⁾	On	Reminds user that verification is due.
	Off	Default setting.

¹⁾ not TB350

4.7 Logs

4.7.1 Selecting a Log Type

Press the Log Type Button in the Action Area to access the Logs.

Application Bar:

Measurement Logs:

Displays all recorded measurements.

Calibration Logs:

Displays all calibration records.

Verification Logs:

Displays all verification records.

¹⁾

GLP Mode Logs:

Displays record of GLP mode being enabled or disabled.

¹⁾ not TB350

4.7.2 Retrieving a Log

1. Select log type

Displays the records of the type of log user would like to view.

2. Export Individual Logs

If a single record is selected, that log is exported in CSV format.
All details associated with log are exported.

3. Select Transmission Method

USB: Exports log to USB drive.

Bluetooth[®]¹⁾: Exports log to AquaLX mobile application.
Not available on TB 350.

Note: an external Lovibond[®] Bluetooth[®] dongle is required for connection to the USB-A port.

¹⁾ not TB 350

EN

4.8 Instrument Settings

Press the Settings Button to access the Settings Mode.

Icon	Description	Option	Function
	Security	Enable	Set up, verify, and change password.
		Disable	Default Setting.
	Power Settings	Screen Brightness	Adjust brightness of the screen.
		Backlight Sleep	Off 20 Seconds 30 Seconds 1 Minute* 5 Minutes
		System Shutoff	Off 1 Minute 2 Minutes 5 Minutes 15 Minutes*
			Note: If USB-C power cable is in use, System Shutoff is overridden to "Off" until instrument is unplugged, then will default back to selected option.

**Bluetooth¹⁾**

Enable

Enables Bluetooth® when hardware is present.

Disable

Default Setting.

**USB Shortcuts¹⁾**

Enable

If enabled, when a USB storage device is connected, the user is presented with a modal which presents an option for “transferring logs” or “backup/restore” settings. By selecting an option, the user is directed to the appropriate function.

Disable

Default Setting; users must access the “Logs” or “Backup/Restore” settings when USB storage device is connected.

**Language**

Select language.

**Date/Time**

Select time displayed and format.

**Backup/Restore**Backup All Settings¹⁾

Option to store settings to local SD Card or save to external USB drive.

Restore All Settings¹⁾

Recovers previously stored settings saved to local SD card or external USB drive.

Restore Factory Settings

Resets all settings to factory defaults. Saved logs are still available.

**Unit Type**

NTU	Select unit displayed.
FNU	Defines the unit type of the displayed measurement value.
mg.l Kaolin	NTU is Default setting.
mg.l PSL	
Degrees	
mNTU ¹⁾	
mFNU ¹⁾	
Custom ¹⁾	Allow user defined unit type and ability to set a factor based on NTU/FNU values.

**GLP¹⁾**

Enable	Establishes protocol for data integrity for laboratories following Good Laboratory Practices.
Disable	Default setting.

**Owner**

Contact details of instrument owner.

**Contact**

Contact information for Lovibond® offices.

**About**

About Lovibond®.

¹⁾ not TB 350

* Default Setting

5 Maintenance

5.1 Recommended maintenance schedule

Task	Do	Do not
1. Clean measurement chamber	1. Ensure the measurement chamber is clean and dry prior to measurement. Use a dry microfiber cloth to remove impurities.	1. Spray or pour any fluid directly into the measurement chamber.
2. Clean instrument display	2. Wipe touchscreen with a dry or slightly damp microfiber cloth.	2. Measure a sample that is not tightly capped.
Clean exterior housing	Clean the outer housing with a mild detergent solution or diluted alcohol.	
Sample Cell Storage & Care	Always rinse vials with turbidity-free water after use.	Use excessive silicone oil that has not been wiped from the vial. A single bead is all that is necessary.
	Fill vial with DI water. Cap vial to prevent formation of water spots.	Leave sample cells uncapped to air dry.
	Wipe outside of vial dry with a microfiber cloth.	Use any vial that has visible imperfections, including pits, scratches, or cracks.
Battery Replacement	It is suggested to remove batteries if instrument is stored for more than 30 days to prevent battery leakage.	Use unspecified type of battery.
Calibration	Calibrate as needed under the following conditions: <ol style="list-style-type: none"> 1. Verification failure. 2. Maintenance of optical or electronic systems are performed. 3. As required by regional regulatory authorities. 	Use expired standards.
	Carefully follow the instructions preparation and use of standards.	

Task	Do	Do not
Verification	Verify instrument performance under the following conditions: <ol style="list-style-type: none"> 1. After performing a calibration. 2. Per internal quality control procedures. 3. As required by regional regulatory authorities. 	Use expired standards.
	Carefully follow the instructions for preparation and use of standards.	
Connectors	Ensure connector plugs are in place when not in a laboratory environment.	Store instrument without connector plugs in place. Transport or support instrument via a connector cable.

**Warning!**

Improper handling of certain reagents can cause damage to your health. In any case follow the safety labels on the packing, the safety instructions of the package insert and available SDS. Protective measures specified there have to be followed exactly.

6 Troubleshooting

6.1 Error messages and what to do

Symptom	Possible Cause	Solution
Overrange	The turbidity of the sample is beyond the measuring range of the instrument.	<ul style="list-style-type: none"> • Dilute the sample • Verify calibration using a high-range standard.
Under range	The turbidity of the sample is below the detection limit of the instrument	<ul style="list-style-type: none"> • Repeat the stray light measurement (zero adjustment) as instructed ensuring the sample lid is closed. • Repeat the calibration after the stray light was determined.
Instrument does not turn on	Batteries are defective Batteries are not inserted correctly	<ul style="list-style-type: none"> • Replace batteries with a complete new set.
“Save” button for a function does not appear, so a log cannot be created.	Required information is missing	Ensure that the user initials and all other required fields are properly addressed.
Instrument will not take a measurement	Lid is not closed	Close the lid.
Low battery notification	Battery is below 20%	Replace or recharge batteries
Instrument settings cannot be accessed by user	Instrument security protocol is in place	Enter the password to unlock settings and calibration functions
Calibration Error	Sample vials were not cleaned prior to calibration	<ul style="list-style-type: none"> • Re-wash sample vials using the procedure in the manual. Store vials as prescribed in the manual. • Check expiration date of standard
Calibration Error	Calibration standards were ran in the wrong order	Repeat calibration with standards run from lowest to highest value.

Symptom	Possible Cause	Solution
Calibration Error	Calibration standards were not adequately mixed.	Mix standards as instructed immediately prior to performing the calibration.
Calibration Error	Calibration Standards prepared to wrong turbidity values	Repeat preparation of calibration standards, making sure dilutions are correct.
Verification Error	Verification falls outside the Error boundaries	<ul style="list-style-type: none"> • Mix the verification standard prior to verification. • Check pass/fail criteria to ensure the criteria is not tighter than the accuracy of the standard and adjust these criteria accordingly. • Check expiration date of standard

7 Accessories & Replacement parts

7.1 List of Accessories

Title		Part Number
Batteries (AA), set of 4	TB350 IR, TB350 WL	1950025
Brush, 11 cm length	TB350 IR, TB350 WL	380230
Cleaning cloth	TB350 IR, TB350 WL	197635
Cuvette stand for 6 round cuvettes Ø 24 mm	TB350 IR, TB350 WL	418951
Power supply TB series	TB350 IR, TB350 WL	19820-170
Rechargeable battery pack, Ni-MH	TB350 IR, TB350 WL	19820-011
Sample cuvettes with black lid, Height 55 mm, ø 24 mm, set of 12	TB350 IR, TB350 WL	197655
Set turbidity standards T-CAL (5, 20, 800, 2000, 4000 NTU)	TB350 WL	194152
Set turbidity standards T-CAL (5, 20, 800, 4000 NTU)	TB350 IR	194154
Silicon oil	TB350 IR, TB350 WL	194295
T-CAL®-Standard, 4000 NTU, 125 ml	TB350 IR, TB350 WL	48012912
T-CAL®-Standard, 4000 NTU, 500 ml	TB350 IR, TB350 WL	48012950
USB-C cable 1 m, USB-C to A	TB350 IR, TB350 WL	19820-081

7.2 Replacement Parts

Description	Part Number
Carrying Case	19820-130
Lid, sample chamber	19820-095
Cover, Battery compartment	19820-098
Plug, Elastomeric, IP67, I/O	19820-084
Insert for 4 AA batteries	19820-012

8 Specifications

Model	Model WM-040-1004 (IR) Model WM-040-1001 (WL)
Light Source	IR LED (860 nm) White LED
Regulatory Compliance	ISO 7027 EPA approved (alternative to US EPA 180.1)
Measuring Principle	Nephelometric (Multipath 90° BLAC)
Measuring Range	0.01 - 4000 NTU (FNU)
Resolution	Auto ranging: 0.01-9.99: 0.01 NTU ; 10.0-99.9: 0.1 NTU ; 100-4000: 1 NTU
Accuracy	± 1.8 % of reading + straylight
Repeatability	< 1 % or ± 0.01 NTU
Stray Light	< 0.014 NTU
Display Units	NTU, FNU, Degrees, mg/L Kaolin, mg/L PSL
Response Time	7 seconds
Reading Modes	Single, Signal Averaging, Fast-Settling
Calibration Options	Full-range calibration from 0 - 4000 NTU
Detector	Silicon Photodiode
Ratio Mode	Yes
Sample Required	10 mL (0.4 oz)
Display	Full-Colour Touchscreen
Display Size	55 x 95 mm (W x H)
Interfaces	<ul style="list-style-type: none"> • USB-A • USB-C
Data Logger	250 measurements, all calibrations, verifications and factory-restore actions
Data Transfer Format	.csv
Power Save Mode	Yes
Auto – OFF	Selectable

Power Supply	<ul style="list-style-type: none"> • 3 NiMH-battery pack (AA) • 4 batteries (AA) • Universal (90 - 240 VAC), via external power supply
Environmental Conditions	<p>Temperature: 0 - 50 °C (32 - 122 °F)</p> <p>Humidity: 0 - 90 % at 30 °C, noncondensing</p> <p>Humidity: 0 - 80 % at 40 °C, noncondensing</p> <p>Humidity: 0 - 70 % at 60 °C, noncondensing</p>
Protection Class	IP 67 (Instrument)
Compliance	<ul style="list-style-type: none"> • CE • CSA • FCC A • RoHS • TÜV • UL • WEE
Languages User Interface	<ul style="list-style-type: none"> • Chinese • Dutch • English • French • German • Japanese • Polish • Portuguese • Spanish • Turkish
Languages Quick Start Guide	<ul style="list-style-type: none"> • Chinese • Dutch • English • French • German • Italian • Portuguese • Russian • Spanish • Turkish

Languages Full User Manual	<ul style="list-style-type: none"> • Chinese • Dutch • English • French • German • Italian • Portuguese • Russian • Spanish • Turkish
Dimensions	155 x 83 x 225 mm
Weight	804 g (898 g with four AA alkaline batteries)
Weight with Packaging	(898 g with four AA alkaline batteries)

**Caution!**

Subject to technical modification!
To ensure maximum accuracy of test results, always use the reagent systems supplied by the instrument manufacturer.

8.1 Specifications - Mains adapter

Type	WR9QA2500USB52NMR6B
Design	Mains adapter plug, switching power supply
Input voltage, frequency	100 – 240 V \pm 10 %, 50 / 60 Hz \pm 5 %
Input current	0.3 A
Primary adapter	Europe, UK, Australia, USA
Protection class	II
Output voltage, frequency	5.2 V, DC
Max. output current	2.5 A
Output protection	Output shut down and auto restart
Environmental conditions	0 °C to 40 °C, with max. 95 % relative air humidity
Energy efficiency class	VI
Safety standard	EN 62368
EMC	EN 55032 FCC Part 15 Class A EN 55024
Dimensions [mm]	31.5 x 41 x 71 (without primary adapter)
Weight, with EU adapter	48 g (without primary adapter blades)

Type **WR9QA2500USB52NMR6B**

Approval, EMC CE, FCC

Permitted altitude meters 5000 m

Pollution level 2

**Caution!**

Subject to technical modification!
To ensure maximum accuracy of test results, always use the reagent systems supplied by the instrument manufacturer.

EN

9 Appendix

9.1 List of Trademarks

Lovibond®, Tintometer® and T-CAL® are registered trademarks of the Tintometer group of companies. All translations and transliterations of Lovibond® and Tintometer® are asserted as trademarks of The Tintometer® Group.

The Bluetooth® mark, figure mark and combination mark are trademarks owned by the Bluetooth SIG and any use by the Tintometer® Group is under license.

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9.1 Hinweis zu Copyright und Warenzeichen
Trademarks of the Tintometer Group
Tintometer Inc.