
# **Dipslides Information leaflet**

Leaflet-UC



## Procedures

Remove the dipslide from the container by pulling the white tab while holding on to the outer cover.

### Sampling Fluids

The dipslide is then immersed into the sample ensuring both sides of the paddle are submerged into the fluid to be tested.

### Sampling Surfaces



The sample should be taken by allowing direct contact between the agar surface and the test material. The flexible paddle can allow both surfaces to come into contact with the surface. Sweep an area approximately twice that of the paddle to give a more accurate result.

## Procedures

Excess fluid sample should be gently shaken from the paddle before it is placed back in to the container.

After sampling, the dipslides should be placed immediately into incubation at 30-35°C for 24-48 hours. For fungi and moulds incubate for up to 5 days.

After incubation compare results of the dipslide to the dipslide comparison chart to quantify the results.

Our dipslides are designed to allow easy enumeration of aerobic bacteria, fungi, yeasts and moulds from fluids and surfaces. Nutrient TTC agar is used for the general cultivation of aerobic organisms and allows for a total viable count (TVC) of a sample.

### **Storage and Product care**

Please store in a cool and dry environment at 10 - 15°C. Large variation in temperature can lead to condensation forming inside the tube and can eventually lead to dehydration of the agar.

### **Disposal**

Ensure to keep vials closed and sterilise before disposal. Alternatively, dispose via incineration. Please ensure compliance with local health and safety procedures. This product has been made in Britain to the highest standards, incubated and quality checked.

**Dipslides**

**GB**

## **Sampling**

### **Fluids**

Remove the dipslide from the container by pulling the white handle while holding on to the outer cover. The sample should be taken by immersing both sides of the paddle into the fluid to be tested. Excess sample should be gently shaken from the paddle before it is replaced in the container.

### **Surfaces**

The sample should be taken by allowing direct contact between the agar surface and the test material. The paddle is flexible and can be bent at the upper end to allow both surfaces to come into intimate contact. Bacterial recovery rate is about 50% so that sweeping an area approximately twice that of the paddle will give a more accurate result.

### **Incubation**

After sampling, the dipslides should be placed immediately into incubation at 30 - 35°C for 24 - 48 hours. For fungi and moulds incubate for up to 5 days.

<b>Code</b>	<b>Type</b>	<b>Description</b>	<b>Growth</b>
56B010110	D001 TTC/TTC	TVC (Aerobic bacteria)	Red
56B010210	D002 TTC/MALT	TVC/Fungi & Moulds	Red/Green, Black, White: Beige (Yeast)
56B010310	D003 TTC/ROSE	TVC/Fungi & Moulds	Red/Green, Black, White: Beige (Yeast)
56B010410	D004 TTC/MAC	TVC/Coliforms	Red/Red or Yellow
56B010510	D005 TTC/E.COLI	TVC/ Coliforms & <i>E.coli</i>	Red/ Red or Yellow(Coliforms)or Blue ( <i>E.coli</i> ) or Cream ( <i>P.aeruginosa</i> )
56B010610	D006 MAC/PDM	Coliforms/ Pseudomonas	Red/Green ( <i>P.aeruginosa</i> ) or Cream ( <i>P.spp</i> )
56B010710	D007 TTC/PDM	TVC/ Pseudomonas	Red/ Green ( <i>P.aeruginosa</i> ) or Cream ( <i>P.spp</i> )
56B010810	D008 SRB Tube Test	Sulphate Reducing Bacteria	Black Diffusion
56B010910	D009 NRB Tube Test	Nitrite Reducing Bacteria	Pink Diffusion
56B011110	D011 R2A/R2A	Aerobic bacteria	Red
56B010010	D012 TSA/MALT	Aerobic bacteria/ Fungi & Moulds	Red/Green, Black, White: Beige (Yeast)

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