Lovibond® Colour Measurement

Tintometer® Group



Filtersatz EP Rot



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

Filtersatz EP Rot

A set of two glass conformance filters for checking PFXi instruments in European Pharmacopoeia Red scale.

Bestell-Nr.: 139410

Technische Daten

Glass Standard Values - R2 - R6

Filters

Lieferumfang

- EP Red 2 filter
- EP Red 6 filter
- Certificate of Conformity