

# Lovibond® Water Testing Tintometer® Group



## DI10 Incubator

For controlled bacterial growth



In order to compare the results of Dipslides and other microbiological tests such as the Coliform / E-coli Test Kit from week to week, the incubation time and temperature must be the same each time the test is performed. This ensures that bacterial growth is always controlled and consistent. The DI10 incubator is therefore an important and necessary step in the analysis of microbiological activity in water. The Lovibond® DI10 incubator provides reliable incubation in the field, in the laboratory or even on the road. For incubation of dipslides we recommend the purchase of the dipslide holder. This provides a secure support for 12 dipslides in the DI10 during the incubation period.

- Robust design
- Holds up to 12 dipslides or 10 quanti-discs
- Excellent temperature stability
- In-car operation

Part Number: 56B000701

### Industry

Chemical Industry | Food and Beverage Industry | Industries Others | Marine Industry | Oil and Gas | Pharmaceutical Industry | Power and Energy

### Technical Data

<b>Voltage</b>	12 V
<b>Display</b>	1 x 16 bit Backlit LCD
<b>Power Supply</b>	12 V 110 V 240 V
<b>Temperature Tolerance</b>	± 0.5 °C
<b>Environmental Conditions</b>	5 °C – 60 °C
<b>Compliance</b>	CE
<b>Dimensions</b>	246 x 215 x 162 mm

# Delivery Scope

- DI10
- Power supply
- Power cable set (UK, EU and USA)
- Cigarette lighter cable & 8 plug-in adapters
- User manual

Title	Part Number
Power cable set	176300
External power supply	182049

## Tintometer GmbH

Lovibond® Water Testing  
Schleefstraße 8-12  
44287 Dortmund  
Tel.: +49 (0)231/94510-0  
sales@lovibond.com  
www.lovibond.com  
Germany

## The Tintometer Limited

Lovibond House  
Sun Rise Way  
Amesbury, SP4 7GR  
Tel.: +44 (0)1980 664800  
sales@lovibond.uk  
www.lovibond.com  
UK

## Tintometer China

9F, SOHO II C.  
No.9 Guanghualu,  
Chaoyang District,  
Beijing, 100020  
Customer Care China Tel.: 4009021628  
Tel.: +86 10 85251111 Ext. 330  
Fax: +86 10 85251001  
chinaoffice@tintometer.com  
www.lovibond.com  
China

## Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,  
Lebu Nilam 2, Bandar Bukit Tinggi,  
Klang, 41200, Selangor D.E  
Tel.: +60 (0)3 3325 2285/6  
Fax: +60 (0)3 3325 2287  
lovibond.asia@tintometer.com  
www.lovibond.com  
Malaysia

## Tintometer Brazil

Caixa Postal: 271  
CEP: 13201-970  
Jundiaí – SP  
Tel.: +55 (11) 3230-6410  
sales@lovibond.us  
www.lovibond.com.br  
Brazil

## Tintometer Inc.

6456 Parkland Drive  
Sarasota, FL 34243  
Tel: 941.756.6410  
Fax: 941.727.9654  
sales@lovibond.us  
www.lovibond.us  
USA

## Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2<sup>nd</sup>, 3<sup>rd</sup> & 4<sup>th</sup> Floor  
Sanathnagar Industrial Estate,  
Hyderabad, 500018  
Telangana  
Tel: +91 (0) 40 23883300  
Toll Free: 1 800 599 3891/ 3892  
indiaoffice@lovibond.in  
www.lovibondwater.in  
India

## Tintometer Spain

Postbox: 24047  
08080 Barcelona  
Tel.: +34 661 606 770  
sales@tintometer.es  
www.lovibond.com  
Spain