

Lovibond® Colour Measurement

Tintometer® Group



Benchtop Stand



For use with the Lovibond® RT300, RT400 & RT500

Product Features and Highlights

- Simplified, expandable accessory approach
- No need for tools to install or remove the instrument
- Installation time approximately 30 seconds
- One design for RT300,-400,-500 series
- Enables vertical or horizontal sample plane orientation
- Removable, dampened sample arm
- Sample support shelf fully adjustable
- Securable to bench or desk top environment
- Open sample preview capability
- Standard adaptors available to measure
 - > Powders > Liquids > Gels > Granules >

The **Lovibond® 400620 Benchtop Stand** is designed to hold the Lovibond® RT300, 400 & 500 Series of Spectrophotometers and samples steady while taking measurements. The fixture ensures consistent and repeatable colour measurements and can be used with other accessory-measuring attachments.

Dimensions: W: 18.42cm; L: 25.4cm; H: 27.31cm

The Tintometer Ltd • Lovibond House • Sun Rise Way • Amesbury • SP4 7GR • UK

Tel: +44 (0)1980 664800 • Fax: +44 (0)1980 625412 • Email: sales@tintometer.com • Website: www.lovibondcolour.com

Lovibond® & Tintometer® are Registered Trade Marks of The Tintometer Limited. All translations and transliterations of Lovibond® & Tintometer® are asserted as Trade Marks of The Tintometer Limited. Registered in England Number 45024. Registered Office: Lovibond House • UK. 920390_V7_08/16