

Lovibond® Colour Measurement



Tintometer® Group



PFX_i - What's in the Box?

- ✓ PFX_i series instrument
- ✓ External power supply
- ✓ Set of 3 power leads (UK, European and United States)
- ✓ USB and LAN Cable
- ✓ Quick start guide
- ✓ CD ROM (Containing Windows control software and manuals)
- ✓ RCMSi Starter Pack
- ✓ Certificate of Conformity
- ✓ Accessory Box: Containing 1 of each item specified in table below for the instrument type

www.lovibondcolour.com

Description																													
	195\1	195\2	195\3	195\4	195\5	195\6	195\7	195\8	195\9	880\L	880\L + Heater	880\AT	880\AT + Heater	880\P	880\P + Heater	880\S	880\S + Heater	880\IP17	880\IP17 + Heater	880\CIE	880\CIE + Heater	950	950 + Heater	950\P	950\P + Heater	995	995 + Heater	995\P	995\P + Heater
5 Volt Halogen Lamp	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Conformance Filter	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
W600/OG/10 Cell	●		●	●		●	●	●	●											●		●					●		●
W600/B/10 Cell																					●		●				●		●
W600/OG/25 Cell				●																									
W600/OG/33 Cell		●												●											●				●
W600/B/33 Cell															●										●				●
W600/OG/40 Cell					●	●																							
W600/OG/50 Cell	●	●			●	●	●	●																					
W600/OG/100 Cell														●		●									●		●		●
W600/B/100 Cell															●		●								●		●		●
W600/OG/1" Cell										●		●						●				●				●		●	
W600/B/1" Cell											●		●						●				●			●		●	
W600/OG/5¼" Cell										●		●										●				●		●	
W600/B/5¼" Cell											●		●									●			●		●		●
W600/OG/6" Cell																		●										●	
W600/B/6" Cell																			●									●	

Genuine Lovibond® Cells are supplied with each instrument. Replacement Lovibond® cells may be purchased by quoting the description (e.g. W600/B/10). The cells are coded as follows:

- W600/** = Type of cell – size/height etc
- OG/** = Optical glass
- B/** = Borosilicate glass for high temperature samples
- 10** = 10 mm path length
- 50** = 50 mm path length
- 5¼"** = 5¼" path length etc...

The Tintometer Ltd
T. +44 (0)1980 664800
F. +44 (0)1980 625412

Tintometer North America
T: +1 (0)941 758 8671
F: +1 (0)941 727 9654
Toll Free: 855-Color12
(855-265-6712)

Tintometer – South East Asia
T: +603 3325 2286
F: +603 3325 2287

Tintometer – China
T: +86 10 8525 1111 Ext 330
F: +86 10 8525 1001