Lovibond® Colour Measurement

Tintometer® Group



Betacaroteno



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The

values achieved and the performance tolerances expected will, of

course, always be reported.

Betacaroteno

A glass conformance filter to check the Model Fx range of instruments accuracy in measuring beta Carotene.

Referencia No: 169520

Datos técnicos

Glass Standard – 45 ppm at 10 mm path length Values

Filters

Volumen de suministro

45ppm filter

Certificate of Conformity