



Water Testing

T-CAL® Turbidity Standards

suitable for all
turbidimeter
brands



Get the full range of formazin standards

Calibration of turbidimeters has never been easier!

Valid measurement is not only an important requirement for you in drinking water and waste water monitoring, it is a matter of course.

Whether for routine turbidity measurements in the laboratory or directly on site: use turbidity standards for all turbidimeters to save additional work.

With Lovibond® you get:

- Primary standards based on formazin for all applications
- Ready to use standards straight from the bottle
- Precise solutions to verify or calibrate your instrument
- Stabilised standards with long shelf life
- The perfect tool to assure your analytical quality
- Expert know-how based on about 140 years of experience

T-CAL[®] our solution for you

ISO & US EPA
compliant

T-CAL[®] long time stable turbidity standards for calibration & verification of any turbidimeter

Get the whole range of formazin based turbidity primary standards for turbidimeters of any manufacturer with our T-CAL[®] series. Ready to use, prepared in a wide variety of packaging sizes.



Single Standard	Quantity	Code				
T-CAL [®] -Standard <0.1 NTU	125 mL	48012012	500 mL	48012050	1000 mL	48012099
T-CAL [®] -Standard 0.3 NTU			500 mL	48011050		
T-CAL [®] -Standard 1.0 NTU	125 mL	48011112	500 mL	48011150	1000 mL	48011199
T-CAL [®] -Standard 5.0 NTU			500 mL	48012250		
T-CAL [®] -Standard 10 NTU	125 mL	48011212	500 mL	48011250		
T-CAL [®] -Standard 20 NTU	125 mL	48012312	500 mL	48012350	1000 mL	48012399
T-CAL [®] -Standard 100 NTU	125 mL	48011512	500 mL	48011550		
T-CAL [®] -Standard 200 NTU	125 mL	48011612	500 mL	48011650		
T-CAL [®] -Standard 800 NTU	125 mL	48011712	500 mL	48011750		
T-CAL [®] -Standard 1000 NTU	125 mL	48011812	500 mL	48011850		
T-CAL [®] -Standard 4000 NTU	125 mL	48012912	500 mL	48012950		



Get the perfect set for your Hach[®] instrument. Choose the appropriate composition to meet your requirements:

Standards in Set	Set Quantity	Code
T-CAL [®] -Standard Set	<0.1 / 20 NTU, each 2x 1000 mL	48019399
T-CAL [®] -Standard Set	<0.1 / 20 / 200 / 1000 / 4000 NTU, each 125 mL	48019712
T-CAL [®] -Standard Set	<0.1 / 20 / 200 / 1000 / 4000 NTU, each 500 mL	48019750
T-CAL [®] -Standard Set	<0.1 / 20 / 100 / 800 NTU, each 125 mL	48019412
T-CAL [®] -Standard Set	<0.1 / 20 / 100 / 800 NTU, each 500 mL	48019450
T-CAL [®] -Standard Set	10 / 20 / 100 / 800 NTU, each 125 mL	48019512
T-CAL [®] -Standard Set	10 / 20 / 100 / 800 NTU, each 500 mL	48019550

HACH[®] is a registered trademark of Hach Company, Loveland, Colorado. The use of the HACH[®] trademark does not imply any affiliation with or approval by Hach Company regarding the formulation, testing or compatibility of these products for use in HACH[®] brand spectrophotometers or other instruments or systems.