

# Lovibond<sup>®</sup> Colour Measurement Tintometer<sup>®</sup> Group



## AOCS-Tintometer Colour 3.0R 28.0Y (5¼")

Part Number: 134280



Reference liquid standards are assigned a nominal value and Lovibond Tintometer always endeavours to match this value as closely as possible. On occasion, it may not be possible to match the requested value exactly. The value achieved and the expected performance tolerances will, of course, always be reported.

### AOCS-Tintometer Colour 3.0R 28.0Y (5¼")

500ml Bottle certified under ISO 9001 quality system with a 12-month shelf life.

Full traceability to internationally recognised standards

All classified as non-hazardous according to EU directives

Each bottle supplied with full certification including MSDS