

Lovibond[®] Colour Measurement Tintometer[®] Group



LC 100 Spectrocolorimeter

Hand-Held Imaging Spectrocolorimeter



- Handheld and Lightweight Spectrocolorimeter providing Fast, Precise and Accurate Colour Measurement
- Up to 350 Sample Measurements can be Stored and Compared to Reference Standards
- Conformance Reports Available in PDF and Excel Formats
- Precise Positioning on Sample with Preview of Measurement Area

Part Number: 401100

Handheld and Portable

The LC 100 offers the ultimate in portability; its compact and lightweight design ensures the instrument can be with the user at all times, offering immediate access to accurate colour measurement. Fitting in the palm of your hand the LC 100 ensures easy one-handed operation anywhere on site for a prolonged period of time.

Accuracy

Despite its lightweight and compact design, the LC 100 does not compromise on accuracy. Using 8 different visible illuminations and 1 UV LED along with the 45°/0 optical geometry and image capture technology, the user can feel confident in the results they are receiving.

Reporting Results

To make reporting as quick and easy as possible, the LC 100 automatically creates conformance reports in either PDF or Excel formats. The user can then share these with a PC or directly to a compatible printer.

Storage

The LC 100 can store up to 350 sample measurements, these samples can then be tagged via the LC 100's on-screen keyboard or built in voice recorder. The LC 100 can also store up to 30 reference standards by date and time, ensuring a high level of quality control when comparing results.

Suitable for Multiple Applications

The versatility of the LC 100 means its suitable for multiple applications and has the ability to report with a display of ΔE for all common methods i.e. CIELAB, CMC, CIE 94, 2000, plus a display of measured and delta differences with $L^*a^*b^*C^*h^\circ$ for standard and sample, graphical $L^*a^*b^*$ plot.

Industry

Chemical Industry | Food and Beverage Industry | Industries Others | Oil and Gas | Pharmaceutical Industry

Application

Chemicals | Coatings | Food and Beverage | Household Products | Medical and Clinical Products | Pharmaceuticals and Cosmetics | Plastics | Sugar Solutions, Syrups and Honey | Textiles

LC 100 Spectrocolorimeter

The LC 100 is a highly flexible and lightweight spectrocolorimeter. This handheld unit offers the user the ability to capture fast and accurate colour measurement throughout the manufacturing process. In just over a second, the LC 100 will provide a sample reading along with a comprehensive report of how close the colour is to a self-programmed standard.

Technical Data

Applicable Standard	ASTM E308
Colour Type	Opaque
Display	4.5 cm Colour TFT
Interfaces	USB
Measuring Geometries	45/0 Image Capture
Internal Storage	30 standards/350 sample results
Measurement Time	1.8 seconds
Storage Temperature	-20° to 60°C (-4°F to 104°F)
Short Term Repeatability	Typical 0.10 ΔE 94 on white (D65/10)
Illuminants	CIE Illuminant A CIE Illuminant D65
Portability	Portable

Delivery Scope

- LC 100 Handheld Spectrocolorimeter
- Carrying Pouch
- UBS Cable for Power Supply and Data Transfer

Tintometer GmbH

Lovibond® Water Testing
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
sales@lovibond.com
www.lovibond.com
Germany

The Tintometer Limited

Lovibond House
Sun Rise Way
Amesbury, SP4 7GR
Tel.: +44 (0)1980 664800
Fax: +44 (0)1980 625412
sales@lovibond.uk
www.lovibond.com
UK

Tintometer China

9F, SOHO II C.
No.9 Guanghualu,
Chaoyang District,
Beijing, 100020
Customer Care China Tel.: 4009021628
Tel.: +86 10 85251111 Ext. 330
Fax: +86 10 85251001
chinaoffice@tintometer.com
www.lovibond.com
China

Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,
Lebuh Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com
www.lovibond.com
Malaysia

Tintometer Brazil

Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@lovibond.us
www.lovibond.com.br
Brazil

Tintometer Inc.

6456 Parkland Drive
Sarasota, FL 34243
Tel: 941.756.6410
Fax: 941.727.9654
sales@lovibond.us
www.lovibond.us
USA

Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2nd, 3rd & 4th Floor
Sanathnagar Industrial Estate,
Hyderabad, 500018
Telangana
Tel: +91 (0) 40 23883300
Toll Free: 1 800 599 3891/ 3892
indiaoffice@lovibond.in
www.lovibondwater.in
India

Tintometer Spain

Postbox: 24047
08080 Barcelona
Tel.: +34 661 606 770
sales@tintometer.es
www.lovibond.com
Spain

Technical changes without notice

Printed in Germany

Lovibond® and Tintometer® are Trademarks of the Tintometer Group of Companies