

# Lovibond® Water Testing Tintometer® Group



## Pooltester Chlorine LR/pH in box

The fast assistants for Pool Testing



- Easy to use
- Futuristic design
- RAPID tablets fast dissolving
- Highest accuracy

Part Number: 151600

### High measuring accuracy

Results are shown with high measurement accuracy. They have a certainty and can counteract if necessary.

### Intuitive

Self-explanatory is the operation of the tester. And if you are ever unsure: A look at the packaging and the intuitively understandable instructions are sufficient to eliminate the ambiguity.

### Ergo mixed

The ergonomically optimized lid ensures a secure grip.

### Industry

Marine Industry | NGO | Pool Domestic | Pool Services

### Application

Disaster Management | Pool Water Control

### Pooltester Chlorine LR/pH in box

The pool tester enables the simultaneous measurement of the water treatment agent and the pH value.

### Measuring Range

Test Name	Measuring Range
Chlorine	0.1 - 3.0 mg/L Cl <sub>2</sub>
pH	6.8 - 8.2 pH

## Delivery Scope

- Pooltester in a sturdy plastic box
- Tablet reagents for 20 tests
- Instruction manual