

Lovibond® Water Testing

Tintometer® Group



Legionella Swab Test Kit

The World's Fastest *Legionella* Test



- Fast, easy and accurate
- Results in as little as 25 minutes!
- Confirms performance of *Legionella* prevention and treatment protocols
- Test Kits available for detecting both waterborne and biofilm *Legionella*

Part Number: 56B006401

Legionella is the causative agent of Legionnaire's disease, a potentially fatal form of pneumonia. Guidelines exist across the world for the control of *Legionella*.

In 2002, the Health and Safety Executive released guidelines for the control of *Legionella*, known as the L8 guidelines. These guidelines detail the key parameters that should be tested on a weekly and monthly basis to minimise the risk of *Legionella* proliferating to harmful levels in open water system.

The Lovibond® *Legionella* "L8" Test Kits contain the recommended tests for both weekly and monthly monitoring.

Legionella Swab Test Kit

Within a water system *Legionella* usually grow within biofilm (slime) found on the insides of pipe work, shower heads, cooling tower packing, or water tanks. Bacteria found in the planktonic phase (i.e. suspended in the water itself) are usually an indicator of colonization of the biofilm.

Since a clean water sample does not equate to a clean system (proliferation of *Legionella* is greater in a biofilm than in water), this kit can be used to pinpoint the contamination in a system or to demonstrate cleanliness of an area after a targeted cleaning.

Use this kit to find *Legionella* in their breeding ground -- this very selective test can detect legionella even in samples with very high background flora. With results in about 25 minutes, the test allows for immediate understanding and response to the water system conditions. Regular use will increase confidence in your *Legionella* control procedures and reduce risk of outbreaks and litigation!

Limit of Detection
200 CFU/l / Swabbed area

No. of Tests
5

Delivery Scope

- 5 Test Strips with Dispensing Pipettes
- 5 Swabs
- Recovery Buffer
- 5 Reaction tubes
- Instruction Manual
- Carrying Case

Tintometer GmbH

Lovibond® Water Testing
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
sales@lovibond.com
www.lovibond.com
Germany

The Tintometer Limited

Lovibond House
Sun Rise Way
Amesbury, SP4 7GR
Tel.: +44 (0)1980 664800
Fax: +44 (0)1980 625412
sales@lovibond.uk
www.lovibond.com
UK

Tintometer China

9F, SOHO II C.
No.9 Guanghualu,
Chaoyang District,
Beijing, 100020
Customer Care China Tel.: 4009021628
Tel.: +86 10 85251111 Ext. 330
Fax: +86 10 85251001
chinaoffice@tintometer.com
www.lovibond.com
China

Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,
Lebu Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com
www.lovibond.com
Malaysia

Tintometer Brazil

Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@lovibond.us
www.lovibond.com.br
Brazil

Tintometer Inc.

6456 Parkland Drive
Sarasota, FL 34243
Tel: 941.756.6410
Fax: 941.727.9654
sales@lovibond.us
www.lovibond.us
USA

Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2nd, 3rd & 4th Floor
Sanathnagar Industrial Estate,
Hyderabad, 500018
Telangana
Tel: +91 (0) 40 23883300
Toll Free: 1 800 599 3891/ 3892
indiaoffice@lovibond.in
www.lovibondwater.in
India

Tintometer Spain

Postbox: 24047
08080 Barcelona
Tel.: +34 661 606 770
sales@tintometer.es
www.lovibond.com
Spain