

Lovibond® Water Testing Tintometer® Group



Chlorine HR Strips

For fast determination with chemical reaction



- Hazard free
- Rapid
- Always to hand
- Semi quantitative results

Part Number: 56S000190

Industry

Industries Others

Application

Disinfection Control

Chlorine HR Strips

Test strips offer quick, clean results without the fuss of reagents or tablets. Their dip and compare method makes them a popular choice for many applications where speed and simplicity are most important.

Measuring Range

Test Name	Measuring Range
Chlorine HR	10 - 200 mg/L Cl ₂