

Lovibond® Colour Measurement Tintometer® Group



Beta-Carotin

Bestell-Nr.: 169520



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

Beta-Carotin

A glass conformance filter to check the Model Fx range of instruments accuracy in measuring beta Carotene.

Technische Daten

Glass Standard Values – 45 ppm bei 10 mm Pfadlänge

Filters 1

Lieferumfang

- 45ppm filter
- Certificate of Conformity