

Lovibond® Colour Measurement Tintometer® Group



Lovibond® RYBN-Farbzahl 2,0R 7,0Y 0,5N (5¼)

Bestell-Nr.: 134100



Reference liquid standards are assigned a nominal value and Lovibond Tintometer always endeavours to match this value as closely as possible. On occasion, it may not be possible to match the requested value exactly. The value achieved and the expected performance tolerances will, of course, always be reported.

Lovibond® RYBN-Farbzahl 2,0R 7,0Y 0,5N (5¼)

500-ml-Flasche, zertifiziert gemäß QS-System nach ISO 9001 mit einer Haltbarkeit von 12 Monaten.

Vollständige Nachverfolgbarkeit gemäß international anerkannten Standards

Jede Flasche wird mit vollständigen Zertifikaten geliefert, inkl. Sicherheitsdatenblatt.