

# Lovibond® Colour Measurement Tintometer® Group



## Flüssigstandard Platin-Cobalt-Farbzahl 5 - 500 ml



Bestell-Nr.: 134140

Reference liquid standards are assigned a nominal value and Lovibond Tintometer always endeavours to match this value as closely as possible. On occasion, it may not be possible to match the requested value exactly. The value achieved and the expected performance tolerances will, of course, always be reported.

### Flüssigstandard Platin-Cobalt-Farbzahl 5 - 500 ml

500-ml-Flasche, zertifiziert gemäß QS-System nach ISO 9001 mit einer Haltbarkeit von 12 Monaten.

Vollständige Nachverfolgbarkeit gemäß international anerkannten Standards

Jede Flasche wird mit vollständigen Zertifikaten geliefert, inkl. Sicherheitsdatenblatt.