

# Lovibond® Water Testing Tintometer® Group



## Legionella Single Swab Test Kit

The World's Fastest *Legionella* Test



- Fast, easy and accurate
- Results in as little as 25 minutes!
- Confirms performance of *Legionella* prevention and treatment protocols
- Test Kits available for detecting both waterborne and biofilm *Legionella*

Part Number: 56B006108

*Legionella* is the causative agent of Legionnaire's disease, a potentially fatal form of pneumonia. Guidelines exist across the world for the control of *Legionella*.

In 2002, the Health and Safety Executive released guidelines for the control of *Legionella*, known as the L8 guidelines. These guidelines detail the key parameters that should be tested on a weekly and monthly basis to minimise the risk of *Legionella* proliferating to harmful levels in open water system.

The Lovibond® *Legionella* "L8" Test Kits contain the recommended tests for both weekly and monthly monitoring.

### *Legionella* Single Swab Test Kit

Within a water system *Legionella* usually grow within biofilm (slime) found on the insides of pipe work, shower heads, cooling tower packing, or water tanks. Bacteria found in the planktonic phase (i.e. suspended in the water itself) are usually an indicator of colonization of the biofilm.

Since a clean water sample does not equate to a clean system (proliferation of *Legionella* is greater in a biofilm than in water), this kit can be used to pinpoint the contamination in a system or to demonstrate cleanliness of an area after a targeted cleaning.

Use this kit to find *Legionella* in their breeding ground -- this very selective test can detect legionella even in samples with very high background flora. With results in about 25 minutes, the test allows for immediate understanding and response to the water system conditions. Regular use will increase confidence in your *Legionella* control procedures and reduce risk of outbreaks and litigation!

**Limit of Detection**  
200 CFU/l / Swabbed area

**No. of Tests**  
1

# Delivery Scope

- Test Strips with Dispensing Pipettes
- Swabs
- Reaction tubes
- Instruction Manual

## Tintometer GmbH

Lovibond® Water Testing  
Schleefstraße 8-12  
44287 Dortmund  
Tel.: +49 (0)231/94510-0  
sales@lovibond.com  
www.lovibond.com  
Germany

## The Tintometer Limited

Lovibond House  
Sun Rise Way  
Amesbury, SP4 7GR  
Tel.: +44 (0)1980 664800  
sales@lovibond.uk  
www.lovibond.com  
UK

## Tintometer China

9F, SOHO II C.  
No.9 Guanghualu,  
Chaoyang District,  
Beijing, 100020  
Customer Care China Tel.: 4009021628  
Tel.: +86 10 85251111 Ext. 330  
Fax: +86 10 85251001  
chinaoffice@tintometer.com  
www.lovibond.com  
China

## Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,  
Lebu Nilam 2, Bandar Bukit Tinggi,  
Klang, 41200, Selangor D.E  
Tel.: +60 (0)3 3325 2285/6  
Fax: +60 (0)3 3325 2287  
lovibond.asia@tintometer.com  
www.lovibond.com  
Malaysia

## Tintometer Brazil

Caixa Postal: 271  
CEP: 13201-970  
Jundiaí – SP  
Tel.: +55 (11) 3230-6410  
sales@lovibond.us  
www.lovibond.com.br  
Brazil

## Tintometer Inc.

6456 Parkland Drive  
Sarasota, FL 34243  
Tel: 941.756.6410  
Fax: 941.727.9654  
sales@lovibond.us  
www.lovibond.us  
USA

## Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2<sup>nd</sup>, 3<sup>rd</sup> & 4<sup>th</sup> Floor  
Sanathnagar Industrial Estate,  
Hyderabad, 500018  
Telangana  
Tel: +91 (0) 40 23883300  
Toll Free: 1 800 599 3891/ 3892  
indiaoffice@lovibond.in  
www.lovibondwater.in  
India

## Tintometer Spain

Postbox: 24047  
08080 Barcelona  
Tel.: +34 661 606 770  
sales@tintometer.es  
www.lovibond.com  
Spain