

# Lovibond® Water Testing

## Tintometer® Group



### *Legionella* Industrial Test Kit - Refill Pack

The World's Fastest *Legionella* Test



*Legionella* is the causative agent of Legionnaire's disease, a potentially fatal form of pneumonia. Guidelines exist across the world for the control of *Legionella*.

In 2002, the Health and Safety Executive released guidelines for the control of *Legionella*, known as the L8 guidelines. These guidelines detail the key parameters that should be tested on a weekly and monthly basis to minimise the risk of *Legionella* proliferating to harmful levels in open water system.

The Lovibond® *Legionella* "L8" Test Kits contain the recommended tests for both weekly and monthly monitoring.

- Fast, easy and accurate
- Results in as little as 30 minutes!
- Confirms performance of *Legionella* prevention and treatment protocols
- Test Kits available for detecting both waterborne and biofilm *Legionella*

Part Number: 56B006104

#### Industry

Chemical Industry | Food and Beverage Industry | Industries Others | Marine Industry | Municipalities | NGO | Oil and Gas | Pharmaceutical Industry | Power and Energy

#### *Legionella* Industrial Test Kit - Refill Pack

The *Legionella* System Test Kit is for detecting *Legionella* in a water sample — it is ideal for on-line and a variety of industrial applications.

This test procedure utilizes hollow fiber filtration modules which attach directly to the provided sample point fittings (both a 1/2" and 3/4" female threaded fittings are provided). The filtered sample increases the sensitivity of the test.

This filtration step requires no pumping — just connect the sample point fitting to the water system (such as a shower head or on-line sample point) and use the system pressure to drive the filtration. The pre-filled syringes contain the required recovery buffer, so there is no need to measure or dilute solutions.

With results in 30-35 minutes (including sample prep-time), the test allows for immediate understanding and response to the water system conditions. Regular use will increase confidence in your *Legionella* control procedures and reduce risk of outbreaks and litigation!

#### Limit of Detection

100 CFU/l

#### No. of Tests

5

# Delivery Scope

- 5 Test Strips with Dispensing Pipettes
- 5 Syringes with Recovery Buffer
- Sample Collection Jar
- 5 single-use filters
- Female threaded sample point fitting
- Instruction Manual

## Tintometer GmbH

Lovibond® Water Testing  
Schleefstraße 8-12  
44287 Dortmund  
Tel.: +49 (0)231/94510-0  
sales@lovibond.com  
www.lovibond.com  
Germany

## The Tintometer Limited

Lovibond House  
Sun Rise Way  
Amesbury, SP4 7GR  
Tel.: +44 (0)1980 664800  
Fax: +44 (0)1980 625412  
sales@lovibond.uk  
www.lovibond.com  
UK

## Tintometer China

9F, SOHO II C.  
No.9 Guanghualu,  
Chaoyang District,  
Beijing, 100020  
Customer Care China Tel.: 4009021628  
Tel.: +86 10 85251111 Ext. 330  
Fax: +86 10 85251001  
chinaoffice@tintometer.com  
www.lovibond.com  
China

## Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,  
Lebu Nilam 2, Bandar Bukit Tinggi,  
Klang, 41200, Selangor D.E  
Tel.: +60 (0)3 3325 2285/6  
Fax: +60 (0)3 3325 2287  
lovibond.asia@tintometer.com  
www.lovibond.com  
Malaysia

## Tintometer Brazil

Caixa Postal: 271  
CEP: 13201-970  
Jundiaí – SP  
Tel.: +55 (11) 3230-6410  
sales@lovibond.us  
www.lovibond.com.br  
Brazil

## Tintometer Inc.

6456 Parkland Drive  
Sarasota, FL 34243  
Tel: 941.756.6410  
Fax: 941.727.9654  
sales@lovibond.us  
www.lovibond.us  
USA

## Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2<sup>nd</sup>, 3<sup>rd</sup> & 4<sup>th</sup> Floor  
Sanathnagar Industrial Estate,  
Hyderabad, 500018  
Telangana  
Tel: +91 (0) 40 23883300  
Toll Free: 1 800 599 3891/ 3892  
indiaoffice@lovibond.in  
www.lovibondwater.in  
India

## Tintometer Spain

Postbox: 24047  
08080 Barcelona  
Tel.: +34 661 606 770  
sales@tintometer.es  
www.lovibond.com  
Spain