

Lovibond® Water Testing Tintometer® Group



MD 200 Chlorine, pH, Cyanuric acid, Alkalinity-m, tablet reagents (OTZ)

Precise results using high-quality interference filters



- Scroll memory
- Automatic switch-off
- Real-Time-Clock and date
- Calibration mode indicator

Part Number: 2860502

Versatile use

The MD 200 photometer can be used in almost all areas of water analysis.

Precise and low maintenance

The high-precision optics, which in addition to long-term stable LEDs as light sources, also offer high-quality interference filters, ensure precise and reproducible measurement results. In addition, these components make the photometer low-maintenance.

User friendly

The devices impress with their ease of use, ergonomic design and compact dimensions, which makes handling the devices easy and safe.

Reagents

Use long-term stable Lovibond® reagent tablets (5-10 years guaranteed minimum durability) or Lovibond® liquid reagents.

Industry

Marine Industry | Municipalities | Pool Domestic | Pool Public | Pool Services

Application

Boiler Water | Cooling Water | Disinfection Control | Drinking Water Treatment | Pool Water Control | Raw Water Treatment | Waste-Water Treatment

MD 200 Chlorine, pH, Cyanuric acid, Alkalinity-m, tablet reagents (OTZ)

Versatile, compact and: The benchtop photometer MD 200 combines high-quality interference filters with long-term stable LEDs as light sources in a high-precision optic. The result: precise measurement results at all times - no matter where you use it.

Measuring Range

| Test Name | Measuring Range | Chemical Method |
|----------------|---|------------------|
| Alkalinity-m T | 5 - 200 mg/L CaCO ₃ | Acid / Indicator |
| Chlorine HR T | 0.1 - 10 mg/L Cl ₂ ^{a)} | DPD |
| Chlorine L | 0.02 - 4.0 mg/L Cl ₂ ^{a)} | DPD |
| Chlorine T | 0.01 - 6.0 mg/L Cl ₂ ^{a)} | DPD |
| CyA T | 10 - 160 mg/L CyA | Melamine |
| pH value L | 6.5 - 8.4 pH | Phenol Red |
| pH-value T | 6.5 - 8.4 pH | Phenol Red |

Technical Data

| | |
|---------------------------------|---|
| Optics | LEDs, interference filters (IF) and photo sensor in transparent sample chamber. Depending on the version, up to 3 different interference filters are used. Wavelength specifications of interference filters: 430 nm $\Delta\lambda = 5$ nm 530 nm $\Delta\lambda = 5$ nm 560 nm $\Delta\lambda = 5$ nm 610 nm $\Delta\lambda = 6$ nm |
| Wavelength Accuracy | ± 1 nm |
| Photometric Range | -2500 - 2500 mAbs |
| Photometric Accuracy | 3 % FS (T = 20 °C – 25 °C) |
| Display | backlit LCD (on keypress) |
| Interfaces | Infrared |
| Operation | Foil keypad |
| Auto – OFF | Yes |
| Calibration | Factory calibration and user calibration. Reset to factory calibration possible. |
| Internal Storage | internal ring memory for 16 data sets |
| Power Supply | 4 batteries (AA) |
| Battery Life Time | approx. 53 hours |
| Clock | Real Time Clock and Date |
| Portability | Portable |
| Environmental Conditions | temperature: 5 - 40 °C rel. humidity: 30 - 90 % (non condensing) |
| Protection Class | IP 68 |
| Compliance | CE |
| Dimensions | 110 x 55 x 190 mm |
| Weight | 455 g |

Delivery Scope

- Instrument in plastic case with 4 batteries (AA)
- 3 Round vials with lids
- Sealing rings
- 1 Stirring rod
- 1 Brush
- 1 Syringe (5 ml)
- 100 tablet reagents DPD No.1
- 100 tablet reagents DPD No.3 Evo
- 100 tablet reagents Phenol Red Photometer
- 100 tablet reagents CyA-Test
- 100 tablet reagents Alka-M Photometer
- Instruction manual
- Manual of Methods
- Certificate (COC)
- Warranty information

| Title | Part Number |
|--|-------------|
| Water sampler with bottle 250 mL and lid, AF 631 | 170500 |
| Batteries (AA), set of 4 | 1950025 |
| Round cuvette 24 mm \varnothing , set of 12 | 197620 |
| Sealing ring for round vials 24 mm \varnothing , set of 12 | 197626 |
| Round cuvette 24 mm \varnothing , set of 5 | 197629 |
| Cleaning cloth | 197635 |
| Round cuvette 16 mm \varnothing , set of 10 | 197665 |
| Adapter for round cuvettes 16 mm \varnothing | 19802190 |
| Battery lid | 19802241 |
| Fixed price service package for MD100/MD110/MD200 | 19802701 |
| Service plan - 3 years for MD100/MD110/MD200 | 19802801 |
| Verification Standard Kit MD 100 / 200 | 215670 |
| Reference Standard Kit Chlorine - 0.2 and 1.0 mg/l (MD 100/110/200) | 275650 |
| Reference Standard Kit Chlorine - 0.5 and 2.0 mg/l (MD 100/110/200) | 275655 |
| Reference Standard Kit Chlorine - 1.0 and 4.0 mg/l (MD 100/110/200) | 275656 |
| Reference Standard Kit pH - 7.45 (MD 100/110/200) | 275670 |
| Stirring rod, 13 cm length | 364100 |
| Stirring rod, 10 cm length | 364109 |
| Stirring rod, 13 cm length, set of 10 | 364120 |
| Stirring rod, 10 cm length, set of 10 | 364130 |
| Membrane filter set for use when preparing samples, 25 membrane filters 0.45 μ m, 2 syringes 20 ml | 366150 |
| 250 mL bottle, AF 631 | 375072 |
| Brush, 11 cm length | 380230 |
| Measuring beaker, 100 ml | 384801 |
| Cuvette stand for 6 round cuvettes \varnothing 24 mm | 418951 |
| Cuvette stand for 10 round cuvettes \varnothing 16 mm | 418957 |
| ValidCheck Chlorine 1,5 mg/l | 48105510 |
| Factory calibration certificate ISO 9001 for MD100/MD110/MD200 | 999750 |
| Additionally for each method | 999855 |

Tintometer GmbH

Lovibond® Water Testing
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
sales@lovibond.com
www.lovibond.com
Germany

The Tintometer Limited

Lovibond House
Sun Rise Way
Amesbury, SP4 7GR
Tel.: +44 (0)1980 664800
Fax: +44 (0)1980 625412
sales@lovibond.uk
www.lovibond.com
UK

Tintometer China

9F, SOHO II C.
No.9 Guanghualu,
Chaoyang District,
Beijing, 100020
Customer Care China Tel.: 4009021628
Tel.: +86 10 85251111 Ext. 330
Fax: +86 10 85251001
chinaoffice@tintometer.com
www.lovibond.com
China

Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,
Lebuh Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com
www.lovibond.com
Malaysia

Tintometer Brazil

Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@lovibond.us
www.lovibond.com.br
Brazil

Tintometer Inc.

6456 Parkland Drive
Sarasota, FL 34243
Tel: 941.756.6410
Fax: 941.727.9654
sales@lovibond.us
www.lovibond.us
USA

Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2nd, 3rd & 4th Floor
Sanathnagar Industrial Estate,
Hyderabad, 500018
Telangana
Tel: +91 (0) 40 23883300
Toll Free: 1 800 599 3891/ 3892
indiaoffice@lovibond.in
www.lovibondwater.in
India

Tintometer Spain

Postbox: 24047
08080 Barcelona
Tel.: +34 661 606 770
sales@tintometer.es
www.lovibond.com
Spain