

# Lovibond® Colour Measurement Tintometer® Group



## LC 100 SV 100 Kit

Hand-Held Imaging Spectrocolorimeter With Viewing Stand



- Handheld and Benchtop Spectrocolorimeter providing Fast, Precise and Accurate Colour Measurement
- Up to 350 Sample Measurements can be Stored and Compared to Reference Standards
- Conformance Reports Available in PDF and Excel Formats
- Precise Positioning on Sample with Preview of Measurement Area

Part Number: 401120

### Handheld and Benchtop

The LC 100/SV 100 offers the ultimate in flexibility. Its compact and lightweight design along with its uniquely designed viewing stand ensures the instrument can be used anywhere within a busy manufacturing plant or laboratory. Fitting in the palm of your hand the lightweight LC 100 unit offers easy one-handed operation for prolonged use. While the SV 100 viewing stand enables awkward samples to be measured using Lovibond® optical cells, giving fast and accurate results.

### Accuracy

Despite the lightweight and compact design, the LC 100/SV 100 does not compromise on accuracy. Using 8 different visible illuminations and 1 UV LED along with the 45°/0 optical geometry and image capture technology, the user can feel confident in the results they are receiving.

### Reporting Results

To make reporting as quick and easy as possible, the LC 100/SV 100 automatically creates conformance reports in either PDF or Excel formats. The user can then share these with a PC or directly to a compatible printer.

### Storage

The LC 100/SV100 can store up to 350 sample measurements, these samples can then be tagged via the LC 100/SV 100's on-screen keyboard or built in voice recorder. The LC 100/SV 100 can also store up to 30 reference standards by date and time, ensuring a high level of quality control when comparing results.

### Suitable for Multiple Applications

The versatility of the LC 100/SV 100 means its suitable for multiple applications and has the ability to report with a display of  $\Delta E$  for all com-

mon methods i.e. CIELAB, CMC, CIE 94, 2000, plus a display of measured and delta differences with  $L^*a^*b^*C^*h^*$  for standard and sample, graphical  $L^*a^*b^*$  plot.

## Industry

Chemical Industry | Food and Beverage Industry | Industries Others | Oil and Gas | Pharmaceutical Industry

## Application

Chemicals | Coatings | Food and Beverage | Household Products | Medical and Clinical Products | Pharmaceuticals and Cosmetics | Plastics | Sugar Solutions, Syrups and Honey | Textiles

## LC 100 SV 100 Kit

The LC 100/SV 100 is a highly flexible and lightweight spectrophotometer. This handheld and benchtop unit offers the user the ability to capture fast and accurate colour measurement throughout the manufacturing process. In just over a second, the LC 100/SV 100 will provide a sample reading along with a comprehensive report of how close the colour is to a self-programmed standard.

## Technical Data

<b>Applicable Standard</b>	ASTM E308
<b>Colour Type</b>	Transparent Translucent Opaque
<b>Display</b>	4.5 cm Colour TFT
<b>Interfaces</b>	USB
<b>Measuring Geometrics</b>	45/0 Image Capture
<b>Internal Storage</b>	30 standards/350 sample results
<b>Measurement Time</b>	1.8 seconds
<b>Storage Temperature</b>	-20° to 60°C (-4°F to 104°F)
<b>Short Term Repeatability</b>	Typical 0.10 $\Delta E$ 94 on white (D65/10)
<b>Illuminants</b>	CIE Illuminant A CIE Illuminant D65
<b>Portability</b>	Portable

## Delivery Scope

- LC 100 Handheld Spectrophotometer
- SV 100 Viewing Stand
- Carrying Pouch
- 10 x Plastic Cells (10ml)
- Cell Adaptor
- UBS Cable for Power Supply and Data Transfer

### Tintometer GmbH

Lovibond® Water Testing  
Schleefstraße 8-12  
44287 Dortmund  
Tel.: +49 (0)231/94510-0  
sales@lovibond.com  
www.lovibond.com  
Germany

### The Tintometer Limited

Lovibond House  
Sun Rise Way  
Amesbury, SP4 7GR  
Tel.: +44 (0)1980 664800  
Fax: +44 (0)1980 625412  
sales@lovibond.uk  
www.lovibond.com  
UK

### Tintometer China

9F, SOHO II C.  
No.9 Guanghualu,  
Chaoyang District,  
Beijing, 100020  
Customer Care China Tel.: 4009021628  
Tel.: +86 10 85251111 Ext. 330  
Fax: +86 10 85251001  
chinaoffice@tintometer.com  
www.lovibond.com  
China

### Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,  
Lebuh Nilam 2, Bandar Bukit Tinggi,  
Klang, 41200, Selangor D.E  
Tel.: +60 (0)3 3325 2285/6  
Fax: +60 (0)3 3325 2287  
lovibond.asia@tintometer.com  
www.lovibond.com  
Malaysia

### Tintometer Brazil

Caixa Postal: 271  
CEP: 13201-970  
Jundiaí – SP  
Tel.: +55 (11) 3230-6410  
sales@lovibond.us  
www.lovibond.com.br  
Brazil

### Tintometer Inc.

6456 Parkland Drive  
Sarasota, FL 34243  
Tel: 941.756.6410  
Fax: 941.727.9654  
sales@lovibond.us  
www.lovibond.us  
USA

### Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2<sup>nd</sup>, 3<sup>rd</sup> & 4<sup>th</sup> Floor  
Sanathnagar Industrial Estate,  
Hyderabad, 500018  
Telangana  
Tel: +91 (0) 40 23883300  
Toll Free: 1 800 599 3891/ 3892  
indiaoffice@lovibond.in  
www.lovibondwater.in  
India

### Tintometer Spain

Postbox: 24047  
08080 Barcelona  
Tel.: +34 661 606 770  
sales@tintometer.es  
www.lovibond.com  
Spain

Technical changes without notice

Printed in Germany

Lovibond® and Tintometer® are Trademarks of the Tintometer Group of Companies