

# Lovibond® Water Testing Tintometer® Group



## pH Strips pH 7-14 Plastic

For fast determination with chemical reaction



- Hazard free
- Rapid
- Always to hand
- Semi quantitative results

Part Number: 56S001190

### Industry

Industries Others

### Application

Cooling Water | Disinfection Control

### pH Strips pH 7-14 Plastic

Test strips offer quick, clean results without the fuss of reagents or tablets. Their dip and compare method makes them a popular choice for many applications where speed and simplicity are most important.

### Measuring Range

Test Name	Measuring Range
pH	7 - 14 pH