

# Lovibond® Water Testing

## Tintometer® Group



### 5in1 Mini

The focus on various parameters simultaneously



- Easy to use
- Futuristic design
- RAPID tablets fast dissolving
- Highest accuracy

Part Number: 151700

#### Inexpensive

The Pool Test Kit is the cost-effective entry-level variant in water analysis.

#### Variable

The different versions of the Pool Test Kit allow you to test different parameters. Possible combinations are, for example, chlorine + pH, bromine + pH, active oxygen + pH or biguanide (PHMB) + pH.

#### Ready for immediate use

The Pool Test Kit includes the tablet reagents for 20 tests each and instruction manual, all packaged in a handy and transparent mini suitcase. So you not only keep track of everything, but also have the test kit at your fingertips at any time.

#### Industry

Pool Domestic | Pool Services

#### Application

Disaster Management | Pool Water Control

#### 5in1 Mini

Perfect for the entry into the water analysis in the private environment. Measure several parameters at the same time.

#### Measuring Range

| Test Name     | Measuring Range                 |
|---------------|---------------------------------|
| Alkalinity-m  | 50 - 300 mg/L CaCO <sub>3</sub> |
| Bromine       | 0.2 - 6.8 mg/L Br <sub>2</sub>  |
| Chlorine      | 0.1 - 3.0 mg/L Cl <sub>2</sub>  |
| Cyanuric acid | 20 - 200 mg/L CyA               |
| pH            | 6.8 - 8.2 pH                    |

# Delivery Scope

- Pool Test Kit in a bubble pack or mini case
- Instruction manual

## **Tintometer GmbH**

Lovibond® Water Testing  
Schleefstraße 8-12  
44287 Dortmund  
Tel.: +49 (0)231/94510-0  
sales@lovibond.com  
www.lovibond.com  
Germany

## **The Tintometer Limited**

Lovibond House  
Sun Rise Way  
Amesbury, SP4 7GR  
Tel.: +44 (0)1980 664800  
sales@lovibond.uk  
www.lovibond.com  
UK

## **Tintometer China**

9F, SOHO II C.  
No.9 Guanghualu,  
Chaoyang District,  
Beijing, 100020  
Customer Care China Tel.: 4009021628  
Tel.: +86 10 85251111 Ext. 330  
Fax: +86 10 85251001  
chinaoffice@tintometer.com  
www.lovibond.com  
China

## **Tintometer South East Asia**

Unit B-3-12, BBT One Boulevard,  
Lebuh Nilam 2, Bandar Bukit Tinggi,  
Klang, 41200, Selangor D.E  
Tel.: +60 (0)3 3325 2285/6  
Fax: +60 (0)3 3325 2287  
lovibond.asia@tintometer.com  
www.lovibond.com  
Malaysia

## **Tintometer Brazil**

Caixa Postal: 271  
CEP: 13201-970  
Jundiaí – SP  
Tel.: +55 (11) 3230-6410  
sales@lovibond.us  
www.lovibond.com.br  
Brazil

## **Tintometer Inc.**

6456 Parkland Drive  
Sarasota, FL 34243  
Tel: 941.756.6410  
Fax: 941.727.9654  
sales@lovibond.us  
www.lovibond.us  
USA

## **Tintometer India Pvt. Ltd.**

Door No: 7-2-C-14, 2<sup>nd</sup>, 3<sup>rd</sup> & 4<sup>th</sup> Floor  
Sanathnagar Industrial Estate,  
Hyderabad, 500018  
Telangana  
Tel: +91 (0) 40 23883300  
Toll Free: 1 800 599 3891/ 3892  
indiaoffice@lovibond.in  
www.lovibondwater.in  
India

## **Tintometer Spain**

Postbox: 24047  
08080 Barcelona  
Tel.: +34 661 606 770  
sales@tintometer.es  
www.lovibond.com  
Spain