

# Lovibond® Water Testing Tintometer® Group



## AF103 Chlor HR Kit 3/2 Iod

More than 400 different test discs available



- More than 400 different test discs available
- Compensation for coloured and turbid samples
- Guaranteed constancy of the coloured glass standards
- Integrated prism

Part Number: 56K014401

### Usable in different applications

The Comparator 2000+ is versatile: it is used in swimming pools as well as in research centers, universities or general drinking water treatment.

### Daylight lamp

So that you can always measure safely and stably under unfavorable visibility conditions, we recommend the use of the (optional) mains or battery-powered daylight lamp. So the lighting conditions always remain the same. No matter where you are measuring.

### Various analyses possible

For the Comparator 2000+, more than 400 color discs and numerous reagents are available. Your combination is guaranteed. And if you have to change them once, it works without opening the comparator.

### Integrated prism

The integrated prism visually brings the glass standards of the color disc and the colored sample into view.

### Industry

Chemical Industry | Food and Beverage Industry | Industries Others | Oil and Gas | Pharmaceutical Industry | Power and Energy

### AF103 Chlor HR Kit 3/2 Iod

The Comparator 2000+ is a versatile colorimetric system for water analysis. Comfortable handling, no compromise on accuracy and reproducibility: The Comparator 2000+ can be used for many applications. Swimming pools, research or drinking water treatment are just a few examples.

### Measuring Range

Test Name	Measuring Range	Chemical Method
Chlorine HR T	5.0 - 250 mg/L Cl <sub>2</sub>	Potassium Iodide