

# Lovibond® Water Testing Tintometer® Group



## HARDNESS-total TEST KIT (°dH)

For the standard situation in water testing



- Easy to use
- Futuristic design
- RAPID tablets fast dissolving
- Highest accuracy

Part Number: 2418411

### HARDNESS-total TEST KIT (°dH)

These test kits are specially developed for testing boiler, cooling and industrial process water. They make use of both colorimetric and titrimetric techniques. Each test kit contains all the necessary chemicals and reagents in liquid or powder form to conduct the tests. The detailed instructions contain a step-by-step explanation of the test procedure. The kits are supplied in a sturdy, compact plastic case. Keenly priced refill reagent packs are available for all test kits.

### Measuring Range

Test Name	Measuring Range
Hardness carbonate (new version)	1 drop = 1 or 0.5 °dH

### Technical Data

Dimensions	90 x 36 x 36 mm
Weight with Packaging	30 g

Title	Part Number
Hardness-Total-Reagent	418563