

Lovibond® Water Testing Tintometer® Group



CHLORIDE TEST KIT CD-2

Determination of Boiler-, Cooling- and Industrial Process Water



- Fast quantitative determination
- For testing boiler, cooling and industrial process water
- Suitable for field and laboratory testing
- Cost-effective use due to keenly priced re-fill packs

Part Number: 2418506

CHLORIDE TEST KIT CD-2

These test kits are specially developed for testing boiler, cooling and industrial process water. They make use of both colorimetric and titrimetric techniques. Each test kit contains all the necessary chemicals and reagents in liquid or powder form to conduct the tests. The detailed instructions contain a step-by-step explanation of the test procedure. The kits are supplied in a sturdy, compact plastic case. Keenly priced refill reagent packs are available for all test kits.

Measuring Range

Test Name	Measuring Range
Chloride HR	1 drop = 50 or 25 Cl ⁻

Technical Data

Dimensions	157 x 135 x 55 mm
Weight with Packaging	185 g

Title	Part Number
Chloride-Refill Reagents for CD-2	2418507