

# Lovibond<sup>®</sup> Colour Measurement Tintometer<sup>®</sup> Group



## Lovibond<sup>®</sup> RYBN Colour 0.8R 2.0Y 0.1N (5¼)

Part Number: 134080



Reference liquid standards are assigned a nominal value and Lovibond Tintometer always endeavours to match this value as closely as possible. On occasion, it may not be possible to match the requested value exactly. The value achieved and the expected performance tolerances will, of course, always be reported.

### Lovibond<sup>®</sup> RYBN Colour 0.8R 2.0Y 0.1N (5¼)

500ml Bottle certified under ISO 9001 quality system with a 12-month shelf life.

Full traceability to internationally recognised standards

Each bottle supplied with full certification including MSDS