

Lovibond® Colour Measurement Tintometer® Group



Juego de filtros EP marrón/amarillo

Referencia No: 139440



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

Juego de filtros EP marrón/amarillo

A set of two glass conformance filters for checking PFXi instruments in European Pharmacopoeia Brown/Yellow scale.

Datos técnicos

Glass Standard Values	- BY3
	- BY5

Filters	2
----------------	---

Volumen de suministro

- EP Brown/Yellow 3 filter
- EP Brown/Yellow 5 filter
- Certificate of Conformity