

# Lovibond® Water Testing Tintometer® Group



## Hydrogen Peroxide Strips

For fast determination with chemical reaction



- Hazard free
- Rapid
- Always to hand
- Semi quantitative results

Part Number: 56S000850

### Industry

Industries Others

### Application

Cooling Water | Disinfection Control

### Hydrogen Peroxide Strips

Test strips offer quick, clean results without the fuss of reagents or tablets. Their dip and compare method makes them a popular choice for many applications where speed and simplicity are most important.

### Measuring Range

**Test Name**

Hydrogen Peroxide test strips

**Measuring Range**

0.5 - 100 mg/L H<sub>2</sub>O<sub>2</sub>