

Lovibond® Water Testing Tintometer® Group



Hardness Total Strips

For fast determination with chemical reaction



- Hazard free
- Rapid
- Always to hand
- Semi quantitative results

Part Number: 56S000950

Industry

Industries Others

Application

Boiler Water | Cooling Water

Hardness Total Strips

Test strips offer quick, clean results without the fuss of reagents or tablets. Their dip and compare method makes them a popular choice for many applications where speed and simplicity are most important.

Measuring Range

Test Name	Measuring Range
Hardness Total-Test Strips	40 - 1000 mg/L CaCO ₃