

Lovibond® Water Testing Tintometer® Group



pH Strips pH 4.5-10 Plastic

For fast determination with chemical reaction



- Hazard free
- Rapid
- Always to hand
- Semi quantitative results

Part Number: 56S001090

Industry

Industries Others

Application

Cooling Water

pH Strips pH 4.5-10 Plastic

Test strips offer quick, clean results without the fuss of reagents or tablets. Their dip and compare method makes them a popular choice for many applications where speed and simplicity are most important.

Measuring Range

Test Name	Measuring Range
pH	4.5 - 10 pH