

# Lovibond® Water Testing Tintometer® Group



## TEST KIT TOTAL- AND HARDNESS CARBONATE

For the standard situation in water testing



- Rapid quantitative analysis
- Ideal for boiler, cooling and process water
- Field and laboratory tests possible
- Economical due to inexpensive refill units

Part Number: 2418412

### TEST KIT TOTAL- AND HARDNESS CARBONATE

These test kits are specially developed for testing boiler, cooling and industrial process water. They make use of both colorimetric and titrimetric techniques. Each test kit contains all the necessary chemicals and reagents in liquid or powder form to conduct the tests. The detailed instructions contain a step-by-step explanation of the test procedure. The kits are supplied in a sturdy, compact plastic case. Keenly priced refill reagent packs are available for all test kits.

#### Measuring Range

Test Name	Measuring Range
Hardness Total (new version) + Carbonate GKH-1	1 drop = 1 or 0.5°dH

#### Technical Data

Dimensions	90 x 35 x 70 mm
Weight with Packaging	80 g