

# Lovibond® Water Testing Tintometer® Group



## total-PHOSPHATE TEST KIT PO-2

Determination of Boiler-, Cooling- and Industrial Process Water



- Fast quantitative determination
- For testing boiler, cooling and industrial process water
- Suitable for field and laboratory testing
- Cost-effective use due to keenly priced re-fill packs

Part Number: 2418523

### total-PHOSPHATE TEST KIT PO-2

These test kits are specially developed for testing boiler, cooling and industrial process water. They make use of both colorimetric and titrimetric techniques. Each test kit contains all the necessary chemicals and reagents in liquid or powder form to conduct the tests. The detailed instructions contain a step-by-step explanation of the test procedure. The kits are supplied in a sturdy, compact plastic case. Keenly priced refill reagent packs are available for all test kits.

#### Note according to EU regulation 2019/1148

This product is regulated by Regulation (EU) 2019/1148: all suspicious transactions, and significant disappearances and thefts should be reported to the relevant national contact point. See [https://ec.europa.eu/info/index\\_en](https://ec.europa.eu/info/index_en)

### Measuring Range

Test Name	Measuring Range
Phosphate (total)	2.5 - 25 mg/L PO <sub>4</sub> <sup>3-</sup>