


Lovibond® Water Testing

Tintometer® Group



Closed System Engineers Test Kit

The expert tool for reliable result



These Test Kits for cooling water monitoring are specially designed for engineers and mobile use. It provides the equipment and reagents required to monitor the correct amount and dosing levels of various important parameters. You can carry everything you need with you and test directly onsite. The Kit is provided in a robust carry case and includes all important equipment from measurement instruments to reagents.

- Complete analysis kit including molybdate, iron and conductivity
- Prevents corrosion and scale
- Accurate inhibitor dosing rates
- Reproducible, accurate results

N. ordine: 56K701600

Industria

Altre industrie | Comuni | Fornitori di energia | Industria alimentare e delle bevande | Industria chimica | Industria farmaceutica | Industria olearia | Navigazione

Applicazione

Acqua di caldaia | Acqua di raffreddamento | Trattamento acqua di scarico

Closed System Engineers Test Kit

Due to its heat absorption capacity, water is the preferred cooling medium in many areas of industry. However, water always contains many dissolved and undissolved substances that can cause trouble even when travelling in closed cooling circuits. Extensive water treatment or chemical post-treatment by experts is always necessary. Problems are mainly caused by deposits, corrosion and microbiological activities. For the experts, a regular check of certain parameters is therefore just as essential as reliable equipment.

Campo di misura

Test Name	Campo di misura	Metodo chimico
Alkalinity M	50 - 2400 mg/L CaCO ₃	
Alluminio T	0 - 0.3 mg/L Al	Eriocromocianina R
Chloride	20 - 12000 mg/L Cl ⁻	
Durezza totale	5 - 600 mg/L CaCO ₃	
Ferro LR T	0.05 - 1 mg/L Fe	PPST
Glicole	% PEG/MEG	
Molybdato HR T	50 - 500 mg/L MoO ₄	Tioglicolato
Nitrito	10 - 2000 mg/L NaNO ₂	
Rame HR T	0 - 5 mg/L Cu	Bichinolina

Fornitura

- Plastic case
- SD40 Multi Pocket Tester
- Checkit Comparator
- 50 - 500 mg/L Molybdate HR+ Checkit Disc
- 0 - 1.0 mg/L Iron LR Checkit Disc
- 0 - 0.3 mg/L Aluminium Checkit Disc
- 0 - 5.0 mg/L Copper HR Free and Total Checkit Disc
- 2 Plastic vials 10 mL with lid
- Aluminium No.1 tablets
- Aluminium No.2 tablets
- Copper No. 1 Test tablets
- Copper No. 2 Test tablets
- Molybdate No.1 tablets
- Molybdate No.2 tablets
- pH Buffer 4 solution
- pH Buffer 7 solution
- pH Buffer 10 solution
- Iron Reagent FE5
- Hardness Total Buffer TH2
- Hardness Total Indicator TH1P
- Hardness LR Titrant TH3
- Hardness HR Titrant TH4
- 1413µS Standard Conductivity Solution
- Alkalinity 4.5 Indicator TA4
- Alkalinity LR Titrant TA3
- Alkalinity HR Titrant PA2/TA2
- Nitrite Indicator N1
- Nitrite Titrant N2
- Chloride Indicator CC1
- Chloride LR Titrant CC2
- Chloride HR Titrant BC2
- 2 Titration jars with caps, plastic, 60 mL
- Syringe plastic, 20 mL
- Wash Bottle, 250 mL
- Syringe Filter Holder 25 mm
- Filter Circle 0.45 µm, 25 mm
- Glycol Refractometer, RHA-218ATC
- Instruction manual

Titolo	N. ordine
CHECKIT® Comparator	145000
Siringa, plastica, 20 mL	56A006501
Vaso di titolazione con coperchio, plastica, 60 mL	56A006701
Siringa di plastica, 5 mL	56A008501
Flacone di lavaggio 250 ml	56A009001
Siringa di plastica 1 ml	56A013501
Cerchio del filtro 0,45 um, 25 mm	56A020050
Provetta 5/10 mL + CAP	56A600401
Rifrattometro per glicole	56M000401

Tintometer GmbH

Lovibond® Water Testing
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
sales@lovibond.com
www.lovibond.com
Germania

The Tintometer Limited

Lovibond House
Sun Rise Way
Amesbury, SP4 7GR
Tel.: +44 (0)1980 664800
sales@lovibond.uk
www.lovibond.com
Regno Unito

Tintometer China

Room 1001, China Life Tower
16 Chaoyangmenwai Avenue,
Beijing, 100020
Customer Care China Tel.: 4009021628
Tel.: +86 10 85251111 App. 330
Fax: +86 10 85251001
chinaoffice@tintometer.com
www.lovibond.com
Cina

Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,
Lebuh Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com
www.lovibond.com
Malesia

Tintometer Brazil

Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@lovibond.us
www.lovibond.com.br
Brasile

Tintometer Inc.

6456 Parkland Drive
Sarasota, FL 34243
Tel: 941.756.6410
Fax: 941.727.9654
sales@lovibond.us
www.lovibond.us
Stati Uniti d'America

Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2nd, 3rd & 4th Floor
Sanathnagar Industrial Estate,
Hyderabad, 500018
Telangana
Tel: +91 (0) 40 23883300
Toll Free: 1 800 599 3891/ 3892
indiaoffice@lovibond.in
www.lovibondwater.in
India

Tintometer Spain

Postbox: 24047
08080 Barcelona
Tel.: +34 661 606 770
sales@tintometer.es
www.lovibond.com
Spagna