

Lovibond® Water Testing Tintometer® Group



Water Safety Kit Basic

Uitgerust voor alle vereisten



Don't wait for solutions in disasters or difficult environments - simply react. With the Lovibond® Water Safety Kit Basic you take everything with you to manage each testing challenge and provide safe potable water. The kit combines microbiological and simple chemical methods of analysis.

Perform microbiological membrane filtration and apply chemical test equipment directly on-site - rapidly & with reliable results.

Its compact design and the optional backpack ensures your mobility in each environment and situation. The Heating and cooling function of the incubator allows you test different sorts of bacteria. Be sure that you don't need any further test kit or test equipment to control the most important parameters you need for safe potable water quality.

- Emergency response kit: Manage each testing challenge
- Portable basic test equipment: Microbiological & chemical analysis in one kit
- Incubator with heating and cooling as standard
- Safe & simple testing according to WHO standards

Bestelnr.: 56K681250

Industrie

Gemeenten | Ngo | Openbaar zwembad | Overige industrieën | Scheepvaart | Voedingsmiddelen- en drankenindustrie

Applikation

Behandeling drinkwater | Controle desinfecteermiddel | Koelwater | Rampenbeheer | Voedingsmiddelen en dranken

Water Safety Kit Basic

Whether disaster or monitoring in difficult environments: You need mobile testing equipment for portable water in each situation. the Water Safety Kit Basic accompanies you with the complete test methods for the most important parameters of both microbiological and physio chemical analysis.

Profit from the advantages of rugged, waterproof cases with compact design and easy to use content.

User

Non Profit Organizations (NGO's) | Medical research and development | Institutes, Universities, Schools | World Health Organization (WHO)

Meetbereik

Test Name	Meetbereik
Aërobe bacteriën (TVC)/colibacteriën	
Aërobe bacteriën (TVC)/E.Coli, coliformen	
Chlorine	0.1 - 3.0 mg/L Cl ₂
Conductivity	0 - 20 mS/cm
pH	6.8 - 8.2 pH
Turbidity	30 - 400 NTU

Technische gegevens

Afmetingen 487 x 229 x 386 mm

Gewicht met verpakking 16 kg

Omschrijving	Bestelnr.
Stroomkabelset	136300
Incubator Battery, Lead/Acid	192345
Flesjes	363590
Universal container with spoon and lid, 30 ml	366239
Plastic beker, gegraduateerd, 100 ml, met deksel	385130
Petrischaaltjes, 10 stuks	400855
Troebelheidstest 5-500 JTU	56A010801
Rugzak voor Water Safety Kits	56A01400
Media Pads and Dispenser (Pack of 200)	56A016330
Maatbeker	56A018501
Membraanfiltratie	56A023801
Pincet	56A604201
Turbidity tube (one piece) - 30-400 NTU	56A607101
Membraanfiltratie	56B000710
D004 TTC/MAC Dipslides (Macconkey)	56B010410
D005 TTC/Chromogeen Dipslides	56B010510
KS16-1413US Standaardgeleidbaarheidsoplossing	56L001665
Membrane Lauryl Sulphate Broth, 38.1 gr.	56P069740
Pump (35 ml stroke)	56X610301
Pipette, disposable, 1.5 ml	BI-K10060
Marker Pen	BI-K14087

Leveromvang

- DI 20 Incubator
- Incubator battery
- 50 Petri dish for DI 20 incubator
- SD 70 Conductivity Meter
- Marker pen
- 2.5 " Bailer with rope
- Filter flask, 500 ml
- Tweezer
- Pump (35 ml)
- Universal container and spoon, 30 ml
- DB 513 Shaker tube with lid, 100 ml
- 150 EZ Pack filters
- Pipette disposable 1.5 ml (EDT222)
- Turbidity Tube
- KS16 1413 µS Conductivity standard solution (65 ml)
- 200 Media pads and dispenser
- Filtration device
- Lead set IEC C13 (3 pin)
- Membrane Lauryl Sulphate Broth (40 g)
- Chlorine-pH Tester
- 250 DPD No.1 Photometer tablet reagents
- 250 DPD No.3 Evo Photometer tablet reagents
- 250 Phenol Red Photometer tablet reagents
- Instruction manual for Water Safety Kit Basic
- Instruction manual for DI 20 Incubator
- Case

Tintometer GmbH

Lovibond® Water Testing
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
sales@lovibond.com
www.lovibond.com
Duitsland

The Tintometer Limited

Lovibond House
Sun Rise Way
Amesbury, SP4 7GR
Tel.: +44 (0)1980 664800
sales@lovibond.uk
www.lovibond.com
Verenigd Koninkrijk

Tintometer China

9F, SOHO II C.
No.9 Guanghualu,
Chaoyang District,
Beijing, 100020
Customer Care China Tel.: 4009021628
Tel.: +86 10 85251111 Ext. 330
Fax: +86 10 85251001
chinaoffice@tintometer.com
www.lovibond.com
China

Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,
Lebu Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com
www.lovibond.com
Maleisië

Tintometer Brazil

Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@lovibond.us
www.lovibond.com.br
Brazilië

Tintometer Inc.

6456 Parkland Drive
Sarasota, FL 34243
Tel: 941.756.6410
Fax: 941.727.9654
sales@lovibond.us
www.lovibond.us
Verenigde Staten

Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2nd, 3rd & 4th Floor
Sanathnagar Industrial Estate,
Hyderabad, 500018
Telangana
Tel: +91 (0) 40 23883300
Toll Free: 1 800 599 3891/ 3892
indiaoffice@lovibond.in
www.lovibondwater.in
India

Tintometer Spain

Postbox: 24047
08080 Barcelona
Tel.: +34 661 606 770
sales@tintometer.es
www.lovibond.com
Spanje