

# Lovibond® Colour Measurement Tintometer® Group



## Lovibond® RYBN-kleurcode 2,0R 7,0Y 0,5N (5¼)

Bestelnr.: 134100



Reference liquid standards are assigned a nominal value and Lovibond Tintometer always endeavours to match this value as closely as possible. On occasion, it may not be possible to match the requested value exactly. The value achieved and the expected performance tolerances will, of course, always be reported.

### Lovibond® RYBN-kleurcode 2,0R 7,0Y 0,5N (5¼)

500-ml-fles, gecertificeerd volgens QS-systeem ISO 9001 met een houdbaarheid van 12 maanden.

Volledige traceerbaarheid volgens internationaal erkende normen

Elke fles wordt geleverd met volledige certificaten, inclusief veiligheidsinformatieblad.