

Lovibond® Colour Measurement Tintometer® Group



Vloeibare standaard platinakobalt-kleurcode 100 - 500 ml

Bestelnr.: 134190



Reference liquid standards are assigned a nominal value and Lovibond Tintometer always endeavours to match this value as closely as possible. On occasion, it may not be possible to match the requested value exactly. The value achieved and the expected performance tolerances will, of course, always be reported.

Vloeibare standaard platinakobalt-kleurcode 100 - 500 ml

500-ml-fles, gecertificeerd volgens QS-systeem ISO 9001 met een houdbaarheid van 12 maanden.

Volledige traceerbaarheid volgens internationaal erkende normen

Elke fles wordt geleverd met volledige certificaten, inclusief veiligheidsinformatieblad.