

# Lovibond® Water Testing Tintometer® Group



## DBNPA Test Kit

For the standard situation in water testing



- For testing cooling water and other industrial applications
- Test equipment with Checkit® Comparator
- Reagents for up to 100 tests
- Easy to use & result evaluation with visual comparison

Номер заказа: 56K701190

### A fast biocide

DBNPA (dibromonitrile propionamide) is a fast-acting biocide that rapidly degrades to non-hazardous by-products, including ammonia, bromide ions, dibromoacetonitrile and dibromoacetic acid. It therefore also meets strict environmental requirements. It is one of the biocides recommended by many guidelines for the treatment of industrial process water.

### Reliably eliminating microorganisms

In water, DBNPA kills microorganisms quickly and acts reliably against bacteria, fungi, yeasts, algae and microbiological slime: it can quickly penetrate the microbial cell membrane, act on certain protein genes and thus lead to cell death.

### Test Kit with comparator

The application of the Test Kit is very easy. After pre-filtration of the water sample, a certain number of reagents is added step by step. Further filtration and the addition of reagent tablets follows. The colour change is used to interpret the results with the help of the colour disc in the comparator.

### Промышленность

Другие отрасли

### Применение

Контроль дезинфицирующих средств

### DBNPA Test Kit

The use of DBNPA in cooling water treatment is widespread because it is a fast and reliable biocide. It is used to combat microorganisms that can cause damage to the water system and material. To ensure that it can work efficiently, regular checks of the DBNPA content should be carried out. This test kit is a reliable companion for these routine tasks.

### Диапазон измерений

Test Name

DBNPA

Диапазон измерений

0 - 6.8 mg/L DBNPA

## Объем поставки

- Plastic case
- CHECKIT® Comparator with standard optics
- CHECKIT® Disc DBNPA
- 30 tablet reagents DPD No.1 rapid
- 30 tablet reagents DPD No.3 rapid
- 2 Plastic tube with lid, 10 mL
- Instruction manual

Заголовок	Номер заказа
Компаратор CHECKIT®	145000
Пластиковая кювета, глубина слоя 13,5 мм, с серой крышкой, комплект из 10 штук	145500
Пластиковая кювета, глубина слоя 13,5 мм, с серой крышкой, комплект из 5 штук	145505
Пластиковая кювета, глубина слоя 13,5 мм, с серой крышкой, комплект из 100 штук	145510
3/150	255020
Шприц, пластик, 20 мл	56A006501

### Tintometer GmbH

Lovibond® Water Testing  
Schleefstraße 8-12  
44287 Dortmund  
Tel.: +49 (0)231/94510-0  
sales@lovibond.com  
www.lovibond.com  
Германия

### The Tintometer Limited

Lovibond House  
Sun Rise Way  
Amesbury, SP4 7GR  
Tel.: +44 (0)1980 664800  
sales@lovibond.uk  
www.lovibond.com  
ВЕЛИКОБРИТАНИЯ

### Tintometer China

Room 1001, China Life Tower  
16 Chaoyangmenwai Avenue,  
Beijing, 100020  
Customer Care China Tel.: 4009021628  
Tel.: +86 10 85251111 App. 330  
Fax: +86 10 85251001  
chinaoffice@tintometer.com  
www.lovibond.com  
Китай

### Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,  
Lebuh Nilam 2, Bandar Bukit Tinggi,  
Klang, 41200, Selangor D.E  
Tel.: +60 (0)3 3325 2285/6  
Fax: +60 (0)3 3325 2287  
lovibond.asia@tintometer.com  
www.lovibond.com  
Малайзия

### Tintometer Brazil

Caixa Postal: 271  
CEP: 13201-970  
Jundiaí – SP  
Tel.: +55 (11) 3230-6410  
sales@lovibond.us  
www.lovibond.com.br  
Бразилия

### Tintometer Inc.

6456 Parkland Drive  
Sarasota, FL 34243  
Tel: 941.756.6410  
Fax: 941.727.9654  
sales@lovibond.us  
www.lovibond.us  
США

### Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2<sup>nd</sup>, 3<sup>rd</sup> & 4<sup>th</sup> Floor  
Sanathnagar Industrial Estate,  
Hyderabad, 500018  
Telangana  
Tel: +91 (0) 40 23883300  
Toll Free: 1 800 599 3891/ 3892  
indiaoffice@lovibond.in  
www.lovibondwater.in  
Индия

### Tintometer Spain

Postbox: 24047  
08080 Barcelona  
Tel.: +34 661 606 770  
sales@tintometer.es  
www.lovibond.com  
Испания