

# Lovibond® Water Testing Tintometer® Group



## Water Safety Kit Microbiology

Mobile microbiological test equipment for demanding situations



- Testing coliforms and E. coli bacteria according to WHO standards
- DI 20 with heating and cooling as standard
- Also available as a duo version

Ürün No: 56K681254

You need certainty and reliability for test results in emergency situations or outback in areas without laboratories. The Water Safety Kit Microbiology contains all the equipment necessary to perform the microbiological analysis component for safe drinking water. Due to the heating and cooling function of the DI 20 incubator you can test different sorts of bacteria. You can perform filtration with sterilisable equipment including sampler, filter flask & vacuum pump on-site - sterilisation equipment not included. A backpack is available as an option for this kit to allow ease of transportation.

### Sanayi

Diğer Sanayiler | Sivil Toplum Kuruluşları

### Uygulama

İçme Suyu Arıtma

### Water Safety Kit Microbiology

Take all test equipment with you in one hand or on your back in every situation. This kit includes all you need to determine bacteria directly on-site and monitor the quality of potable water even when you don't have any laboratory nearby. The DI 20 incubator with its unique capability of heating and cooling as standard ensures reliable results.

### User

Non Profit Organizations (NGO's) | Medical research and development | Institutes, Universities, Schools | World Health Organization (WHO)

### Teknik Veri

Boyutlar	487 x 229 x 386 mm
Ambalajlı Ağırlık	14.75 kg

# Teslimat Kapsam Metni

- Marker Pen
- 2.5" Bailer with Rope
- Filter Flask, 500 ml
- Tweezers
- Pump, 35 ml stroke
- 5 Universal Container and spoon, 30 ml
- DB 513 Shaker Tube with lid, 100 ml
- 150 EZ Pack Filters
- 5 Pipette Disposable, 1.5 ml (EDT222)
- DI 20 Incubator
- 200 Media Pads and dispenser
- Incubator Battery Lead/Acid
- Filtration Device
- 50 Petri Dish for DI 20 Incubator
- Lead Set IEC C13, (3 pin)
- Membrane Lauryl Sulphate Broth, 40 g
- Instruction manual for Water Safety Kit Microbiology
- Instruction manual for DI 20 Incubator
- Case

Başlık	Ürün No
Güç kablosu seti	136300
Incubator Battery, Lead/Acid	192345
Şişeler	363590
Universal container with spoon and lid, 30 ml	366239
Plastik kap, 100 ml derecelendirmeli, kapaklı	385130
Petri kapları, 10 adet.	400855
Su Güvenlik Kitleri için sırt çantası	56A01400
Media Pads and Dispenser (Pack of 200)	56A016330
Ölçü kabı	56A018501
Diyaframlı filtreleme	56A023801
Cımbız	56A604201
Diyaframlı filtreleme	56B000710
D004 TTC/MAC daldırma kızakları (Macconkey)	56B010410
D005 TTC/Kromojenik daldırma kızakları	56B010510
Membrane Lauryl Sulphate Broth, 38.1 gr.	56P069740
Pump (35 ml stroke)	56X610301
Pipette, disposable, 1.5 ml	BI-K10060
Marker Pen	BI-K14087

## Tintometer GmbH

Lovibond® Water Testing  
Schleefstraße 8-12  
44287 Dortmund  
Tel.: +49 (0)231/94510-0  
sales@lovibond.com  
www.lovibond.com  
Almanya

## The Tintometer Limited

Lovibond House  
Sun Rise Way  
Amesbury, SP4 7GR  
Tel.: +44 (0)1980 664800  
Fax: +44 (0)1980 625412  
sales@lovibond.uk  
www.lovibond.com  
Birleşik Krallık

## Tintometer China

Room 1001, China Life Tower  
16 Chaoyangmenwai Avenue,  
Beijing, 100020  
Customer Care China Tel.: 4009021628  
Tel.: +86 10 85251111 App. 330  
Fax: +86 10 85251001  
chinaoffice@tintometer.com  
www.lovibond.com  
Çin

## Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,  
Lebuh Nilam 2, Bandar Bukit Tinggi,  
Klang, 41200, Selangor D.E  
Tel.: +60 (0)3 3325 2285/6  
Fax: +60 (0)3 3325 2287  
lovibond.asia@tintometer.com  
www.lovibond.com  
Malezya

## Tintometer Brazil

Caixa Postal: 271  
CEP: 13201-970  
Jundiaí – SP  
Tel.: +55 (11) 3230-6410  
sales@lovibond.us  
www.lovibond.com.br  
Brezilya

## Tintometer Inc.

6456 Parkland Drive  
Sarasota, FL 34243  
Tel: 941.756.6410  
Fax: 941.727.9654  
sales@lovibond.us  
www.lovibond.us  
Amerika Birleşik Devletleri

## Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2<sup>nd</sup>, 3<sup>rd</sup> & 4<sup>th</sup> Floor  
Sanathnagar Industrial Estate,  
Hyderabad, 500018  
Telangana  
Tel: +91 (0) 40 23883300  
Toll Free: 1 800 599 3891/ 3892  
indiaoffice@lovibond.in  
www.lovibondwater.in  
Hindistan

## Tintometer Spain

Postbox: 24047  
08080 Barcelona  
Tel.: +34 661 606 770  
sales@tintometer.es  
www.lovibond.com  
İspanya