

# Lovibond® Colour Measurement Tintometer® Group



## AOCS Tintometer renk endeksi 2,5R 20,0Y (5¼")

Ürün No: 134270



Reference liquid standards are assigned a nominal value and Lovibond Tintometer always endeavours to match this value as closely as possible. On occasion, it may not be possible to match the requested value exactly. The value achieved and the expected performance tolerances will, of course, always be reported.

AOCS Tintometer renk endeksi 2,5R 20,0Y (5¼")  
500 ml şişe, ISO 9001 uyarınca 12 aylık kullanım ömrüyle QE sistemi-  
ne göre sertifikalanmıştır.

Uluslararası kabul görmüş standartlara göre tam takip edilebilirlik  
Her şişe, güvenlik bilgi formu dahil tam sertifikalarla birlikte teslim edi-  
lir.