

# Lovibond® Water Testing

## Tintometer® Group



### Pooltester Multipooltester 5 in 1 in plastic case, Alkalinity-m, Calcium Hardness, Chlorine, Cyanuric acid, pH

The fast assistants for Pool Testing



- Easy to use
- Futuristic design
- RAPID tablets fast dissolving
- Highest accuracy

Part Number: 151900

#### High measuring accuracy

Results are shown with high measurement accuracy. They have a certainty and can counteract if necessary.

#### Intuitive

Self-explanatory is the operation of the tester. And if you are ever unsure:

A look at the packaging and the intuitively understandable instructions are sufficient to eliminate the ambiguity.

#### Ergo mixed

The ergonomically optimized lid ensures a secure grip.

#### Industry

Pool Domestic | Pool Services

#### Application

Disaster Management | Pool Water Control

#### Pooltester Multipooltester 5 in 1 in plastic case, Alkalinity-m, Calcium Hardness, Chlorine, Cyanuric acid, pH

Water care is a must, even for the smallest pool. No pool owner can get a reliable overview of whether the pH value and chlorine content in their water paradise are in order more easily and quickly than with the pooltester from Lovibond®. Simply immerse the water tester or fill a water sample into the chambers, add a fast-dissolving test tablet and see with your own eyes where there is a need for action. The water changes color and the comparison scale makes it easy to read off the values. From pH to chlorine, bromine, active oxygen to copper and calcium hardness, the pooltester offers all the important parameters that every pool owner needs to keep an eye on, depending on the version.

## Measuring Range

Test Name	Measuring Range
Alkalinity-m	20 mg/L CaCO <sub>3</sub>
Calcium hardness	20 - 800 mg/L CaCO <sub>3</sub>
Chlorine	0.1 - 3.0 mg/L Cl <sub>2</sub>
Cyanuric acid	20 - 200 mg/L CyA
pH	6.8 - 8.2 pH

## Delivery Scope

- Pooltester Chlorine - pH LR in a robust plastic case
- Cyanuric acid tube
- Dilution / shaker tube (100 ml)
- Dilution / shaker tube (30 ml)
- Cleaning brush
- Stirring rod
- 20 tablet reagents each DPD No. 1 Rapid, DPD No. 3 Rapid, Phenolred Rapid
- 10 tablet reagents each CyA-Test, Alk-Test, CAL-Test
- Instruction manual

### Tintometer GmbH

Lovibond® Water Testing  
Schleefstraße 8-12  
44287 Dortmund  
Tel.: +49 (0)231/94510-0  
sales@lovibond.com  
www.lovibond.com  
Germany

### The Tintometer Limited

Lovibond House  
Sun Rise Way  
Amesbury, SP4 7GR  
Tel.: +44 (0)1980 664800  
Fax: +44 (0)1980 625412  
sales@lovibond.uk  
www.lovibond.com  
UK

### Tintometer China

9F, SOHO II C.  
No.9 Guanghualu,  
Chaoyang District,  
Beijing, 100020  
Customer Care China Tel.: 4009021628  
Tel.: +86 10 85251111 Ext. 330  
Fax: +86 10 85251001  
chinaoffice@tintometer.com  
www.lovibond.com  
China

### Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,  
Lebu Nilam 2, Bandar Bukit Tinggi,  
Klang, 41200, Selangor D.E  
Tel.: +60 (0)3 3325 2285/6  
Fax: +60 (0)3 3325 2287  
lovibond.asia@tintometer.com  
www.lovibond.com  
Malaysia

### Tintometer Brazil

Caixa Postal: 271  
CEP: 13201-970  
Jundiaí – SP  
Tel.: +55 (11) 3230-6410  
sales@lovibond.us  
www.lovibond.com.br  
Brazil

### Tintometer Inc.

6456 Parkland Drive  
Sarasota, FL 34243  
Tel: 941.756.6410  
Fax: 941.727.9654  
sales@lovibond.us  
www.lovibond.us  
USA

### Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2<sup>nd</sup>, 3<sup>rd</sup> & 4<sup>th</sup> Floor  
Sanathnagar Industrial Estate,  
Hyderabad, 500018  
Telangana  
Tel: +91 (0) 40 23883300  
Toll Free: 1 800 599 3891/ 3892  
indiaoffice@lovibond.in  
www.lovibondwater.in  
India

### Tintometer Spain

Postbox: 24047  
08080 Barcelona  
Tel.: +34 661 606 770  
sales@tintometer.es  
www.lovibond.com  
Spain

Technical changes without notice

Printed in Germany

Lovibond® and Tintometer® are Trademarks of the Tintometer Group of Companies