

Lovibond[®] Colour Measurement Tintometer[®] Group



Liquid Standard Platinum-Cobalt Color 5 - 500 ml

Part Number: 134140



Reference liquid standards are assigned a nominal value and Lovibond Tintometer always endeavors to match this value as closely as possible. On occasion, it may not be possible to match the requested value exactly. The value achieved and the expected performance tolerances will, of course, always be reported.

Liquid Standard Platinum-Cobalt Color 5 - 500 ml
500ml Bottle certified under ISO 9001 quality system with a 12-month shelf life.
Full traceability to internationally recognized standards
Each bottle supplied with full certification including MSDS