

# Lovibond® Colour Measurement Tintometer® Group



## Índice de color RYBN de Lovibond® 2,5R 14,0Y 0,7N (5¼)



Referencia No: 134120

Reference liquid standards are assigned a nominal value and Lovibond Tintometer always endeavours to match this value as closely as possible. On occasion, it may not be possible to match the requested value exactly. The value achieved and the expected performance tolerances will, of course, always be reported.

## Índice de color RYBN de Lovibond® 2,5R 14,0Y 0,7N (5¼)

Frasco de 500 ml, certificado según el sistema QS de conformidad con la norma ISO 9001 con una durabilidad de 12 meses.

Trazabilidad absoluta de conformidad con los estándares reconocidos a nivel internacional

Cada frasco se suministra con todos los certificados, incl. hoja de datos de seguridad.