

# Lovibond® Colour Measurement Tintometer® Group



## Rojo, amarillo, neutro de Lovibond®

Referencia No: 139592



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

### Rojo, amarillo, neutro de Lovibond®

A single glass conformance filter for checking PFXi instruments in Lovibond RYBN scale. This filter is supplied as part of the PFXi 880/F package.

#### Datos técnicos

Glass Standard Values	Approximate Value
	- 0.5R
	- 1.3Y
	- 0.1N

---

Filters	1
---------	---

### Volumen de suministro

- Lovibond 0.5R, 1.3Y, 0.1N filter
- Certificate of conformity